

Product datasheet for **RG222044**

SMARCD3 (NM_003078) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SMARCD3 (NM_003078) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SMARCD3
Synonyms: BAF60C; CRACD3; Rsc6p
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG222044 representing NM_003078
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTCCAGGTCTTCAGCACCCACCCACCGTGGTACAGCGCCCCGGGATGCCGTCTGGAGCCCGGATGC
CCCACCAGGGGGCGCCATGGGCCCCCGGGCTCCCCGTACATGGGCAGCCCCGCGTGCAGCCCGGCT
GGCCCCCGGGCATGGAGCCCCCGCAAGCGAGCAGCGCCCCCGCCGGGCAGAGCCAGGCACAGAGC
CAGGGCCAGCCGGTGCCACCGCCCCCGCGGGAGCCGCAGTGCCAAGAGGAGGAAGATGGCTGACAAA
TCCTCCCTCAAAGGATTCGGGAGCTGGTCCCCGAGTCCCAGGCTTACATGGACCTTTGGCATTGAGAG
GAAACTGGATCAAACCATCATGCGGAAGCGGGTGGACATCCAGGAGGCTCTGAAGAGGCCATGAAGCAA
AAGCGGAAGCTGCGACTCTATATCTCCAACACTTTTAAACCTGCGAAGCCTGATGCTGAGGATTCGACG
GCAGCATTGCCTCCTGGGAGCTACGGGTGGAGGGGAAGCTCCTGGATGATCCCAGCAAACAGAAGCGGAA
GTTCTCTTTCTTCAAGAGTTTGGTCATCGAGCTGGACAAAGATCTTTATGGCCCTGACAACACCTC
GTTGAGTGGCATCGACACCCACGACCCAGGAGACGGACGGCTTCCAGGTGAAACGGCCTGGGGACCTGA
GTGTGCGCTGCACGCTGCTCCTCATGCTGGACTACCAGCCTCCCCAGTTCAAACTGGATCCCCGCATAGC
CCGGCTGCTGGGGCTGCACACACAGAGCCGCTCAGCCATTGTCCAGGCCCTGTGGCAGTATGTGAAGACC
AACAGGCTGCAGGACTCCCATGACAAGGAATACATCAATGGGGACAAGTATTTCCAGCAGGTAACCTCTC
TCCTGAACCCTGGGATCTGGGGAGCAGGCGGAGTTCTCCAGCAGGGGAAGTACCTTACCAAGGTAAC
CCACCCCTCCCTTCTAGATTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG222044 representing NM_003078
Red=Cloning site Green=Tags(s)

MTPGLQHPPTVVQRPGMPGARMPHQGAPMGPPGSPYMGSPA VRPGLAPAGMEPARKRAAPPPGQSQAQSQ
 QGQPVPTAPARSRSARRKMKADKILPQRIRELVPESQAYMDLLAFERKLDQTIMRKRVDIQEALKRPMKQ
 KRKLRLYISNTFNPAPKDAEDSDGSIASWELRVEGKLLDDPSKQKRKFSFFKSLVIELDKDLYGPDNHL
 VEWHRTPTTQETDGFQVKRPGDLSVRCTLLMLDYQPPQFKLDPRLARLLGLHTQSRSAIVQALWQYVKT
 NRLQDSDHKEYINGDKYFQQVTSLLNPGIWGAGVLLQQGKYLYQGNYHPSLLDS


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_003078

ORF Size: 1410 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003078.3 , NP_003069.2
RefSeq Size:	1728 bp
RefSeq ORF:	1413 bp
Locus ID:	6604
UniProt ID:	Q6STE5
Cytogenetics:	7q36.1
Domains:	SWIB
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Multiple alternatively spliced transcript variants have been found for this gene, but the biological validity of some variants has not been determined. [provided by RefSeq, Jul 2008]