

Product datasheet for **RG222225**

OR10T2 (NM_001004475) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OR10T2 (NM_001004475) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: OR10T2
Synonyms: OR1-3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG222225 representing NM_001004475
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGAGGTTTCAACAAAACCACTGTGGTTACACAGTTCATCCTGGTGGGTTTCTCCAGCCTGGGGAGC
TCCAGTGCTGCTTTTTGTCATCTTTCTTCTCTATACTTGACAATCCTGGTGGCCAATGTGACCATCAT
GGCCGTTATTCGCTTCAGCTGGACTCTCCACACTCCCATGTATGGCTTTCTATTCATCCTTTTCAATTTCT
GAGTCTGCTACACTTTTGTATCATCCCTCAGCTGCTGGTCCACCTGCTCTCAGACACCAAGACCATCT
CCTTCATGGCCTGTGCCACCCAGCTGTTCTTTTCTGGCTTTGCTTGACCAACTGCCTCCTCATTGC
TGTGATGGGATATGATCGCTATGTAGCAATTTGTACCCTCTGAGGTACACACTCATATAAACAAGG
CTGGGTTGGAGTTGATTTCTCTCTCAGGACCCACAGTTTCTTTATTGCTTTGGTGGCCACCAACCTCA
TTTGTGACATGCGTTTTTGTGGCCCAACAGGGTTAACTACTTTTCTGTGACATGGCACCTGTTATCAA
GTTAGCCTGCACTGACACCCATGTGAAAGAGCTGGCTTTATTTAGCCTCAGCATCCTGGTAATTATGGT
CCTTTTCTGTTAATTCTCATATCCTATGGCTTCATAGTTAACACCATCCTGAAGATCCCCCTCAGCTGAGG
GCAAGAAGGCCTTTGTACCTGTGCCTCACATCTCACTGTGGTCTTTGTCCACTATGGCTGTGCCTCTAT
CATCTATCTGCGCCCAAGTCCAAGTCTGCCTCAGACAAGGATCAGTTGGTGGCAGTGACCTACACAGTG
GTTACTCCCTTACTTAATCCTCTTGTCTACAGTCTGAGGAACAAGAGGTAACAACTGCATTGAAAAGAG
TTCTTGAATGCCTGTGGCAACCAAGATGAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG222225 representing NM_001004475
Red=Cloning site Green=Tags(s)

MRGFNKTTVVTQFILVGFSSLGELQLLLFVIFLLLYLTLVANVTIMAVIRFSWTLHTPMYGFLLIFSFS
 ESCYTFV IIPQLLVHLLSDTKTISFMACATQLFFFLGFACTNCLLIAMGYDRYVAICHPLRYTLIINKR
 LGLELISLSGATGFFIALVATNLICDMRF CGPNRVNHYFCMAPVIKLACTDTHVKELALFSL SILVIMV
 PFLLLILISYGFIVNTILKIPSAEGKKAFTVCASHLTVV FVHYGCASIIYLRPKSKSASDKDQLVAVTYTV
 VTPLLNLVYSLRNKEVKTALKRVLGMPVATKMS



TRTRPLE - GFP Tag - V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001004475

ORF Size: 942 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004475.1 , NP_001004475.1
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	128360
UniProt ID:	Q8NGX3
Cytogenetics:	1q23.1
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]