

## Product datasheet for **RG222298**

### OR56B4 (NM\_001005181) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OR56B4 (NM_001005181) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	OR56B4
Synonyms:	OR11-67
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG222298 representing NM_001005181 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATACCTCCACCAGTGTACCTATGACTCCAGCCTCCAGATTTCCAGTTCATCCTGATGGGATTAC  
CAGGCATTCATGAGTGGCAGCACTGGCTCTCCCTGCCCTGACTCTGCTCTACCTCTTAGCTCTGGTGC  
CAACCTCCTCATATAATCACCATTCAACATGAGACCGTGCTACATGAACCCATGTACCATTGCTGGGC  
ATATTAGCAGTGGTGGACATTGGCCTGGCCACCACCATCATGCCAAGATCCTGGCCATCTTCTGGTTTG  
ATGCCAAGGCCATTAGCCTCCCATGTGTTTTGCTCAGATCTATGCCATCCACTGCTTCTTCTGCATAGA  
GTCAGGCATCTTTCTCTGCATGGCAGTAGACAGATACATAGCCATCTGTGCCCTCTCAGTACCCCTCC  
ATAGTCACTAAAGCTTTTGTCTTCAAAGCCACAGGGTTCATCATGCTCAGGAATGGCCTGTTGACCATCC  
CAGTGCCTATACTGGCTGCCAGAGACTACTGTTCCAGGAATGAAATCGAGCACTGCCTCTGCTCTAA  
CTTGGGGGTTATCAGCCTGGCTTGATGACATCACTGTGAACAAATTTACCAACTGATGCTAGCATGG  
GTCTTGGTTGGGAGTGATGGCTCTGGTATTTCTTCTATGCTGTAATCCTTCACTCTGTGCTGAGGC  
TGAACCTCAGCAGAAGCAATGTCCAAGGCTCTGAGCACTGTAGCTCCACCTCATCTCATCCTCTTCCA  
CACAGGTATCATTGTGCTGTCTGTACACACCTTGACAGAGAAAAAGATTCCCTTATTCTGTGTTCCCT  
AATGTGCTGCACAATGTCATCCCCCTGCACTCAACCCCTGGCCTGTGCACTCAGGATGCACAACTCA  
GACTGGGCTTTTCCAGAGACTGCTTGGACTGGGTCAGGACGTGTCCAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG22298 representing NM\_001005181  
Red=Cloning site Green=Tags(s)

MDTSTSVTYDSSLQISQFILMGLPGIHEWQHWSLPLTLLYLLALGANLLIITIQHETVLHEPMYHLLG  
 ILAVVDIGLATTIMPKILAIWFDAKAIISLPMCFAQIYAIHCFECIESGIFLCMAVDRYIAICRPLQYPS  
 IVTKAFVFKATGFIMLRNGLLTIPVPILAAQRHYCSRNEIEHCLCSNLGVISLACDDITVKNFYQLMLAW  
 VLVGSDMALVFSSYAVILHSVLRRLNSAEAMSKALSTCSSHLILILFHTGIIVLSVTHLAEKKIPLIPVFL  
 NVLHNVIPPALNPLACALRMHKLRLGFQRLGLGQDVSK


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001005181

**ORF Size:** 957 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001005181.1</a> , <a href="#">NP_001005181.1</a>
<b>RefSeq Size:</b>	960 bp
<b>RefSeq ORF:</b>	960 bp
<b>Locus ID:</b>	196335
<b>UniProt ID:</b>	<a href="#">Q8NH76</a>
<b>Cytogenetics:</b>	11p15.4
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Olfactory transduction
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]