

Product datasheet for **RG222351**

CAMK2G (NM_172170) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CAMK2G (NM_172170) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CAMK2G
Synonyms:	CAMK; CAMK-II; CAMKG; MRD59
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG222351 representing NM_172170
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCACCACCGCCACCTGCACCCGTTTCACCGACGACTACCAGCTCTTCGAGGAGCTTGCAAGGGTG
 CTTTCTCTGTGGTCCGCAGGTGTGTGAAGAAAACCTCCACGCAGGAGTACGCAGCAAAAATCATCAATAC
 CAAGAAATTGTCTGCCCGGGATCACCAGAACTAGAACGTGAGGCTCGGATATGTCGACTTCTGAAACAT
 CCAAACATCGTGCCTCCATGACAGTATTTCTGAAGAAGGGTTTCACTACCTCGTGTGGACCTTGTTA
 CCGGCGGGGAGCTGTTTGAAGACATTGTGGCCAGAGAGTACTACAGTGAAGCAGATGCCAGCCACTGTAT
 ACATCAGATTCTGGAGAGTGTTAACCACATCCACCAGCATGACATCGTCCACAGGGACCTGAAGCCTGAG
 AACCTGCTGTGGCGAGTAAATGCAAGGGTGCCCGCTCAAGCTGGCTGATTTTGGCCTAGCCATCGAAG
 TACAGGGAGAGCAGCAGGCTTGGTTTGGTTTTGCTGGCACCCAGGTTACTTGTCCCTGAGGTCTTGAG
 GAAAGATCCCTATGGAACCTGTGGATATCTGGCCTCGGGGTCATCCTGTATATCCTCCTGGTGGGC
 TATCCTCCCTTCTGGGATGAGGATCAGCACAAGCTGTATCAGCAGATCAAGGCTGGAGCCTATGATTTCC
 CATCACCAGAATGGGACACGGTAACTCCTGAAGCCAAGAATTGATCAACCAGATGCTGACCATAAACCC
 AGCAAAGCGCATCACGGCTGACCAGGCTCTCAAGCACCCGTGGTCTGTCAACGATCCACGGTGGCATCC
 ATGATGCATCGTCAGGAGACTGTGGAGTGTTCGCAAGTTCAATGCCCGGAGAAAACCTGAAGGGTGCCA
 TCCTCACGACCATGCTTGTCTCCAGGAACCTTCTCAGCTGCCAAAAGCCTATTGAACAAGAACTCGGATGG
 CGGTGTCAAGCCACAGAGCAACAACAAAACAGTCTCGTAAGCCAGCCCAAGAGCCCGCCCTTGCGAG
 ACGGCCATGGAGCCACAAACCACTGTGGTACACAACGCTACAGATGGGATCAAGGGCTCCACAGAGAGCT
 GCAACACCACACAAGATGAGGACCTCAAAGTGCAGAAACAGGAGATCATTAAAGATTACAGAACAGCT
 GATTGAAGCCATCAACAATGGGGACTTTGAGGCCTACACGAAGATTTGTGATCCAGGCCTCACTTCTTTT
 GAGCCTGAGGCCCTTGGTAACCTCGTGGAGGGATGGATTTCCATAAGTTTTACTTTGAGAACTCCTGT
 CCAAGAACAGCAAGCCTATCCATACCACCATCCTAAACCCACACGTCACCGTATTGGGAGGACGCAGC
 GTGCATCGCTACATCCGCCTCACCCAGTACATCAACGGGAGGGTGGCCTCGCACCAGCCAGTCAGAA
 GAGACCCGGTCTGGCACCGTGGGATGGCAAGTGGCTCAATGTCCACTACTGCTCAGGGGCCCTG
 CCGCACCGCTGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG222351 representing NM_172170
 Red=Cloning site Green=Tags(s)

MATTATCTRFDDYQLFEELGKGAFSVVRRVCVKKTSTQEYAAKIINTKKLSARDHQKLEREARICRLLKH
 PNIVRLHDSISEEGFHLYVFDLVTGGELFEDIVAREYYSEADASHCIHQILESVNHIHQHDIVHRDLKPE
 NLLLASKCKGAAVKLADFLAIEVQGEQAWFGFAGTPGYLSPEVLRKDPYGKPVDIWACGVILYILLVG
 YPPFWDEDQHKLYQKIKAGAYDFPSPEWDTVTPEAKNLINQMLTINPAKRITADQALKHPWVCQRSTVAS
 MMHRQETVECLRKFNARRKLGAILTTMLVSRNFSAAKSLLNKNSDGGVQPQSNKNSLVSPAQEPAPLQ
 TAMEPQTTVVHNATDGIKGSTESCNTTTEDEDLKVRKQEI IKITEQLIEAINNGDFEAYTKICDPGLTSF
 EPEALGNLVEGMDFHKFFENLLSKNSKPIHTTILNPHVHVGEDAACIAYIRLTQYINGQGRPRTSQSE
 ETRVWHRRDGKWLNVVHYHCSGAPAAPLQ

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_172170

ORF Size: 1554 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172170.2](#), [NP_751910.1](#)

RefSeq Size: 3712 bp

RefSeq ORF: 1557 bp


Locus ID: 818

Cytogenetics: 10q22.2

Protein Families: Druggable Genome, Protein Kinase

- Protein Pathways:** Calcium signaling pathway, ErbB signaling pathway, Glioma, GnRH signaling pathway, Long-term potentiation, Melanogenesis, Neurotrophin signaling pathway, Olfactory transduction, Oocyte meiosis, Wnt signaling pathway
- Gene Summary:** The product of this gene is one of the four subunits of an enzyme which belongs to the serine/threonine protein kinase family, and to the Ca²⁺/calmodulin-dependent protein kinase subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. In mammalian cells the enzyme is composed of four different chains: alpha, beta, gamma, and delta. The product of this gene is a gamma chain. Many alternatively spliced transcripts encoding different isoforms have been described but the full-length nature of all the variants has not been determined.[provided by RefSeq, Mar 2011]

Product images:



Circular map for RG222351