

Product datasheet for **RG222388**

TEX55 (NM_152539) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TEX55 (NM_152539) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TEX55
Synonyms:	C3orf30; TSCPA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG222388 representing NM_152539
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGAAGAGCCTCCGCAAGAGGCTCTGGCTGAACCTTGAACATGAAAGCCCAGCCGCTCCCTCAAGTG
 CTGGCCACACTAAGGGCCAGGAAGAAGACGACCAGAAGAACAGGCCGAAAGGAAGGCAGATAACCACAC
 TGCTCACAGAATAGCTGACCAGACTGCCCTAAGAGTGCCTAGCCAGGCTGAATCCAGCATATTTAGCCAA
 GCTACCAACGGAGTAGCTGAACAAAATGGGCATAGTACACCTGGTCAGGCTGGCCGCAGAGCATCCAACC
 CTGCTGATGTTTCTGACCTTAGAGCAGATGATCAGGTTAATCAAACACCGTCTGAACAGACTAAAGGCAA
 GGCATCTAGCCAAGCTAATAATGTACAGCATGAACAGAGTGATGGTCAGGTGTCTGGCCTGACGGAGGAA
 AGAACTGCTGAACAGACTGAACGAAGATTACCTACCCAGGCTGAGAGAAGAACTTCTGGGCAGATTGATG
 GTAGACTGGCTATGCCATCTGACCAGAGAGGTTCCAGACAGACCGACCACAGAATGGCAGGCCAGTCTGA
 GAGAAGAGCTTCCGAGCAGATGGACCGCAGAATGTCTGGCGAGGCTGAGCGAAGAACTTCTGAGCAGATT
 ACACACAGATTATCCAACTATCTGAGAGAAGACCTTCTGTGCAGATTGACAGTGGGTATCCGTCCTCAT
 CTGACCAAAGTCCTTCTGTACAGATTGACAGTGGATCGTCCGTACCATCTGACCAAAGACCTTCCGTACA
 GATTGACCGCAGAATGTGAGGAAAGTTAGGAGAAGAAGTTCTGAGAAGACTGACTACAGATTGGCTGGC
 CTGGCTGACCCAGGAAGTTCTGAGCAGACTGACCTCAGATTGTATGGCCTCGTTGACCACAAAACATCTG
 TAAAGACTCACCACCAAGTGTACGGCCAAGCCACTGAAGTGAACACCCAGGCTATTGACCAAGCTCA
 TAGTAATGCTGATCAACCTCCAGTTGACAATGCTCACTACACTGAATCTGACCAGACTGACCACTTAGCA
 GACAGACAAGCTAATCATAAAGACCAGCTGTCTTACTATGAAACACGTGGCCAGTCTGAAGACAGAATAT
 TTCCCCAGTTAGGCAACAGCAAAGAGGACAAAGAGGCTGACTACAGAGTACAACCTGCAAAATTTGAGGA
 TAGCCAAGTAGACCTCAATTCCAAGCCTTCAGTTGAAATGGAACTCAGAATGCAACCACTATCCCACCC
 TACAACCCAGTTGATGCCAGATTACCAGTAACCTTCCAAGCAAAAGACCAAGCTCTTTTCCAAGACTCC
 CCTCCATCTTATCCAAATTGAACTATACCAGCAGTCAAGAAAAAACTCAAGCCATAGTAACCAAACTCTGA
 TGAATTTTCAGAAATTGACCAAGGAAAGGGTTATCATATACGCAATCAAACCTTATAGAAGGTTCCCTTCT
 ATAGTTTATGAAGATCCTTACCAAGTTTCACTCCAATACATGGAAAAACACCACATTCTGCAAAATATTC
 AGCAGATTACTGAAAACTTAGTCTATGAAAAGCCAGAGGACCCCTGAATTTTATGCTGTGCCAGGTA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG222388 representing NM_152539
 Red=Cloning site Green=Tags(s)


MEEPPQEALAEPLKHESPAAPSSAGHTKGQEEDDQKNQAERKADNHTAHRIADQTALRVPSQAESSIFSQ
 ATNGVAEQNGHSTPGQAGRRASNPADVSDLRADDQVNPQTPSEQTKGKASSQANNVQHEQSDGQVSGLTEE
 RTAEQTERRLPQAEERRTSGQIDGRLAMPDQGRSRQTDHRMAGQSERRASEQMDRRMSGAEERRTSEI
 THRLSKLSERRPSVQIDSGSSVPSDQSPSVQIDSGSSVPSDQRPVQIDRRMSGKVRRRSSEKTDYRLAG
 LADPGTSEQTDLRLYGLVDHKTSVKTHHQVYQATELAEHQAIDQAHSNADQPPVDNAHYTESDQTDHLA
 DRQANHDKQLSYETRGQSEDRIFPQLGNSKEDKEADYRVQPKFEDSQVDLNSKPSVEMETQNATTIPP
 YNPVDARFTSNFQAKDQALFPRLPISLSKLNYSSEKQTAIVTKSDEFSEIDQKGKGYHIRNQTYRRFPS
 IYVEDPYQVSLQYMEKHILQIFQQITENLVYKPEDPLNFMLCQV

TRTRPLE – GFP Tag – V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:



ACCN:	NM_152539
ORF Size:	1608 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152539.1 , NP_689752.1
RefSeq Size:	1814 bp
RefSeq ORF:	1611 bp
Locus ID:	152405
UniProt ID:	Q96M34
Cytogenetics:	3q13.32

Product images:



Circular map for RG222388