

Product datasheet for **RG222723**

LIMD1 (NM_014240) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIMD1 (NM_014240) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LIMD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG222723 representing NM_014240
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGATAAGTATGACGACCTGGGCTGGAGGCCAGTAAATTCATCGAGGACCTGAACATGTATGAGGCCT
 CTAAGGATGGGCTCTTCGAGTGGACAAGGGTGCAGGCAACAACCCGAGTTTGAGGAACTCGCAGGGT
 GTTCGCCACCAAGATGGCCAAAATCCACCTCCAGCAGCAGCAGCAGCAGCTCCTGCAGGAGGAGACTCTG
 CCCAGGGGAGTAGAGGCCCTGTCAATGGAGGGGGCCGCTGGGCCACAGGCCGTTGGGAAGTTGTGG
 GCAGCAAGCTGACTGTGGATGGTGTGCCAAGCCTCCTTGTGCTCGACAGGGGCACCTGGGCACT
 CACCACCTCGTGTGGCAGCCCCGTACCCACCGCAGGAGCAGAGATCCAGGCCATACCTGCATGGC
 ACGAGGATGGCAGCCAGGACTGTGGTTCCAGGGAGAGCCTGGCGACTTCTGAGATGTCTGCTTTCCACC
 AGCCAGGGCCCTGTGAGGATCCTTCTGCCTCACTCATGGAGACTATTATGACAACCTCTCCTTGGCAAG
 CCCAAAGTGGGGTGACAAACCAGGAGTGTCCCCAGCATCGGCCTGAGTGTAGGGAGTGGGTGGCCTAGC
 TCCCCGGGGAGTGACCCACCACTGCCAAACCTGCGGGGACCATCCCTAAATACCGACAGCTCTCCC
 TGAGCTCCAGCAGGTCTTCTGAGGGTAGCCTCGGTGGTGAGAATAGTGGCATTGGTGGCCGAGCAGCGA
 GAAGCCAAACAGGCCTTTGGTCCACTGCCTCCTCCAGCGGGTGAAGCCTGGCCTGCCTTCCCAAACCTTG
 GAGAACGAGCACCAGCTGTGGGGCCTGTTCCAGCCAGGACTCCTTCTGTGTGACACCCTTGGCCCTGA
 GCTGCCCCAGGCAAGGAGGTCTTCAAGATCAAACCTCGGGGCTGGGGGTGAGGTTTCAGGTGTGATGTC
 CAAACCAATGTGGACCCCAACCTGGTTCAGGATGGGCCAAATCTTACCTTCCAGTTCTGCCCCG
 TCATCCTCGCAGCTGGCCTGGACGGTTCACAGCAGGGTGGGTCCCTGGGCTGGGGCGAAGCCTGGCT
 GCACAGACCTTGGCACTGGTCCCAAGCTCAGCCCACCACTTGTCCATCCAGTGATGCCACCCTGCC
 TGAGTTATCTTGTAAAGAGGGTCCCCTGGGCTGGTCTTCTGATGGTAGCCTGGGATCTGTGCTCCTGGAC
 AGCCCCAGCTCCCTAGGTAAGGCTGCCTGCCAGCCCCTCGTCCCAGGTCTGAGCTGAGACCCTCTG
 CTGCTGAGTTGAAATTAGAAGCCCTCACCAACGTCTGGAGCGAGAGATGGATGCTACCCGAAGGCTGA
 TTACTTTGGAGCCTGTGTGAAATGCAGCAAAGGGTGTGGGGCTGGCCAGGCCTGTCAGGCCATGGGG
 AACCTCTACCATGACACATGCTTACCTGTGCAGCTTGCAGCCGAAGCTGAGAGGAAAAGCCTTTTATT
 TTGTCAACGGCAAAGTGTGGTGAAGAAGACTTCTGTACTCTGGTTTCCAGCAGTCGGCTGACAGGTG
 TTTTCTTTGTGGACATCTGATCATGGACATGATCTGCAAGCCCTGGGGAAGCTTACCACCCCGGCTGT
 TTCCGCTGTGTCTGTATGAGTGTGGATGGGGTGCCTTACCGTGGACTCAGAGAACAAGATCT
 ACTGTGTCCGAGATTACCACAAGGTGCTGGCCCCAAGTGTGCAGCCTGTGGGCTTCCCATCCTTCCACC
 TGAGGGCTCAGATGAGACCATCCGTGTCGTGTCCATGGACAGAGACTACCACGTGGAGTGTACCACTGC
 GAGGACTGTGGTCTGGAGCTCAATGATGAAGATGGCCACCGCTGTTATCCGCTGGAGGACCACCTGTTCT
 GTCACCTGCCACGTGAAGAGGCTGGAGAAGAGACCCTCATCTACAGCCCTTACCAGCACCCTT

ACGCGTACGCGGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG222723 representing NM_014240
 Red=Cloning site Green=Tags(s)


MDKYDDLGLEASKFIEDLNMYEASKDGLFRVDKGAGNNPEFEETRRVFATKMAIHLQQQQQQLQEETL
 PRGSRGPVNGGRLGPQARWEVVGSKL TVDGAAPPLAASTGAPGAVTTLAAGQPPYPPEQSRPYLHG
 TRHGSQDCGSRESLATSEMSAFHQPGCEDPSC LTHGDYYDNL SLASPKWGD KPGVSPSIGLSVGSGWPS
 SPGSDPPLPKPCGDHPLNHRQLSLSSRSSEGLGGQNSGIGGRSSEKPTGLWSTASSQRVSPGLPSPNL
 ENGAPAVGPVQPRTPSVSAPLALSCPRQGLPRNSGLGGEVSGVMSKPNVDPQPFWDGPKSYLSSSAP
 SSSPAGLDGSQQGAVPGLGPKPGCTDLGTGPKLSPTSLVHPVMSTLPELSCKEGLGWSSDGLSGSVLLD
 SPSSPRVRLPCQPLVPGPELRPSAAELKLEAL TQRLEREMDAHPKADYFGACVKCSKGVFGAGQACQAMG
 NLYHDTCTCAACSRKLRGKAFYFVNGKVFCEEDFLYSGFQQSADRCFLCGHLIMDMILQALGKSYHPGC
 FRCVICNECLDGVPTVDSENKIYCVRDYHKVLAPKCAACGLPILPPEGSDEIRVVSMRDYHVECYHC
 EDCGLELNDEDGHRCPLEDHLFCHSCHVKRLEKRPSSSTALHQHHF

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_014240

ORF Size: 2028 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_014240.1, NP_055055.1</u>
RefSeq Size:	5067 bp
RefSeq ORF:	2031 bp
Locus ID:	8994
UniProt ID:	<u>Q9UGP4</u>
Cytogenetics:	3p21.31
Gene Summary:	Adapter or scaffold protein which participates in the assembly of numerous protein complexes and is involved in several cellular processes such as cell fate determination, cytoskeletal organization, repression of gene transcription, cell-cell adhesion, cell differentiation, proliferation and migration. Positively regulates microRNA (miRNA)-mediated gene silencing and is essential for P-body formation and integrity. Acts as a hypoxic regulator by bridging an association between the prolyl hydroxylases and VHL enabling efficient degradation of HIF1A. Acts as a transcriptional corepressor for SNAI1- and SNAI2/SLUG-dependent repression of E-cadherin transcription. Negatively regulates the Hippo signaling pathway and antagonizes phosphorylation of YAP1. Inhibits E2F-mediated transcription, and suppresses the expression of the majority of genes with E2F1-responsive elements. Regulates osteoblast development, function, differentiation and stress osteoclastogenesis. Enhances the ability of TRAF6 to activate adapter protein complex 1 (AP-1) and negatively regulates the canonical Wnt receptor signaling pathway in osteoblasts. May act as a tumor suppressor by inhibiting cell proliferation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG222723