

Product datasheet for **RG222820**

TLR4 (NM_138554) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TLR4 (NM_138554) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TLR4
Synonyms:	ARMD10; CD284; TLR-4; TOLL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG222820 representing NM_138554
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCTCGCGCCTGGCTGGGACTCTGATCCCAGCCATGGCCTTCTCTCCTGCGTGAGACCAGAAA
GCTGGGAGCCCTGCGTGGAGGTGGTTCTAATATTACTTATCAATGCATGGAGCTGAATTTCTACAAAAT
CCCCGACAACCTCCCCTTCTCAACCAAGAACCTGGACCTGAGCTTTAATCCCCTGAGGCATTTAGGCAGC
TATAGCTTCTTCAGTTTCCAGAACTGCAGGTGCTGGATTTATCCAGGTGTGAAATCCAGACAATTGAAG
ATGGGGCATATCAGAGCCTAAGCCACCTCTCTACCTTAATATTGACAGGAAACCCATCCAGAGTTAGC
CCTGGGAGCCTTTTCTGGACTATCAAGTTTACAGAAGCTGGTGGCTGTGGAGACAAATCTAGCATCTCA
GAGAACTCCCCATTGGACATCTCAAACCTTTGAAAGAAGCTTAAATGTGGCTCACAATCTTATCCAATCTT
TCAAATTACCTGAGTATTTTCTAATCTGACCAATCTAGAGCACTTGGACCTTCCAGCAACAAGATTCA
AAGTATTTATTGCACAGACTTGGGGTTCTACATCAAATGCCCTACTCAATCTCTCTTTAGACCTGTCC
CTGAACCCATGAACCTTATCCAACCAGGTGCATTTAAAGAAATTAGGCTTCATAAGCTGACTTTAAGAA
ATAATTTTGATAGTTTAAATGTAATGAAAACCTGTATTCAAGGTCTGGCTGGTTTGAAGTCCATCGTTT
GGTTCTGGGAGAAATTTAGAAATGAAGGAAACTTGGAAAAGTTTGACAAATCTGCTCTAGAGGGCCTGTGC
AATTTGACCATTGAAGAATCCGATTAGCATACTTAGACTACTACCTCGATGATATTATTGACTTATTTA
ATTGTTTGACAAATGTTTCTTCAATTTCCCTGGTGAAGTGTGACTATTGAAAGGGTAAAAGACTTTTCTTA
TAATTTCCGATGGCAACATTTAGAATTAGTTAACTGTAATTTGGACAGTTTCCACATTGAAACTCAAA
TCTCTCAAAGGCTTACTTTCACCTCCAACAAAGGTGGGAATGCTTTTTCAGAAGTTGATCTACCAAGCC
TTGAGTTTCTAGATCTCAGTAGAAATGGCTTGAAGTTTCAAAGTTGCTGTTCTCAAAGTATTTTGGGCTTA
AACAGCCTAAAGTATTTAGATCTGAGCTTCAATGGTGTATTACCATGAGTTCAAACTCTTGGGCTTA
GAACAAC TAGAACATCTGGATTTCCAGCATTCCAATTTGAAACAATGAGTGAGTTTTCAGTATTCTAT
CACTCAGAAAACCTCATTTACCTTGACATTTCTCATACTCACACCAGAGTTGCTTTCAATGGCATCTTCAA
TGGCTTGTCCAGTCTCGAAGTCTTGAAGTGGCTGGCAATTTCTTCCAGGAAAACCTCCTCCAGATATC
TTCACAGAGCTGAGAACTTGACCTTCTGGACCTCTCTCAGTGTCAACTGGAGCAGTTGTCTCCAACAG
CATTTAACTCACTCTCCAGTCTTCAAGTACTAAATATGAGCCACAACAACCTCTTTTTCATTGGATACGTT
TCCTTATAAGTGTCTGAACCTCCCTCCAGGTTCTTGATTACAGTCTCAATCACATAATGACTTCCAAAAA
CAGGAACTACAGCATTTTCCAAGTACTAGCTTTCTTAAATCTTACTCAGAATGACTTTGCTTGTACTT
GTGAACACCAGAGTTTCTGCAATGGATCAAGGACCAGAGGCAGCTCTTGGTGGAAAGTTGAACGAATGGA
ATGTGCAACACCTTCAGATAAGCAGGGCATGCCTGTGCTGAGTTTGAATATCACCTGTCAGATGAATAAG
ACCATCATTGGTGTGTCGGTCTCAGTGTGCTTGTAGTATCTGTTGTAGCAGTTCTGGTCTATAAGTTCT
ATTTTACCTGATGCTTCTTGTGGCTGCATAAAGTATGGTAGAGGTGAAAACATCTATGATGCCTTTGT
TATCTACTCAAGCCAGGATGAGGACTGGGTAAGGAATGAGCTAGTAAAGAATTTAGAAGAAGGGGTGCT
CCATTTCACTCTGCCTTCACTACAGAGACTTTATCCCGGTGTGGCCATTGCTGCCAACATCATCCATG
AAGGTTCCATAAAAAGCCGAAAGGTGATTGTTGTGGTGTCCAGCACTTCCATCCAGAGCCGCTGGTGTAT
CTTTGAATATGAGATTGCTCAGACCTGGCAGTTTCTGAGCAGTCTGCTGGTATCATCTTCAATGTCTG
CAGAAGGTGGAGAAGACCCTGCTCAGGCAGCAGGTGGAGCTGTACCGCCTTCTCAGCAGGAACACTTACC
TGGAGTGGGAGGACAGTGTCTGGGCGGCACATCTTCTGGAGACGACTCAGAAAAGCCCTGCTGGATGG
TAAATCATGGAATCCAGAAGGAACAGTGGGTACAGGATGCAATGGCAGGAAGCAACATCTATC

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG222820 representing NM_138554
 Red=Cloning site Green=Tags(s)


```

MSASRLAGTLIPAMAFLSCVRPESWEPCVEVVPNITYQCMELNFYKIPDNLPFSTKNLDLSFNPLRHLGS
YSFFSFPELQVLDL SRCEIQTI EDGAYQSLSHLSTLIL TGNPIQSLALGAFSGLSSLQKL VAVETNLASL
ENFPIGHLKTLKELNVAHNL IQSFKLPEYFSNL TNLEHLDLSSNKIQSIYCTDLRVLHQMPLLNLSL DLS
LNPMNFIQPGAFKEIRLHKL TLRNNFDSL NVMKTCIQGLAGLEVHRLVLGEFRNEGNLEKFDKSALEGLC
NLTIEEFRLAYLDYYLDDIIDLFNCLTNVSSFSLVSVTI ERVKDFSYNFGWQHLELVNCKFGQFP TLK LK
SLKRLTFTSNKGGNAFSEVDLP SLEFLDL SRNGLSFKGCCSQSDFGTTSLKYLDLSFN G VITMSSNFLGL
EQLEHLDFQHSNLKQMSEFSVFLSLRNLIYLDISHTHTRVAFNGIFNGLSSLEVLK MAGNSFQENFLPDI
FTEL RNLTFLDLSQCQLEQLSPTAFNSLSSLQVLNMSHNNFFSLDTPYKCLNSLQVLDYSLNHIMTSKK
QELQHFPSSLAFLNLTQNDFACTCEHQSF LQWIKDQRQLLVEVERMECATPSDKQGM PVL SLNITCQM NK
TIIIGVSVLSVLVVSVAVLVYKFYFHLMLLAGCIKYGRGENIYDAFVIYSSQDE DWVRNELVKNLEEGVP
PFQLCLHYRDFIPGVAIAANI IHEGFHKSRKVI VVVSQHF IQSRWCIFEY EIAQTWQFLSSRAGIIFIVL
QKVEKTL LRQQVELYRLLSRNTYLEWEDSVLGRHIFWRRLRKALLDGKSWNPEGTVGTGCNWQEATSI
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-Mlul

Cloning Scheme:



ACCN: NM_138554

ORF Size: 2514 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_138554.2](#), [NP_612564.1](#)

RefSeq Size: 5503 bp

RefSeq ORF: 2520 bp

Locus ID: 7099

UniProt ID: [O00206](#)


Cytogenetics: 9q33.1

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Pathogenic Escherichia coli infection, Toll-like receptor signaling pathway

Gene Summary: The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. In silico studies have found a particularly strong binding of surface TLR4 with the spike protein of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), the causative agent of Coronavirus disease-2019 (COVID-19). This receptor has also been implicated in signal transduction events induced by lipopolysaccharide (LPS) found in most gram-negative bacteria. Mutations in this gene have been associated with differences in LPS responsiveness, and with susceptibility to age-related macular degeneration. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2020]

Product images:



Circular map for RG222820