

Product datasheet for **RG222864**

LEMD1 (NM_001001552) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: LEMD1 (NM_001001552) Human Tagged ORF Clone

Tag: TurboGFP

Symbol: LEMD1

Synonyms: CT50; LEMP-1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG222864 representing NM_001001552
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGGATGTGAAGTGTCTGAGTGACTGTAAATTGCAGAACCACTTGAGAAGCTTGGATTTTCACCTG
GCCCAATACTACCTTCCACCAGAAAGTTGTATGAAAAAAGTTAGTACAGTTGTTGGTCTCACCTCCCTG
TGACCACCTGTGATGAATGGACCCAGAGAGCTGGATGGAGCGCAGGACAGTGATGACAGCGAAGGTGGG
CTGCAAGAGCACCAAGCACCAGAATCACATATGGGACTATCACCAAGAGAGAGACTACTGCGCGGAAGA
CCAGACTATCGAGAGCTGGAGAGAAGAAGTTTCCAGTGGGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG222864 representing NM_001001552
Red=Cloning site Green=Tags(s)

MVDVKCLSDCKLQNQLEKLGFSGPILPSTRKLYEKKLVQLLVSPPCAPPVMNGPRELDGAQDSDSEGG
LQEHQAPESHMGLSPKRETTARKTRLSRAGEKKVSQWA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001001552

ORF Size: 324 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001001552.5](#)

RefSeq Size: 779 bp

RefSeq ORF: 327 bp

Locus ID: 93273
 UniProt ID: [Q68G75](#)
 Cytogenetics: 1q32.1

Product images:



Circular map for RG222864