

Product datasheet for **RG222888**

EME1 (NM_152463) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EME1 (NM_152463) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EME1
Synonyms:	MMS4L; SLX2A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG222888 representing NM_152463
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTCTAAAGAAGTCATCACCTCACTGGATTCTGGTGATAGTGACTCTGAGGAGTTGCCAACATTTG
 CCTTTCTGAAGAAGGAACCATCTTCAACAAAGAGGAGACAGCCTGAAAGGGAAGAGAAGATTGTAGTGTT
 TGACATCTCAGATTGTGAAGCCTCCTGTCTCCAGCACAGAGTTATTTTACCACCTGTCCCAGAAATA
 GCTGAAACTGTCACACAAACACAGCCAGTCAGGTTGCTAAGCAGTGAAAGTGAAGATGAAGAAGAATTTA
 TTCCTCTGGCTCAAAGGCTTACATGTAAGTTTCTGACCCACAAGCAACTGAGCCCTGAGGACTCTAGCTC
 CCCAGTTAAAAGTGTGGATCATCAAATAATGAAGTGCATCATGTGACTGGAAAAAGCCCTTTCCA
 AAGATCCCTGAAGTTCCCTCCATGATACCCAGAGAGGAGTGCAGCAGATAACAAGGACCTGATCTTAG
 ATCCATGCTGTCAGCTTCCAGCCTACCTGTCTACCTGCCCTGGCCAGAGCAGCAGCTTGGCAGTAACCAA
 AACAAATTCTGACATCCTTCCACCCAGAAGAAAACCAAGCCGAGTCAGAAGGTCCAGGGAAGAGGCTCA
 CACGGATGCCGGCAGCAGAGACAAGCAAGGCAGAAGGAAAGCACCTGAGAAGACAGGAAAGAAAGAAATG
 CAGCACTGGTTACCAGGATGAAAGCCCAGAGGCCAGAGGAATGCTTAAAACACATCATTGTAGTGCTGGA
 TCCAGTGCTTTACAGATGGAAGGTGGGGGCCAGCTCCTAGGAGCACTGCAGACCATGGAGTGCCGCTGT
 GTGATTGAGGCGCAGGCTGTGCCCTGTCAGTGTCACTTGGAGGAGAAGGGCTGGGCCGCTGAGGACAGAG
 AGGACTGGGTGGAGGAGCCAACAGTACTGGTGTGCTCCGGGCAGAGGCATTTGTGTCCATGATCGACAA
 TGGAAAGCAGGGAAGCCTGGACAGCACTATGAAAGGGAAGGAAACGCTTCAGGGCTTTGTAAGTACATC
 ACAGCAAAGACAGCAGGGAAGCTCTGTCACTGGTGATTGTGGATCAGGAGAAATGCTTCAGTGTCTAGA
 ATCCTCCAAGAAGAGGGAACAGGGAGCAAATAAACAGACCAAGAAGCAGCAGCAGACAGCAACAGAGGC
 CAGCATAGGGTCCATGGTATCCAGGGTAGACGCTGAAGAGGCATTGGTGGATCTGCAGCTACACACAGAA
 GCCCAGGCTCAAATTGTGCAGAGCTGAAAGAGCTGGCCGACTTCACATGCGCATTACACAAAGGCTGTGG
 CTGAGGCGCCCTTCAAGAAGCTCCGAGATGAAACTACCTTCTCCTTCTGTCTGGAGAGTGACTGGGCTGG
 AGGGGTGAAGGTGGACCTTGCTGGCAGGGGACTCGCACTAGTCTGGAGGAGACAGATTACAGCAGTGAAC
 CGAGTCAGCCTGGAAATGGCCAGTGCAGTTGTGAATGCCATCCCTCCCCACAGCTCCTGGTACAGGCTT
 ATCAGCAGTGTTCGATAAAGAACGCCAGAATTTGCTCGCAGACATACAGGTGCGCCGTGGGGAAGG
 TGTGACATCCACTTCTCGCCGATTGGACCAGAATATCCAGGCGTATCTACCTTCAGATGACCCTTTA
 CAGCCACATCTCTTTAGATAGTGCTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG222888 representing NM_152463
 Red=Cloning site Green=Tags(s)

MALKKSSPSLDSGDSSEELPTFAFLKKEPSSTKRRQPEREEKIVVVDISDCEASCPPAPELFSPPVPEI
 AETVTQTQPVRLLSSESEDEEEFIPLAQRLTCKFLTHKQLSPEDSSSPVKSVDLHQNNEGASCDWKKPFP
 KIPEVPLHDTPERSAADNKDLILDPCQLPAYLSTCPGQSSSLAVTKNSDILPPQKTKPSQKVQGRGS
 HGCRQQRQARQKESTLRRQERKNAALVTRMKAQRPEECLKHIIIVLDPVLLQMEGGGQLLQALQTMERC
 VIEAQAVPCSVTWRRRAGPSEDREDWVEEPTVLLRAEAFVSMIDNGKQGLDSTMKGKETLQGFVTDI
 TAKTAGKALSLVIDQEKCFSAQNPPRRGKQGANKQTKKQQRQPEASIGSMVSRVDAEEALVDLQLHTE
 AQAQIVQSWKELADFTCAFTKAVAEAPFKLRDETTFSFCLESDWAGGVKVDLAGRGLALVWRRQIQQLN
 RVSLEMASAVVNAYSPQLLVQAYQQCFSDKERQNLADIQVRRGEGVTSTSRRIQPELRRRIYLQMTTL
 QPHLSLDSAD

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152463.4
RefSeq Size:	2332 bp
RefSeq ORF:	1713 bp
Locus ID:	146956
UniProt ID:	Q96AY2
Cytogenetics:	17q21.33
Protein Pathways:	Homologous recombination
Gene Summary:	This gene encodes a protein that complexes with methyl methanesulfonate-sensitive UV-sensitive 81 protein to form an endonuclease complex. The encoded protein interacts with specific DNA structures including nicked Holliday junctions, 3'-flap structures and aberrant replication fork structures. This protein may be involved in repairing DNA damage and in maintaining genomic stability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]