

Product datasheet for **RG222924**

CCR10 (NM_016602) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CCR10 (NM_016602) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CCR10
Synonyms: GPR2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG222924 representing NM_016602
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGACGGAGGCCACAGAGCAGGTTTCTGGGGCCATTACTCTGGGGATGAAGAGGACGCATACTCGG
 CTGAGCCACTGCCGGAGCTTTGCTACAAGGCCGATGTCCAGGCCCTCAGCCGGGCCTTCAAACCCAGTGT
 CTCCTGACCGTGGCTGCGCTGGTCTGGCCGGCAATGGCCTGGTCTGGCCACCCACCTGGCAGCCCGA
 CGCGCAGCGCGCTCGCCACCTCTGCCACCTGCTCCAGCTGGCCCTGGCCGACCTTTGCTGGCCCTGA
 CTCTGCCCTTCGCGGCAGCAGGGGCTCTCAGGGCTGGAGTCTGGGAAGTGCCACCTGCCGACCATCTC
 TGGCCTCTACTCGGCCTCCTTCCAGCCGGCTTCTCTTCTGGCCTGTATCAGCGCCGACCGCTACGTG
 GCCATCGCGGAGCGCTCCCAGCCGGCCGCGGCCCTCCACTCCCAGCCGCGCACACTTGGTCTCCGTCA
 TCGTGTGGCTGTGCTCACTGCTCCTGGCGCTGCCTGCGCTGCTCTCAGCCAGGATGGGAGCGGGAAGG
 CCAACGACGCTGTGCGCTCATCTTCCCCGAGGGCCTCACGCAGACGGTGAAGGGGGCAGCGCCGTGGCG
 CAGGTGGCCCTGGGCTTCGCGCTGCCGCTGGGCGTCATGGTAGCCTGCTACGCGCTTCTGGGCCGACGC
 TGCTGGCCGCCAGGGGGCCGAGCGCCGGCGTGCCTGCGCGTGTGGTGGCTCTGGTGGCGGCCTTCGT
 GGTGCTGCAGTGCCTACAGCCTCGCCCTGCTGCTGGATACTGCCGATCTACTGGCTGCGCGGAGCGG
 AGCTGCCCTGCCAGCAAACGCAAGGATGTCGACTGCTGGTACCAGCGGCTTGGCCCTCGCCCGCTGTG
 GCCTCAATCCCCTTCTCTACGCCTTCTGGGCTGCGCTTCCGCCAGGACCTGCGGAGGCTGCTACGGGG
 TGGGAGCTCGCCCTCAGGGCCTCAAACCCGCGCGGCTGCCCGCCGCGCCCGCCTTTCTCTGCTCA
 GCTCCACGGAGACCCACAGTCTCTCTGGGACAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG222924 representing NM_016602
 Red=Cloning site Green=Tags(s)

MGTEATEQVSWGHYSGDEEDAYSAEPLPELCYKADVQAFSRAFQPSVSLTVAALGLAGNGLVLATHLAAR
 RAARSPTSALLQLALADLLLAL TLPFAAAGALQGWSLGSATCRTISGLYSASFHAGFLFLACISADRYV
 AIARALPAGPRPSTPGRAHLVSVIVWLLSLLLALPALLFSQDGGQREGQRRCLIFPEGLTQTVKGASAVA
 QVALGFALPLGVMVACYALLGRLLAARGPERRRALRVVVALVAAFVVLQLPYSLALLLDTADLLAARER
 SCPASKRKDVALLVTSGLALARCLNPVLYAFLGLRFRQDLRLLRGGSSPSGPQPRRGCPRRRLSSCS
 APTETHLSWDN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_016602

ORF Size: 1086 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016602.1 , NP_057686.1
RefSeq Size:	1244 bp
RefSeq ORF:	1089 bp
Locus ID:	2826
UniProt ID:	P46092
Cytogenetics:	17q21.2
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Chemokine signaling pathway, Cytokine-cytokine receptor interaction
Gene Summary:	<p>Chemokines are a group of small (approximately 8 to 14 kD), mostly basic, structurally related molecules that regulate cell trafficking of various types of leukocytes through interactions with a subset of 7-transmembrane, G protein-coupled receptors. Chemokines also play fundamental roles in the development, homeostasis, and function of the immune system, and they have effects on cells of the central nervous system as well as on endothelial cells involved in angiogenesis or angiostasis. Chemokines are divided into 2 major subfamilies, CXC and CC, based on the arrangement of the first 2 of the 4 conserved cysteine residues; the 2 cysteines are separated by a single amino acid in CXC chemokines and are adjacent in CC chemokines. CCR10 is the receptor for CCL27 (SCYA27; MIM 604833); CCR10-CCL27 interactions are involved in T cell-mediated skin inflammation (Homey et al., 2002 [PubMed 11821900]).[supplied by OMIM, Mar 2008]</p>

Product images:



Circular map for RG222924