

Product datasheet for **RG222984**

JPH1 (NM_020647) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	JPH1 (NM_020647) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	JPH1
Synonyms:	CMT2K; JP-1; JP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG222984 representing NM_020647
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGACGGGCGGAAGGTTTCGACTTCGACGATGGCGGCACCTACTGCGGCGGCTGGGAGGAGGCAAGGCGC
 ACGGGCATGGCATCTGCACGGGGCCCAAGGGCCAGGGCGAGTACTCGGGCTCCTGGTCGCACGGCTTCGA
 GGTGGTCGGAGGCTACACCTGGCCAGCGGCAACACCTACCAGGGCTACTGGGCGCAGGGCAAGCGGCAC
 GGGCTGGGGTGGAGACGAAGGGCAAGTGGATGTACCGGGGGAGTGGTCACATGGTTTCAAGGGCGCT
 ACGGGTCCGGCAGAGCCTGTGCACCCCGCTCGCTACGAGGGTACCTGGAGTAACGGGTGCAAGACGG
 GTACGGCGTGGAGACCTACGGGACGGAGGTACCTACCAGGGCCAGTGGCCGGAGGCATGCGGCATGGC
 TACGGCGTGCAGAGCGTGCCTACGGCATGGCCACGGTGATCCGCTACCGCTGCGTACCTCGCTGG
 CCTCGCTGCGCAGCGAGCAGCAATGGCAGCGTGCTCCACGACGCCGAGCCGCCGCGACAGCCCGGC
 CGGCACCCGCGGGGTTTCGTGCTCAACTTCCACGACGCGTGAAGTACGAGGGAAGAAGGGCGGC
 CTCTTCCGGAGGGCTCCCTTCTTGAAGCATGAACTTCGCAAGTCCGAATCCAAGTCTTCCATCTCGA
 GCAAGCGCAGCTCTGTCCGACGACGCGGCCATGAGCAGAATTAGTTCAGCGATGCCAACTCCACGAT
 CAGCTTTGGCGATGTAGATTGTGATTTTGCCCGGTGGAAGACCAGTGGACGCCACCACCACGAAACC
 TACATGGGCGAGTGGAAGAACGACAAGCGCAACGGCTTCGGCGTTAGCGAGCGCTCCAATGGCATGAAGT
 ATGAAGGGGAGTGGGCAATAACAAGAGGCATGGATATGGCTGTACCGTGTTCCTGACGGCTCCAAAGA
 AGAGGGAAAAATACAAAAATAATTTCTGGTCCGTGGGATAAGGAAGCAGCTTATACCAATAAGACATACA
 AAACTAGGGAGAAGGTGGACAGAGCAATTGAAGGCGCCAAAGGGCAGTGCCATGGCCAGAACCAAAG
 TGGAAATAGCAAATTCAGGACTGCACATGCCAGAGCGAAGGCCGATGCCGCCGACCAGGCCGCGTGGC
 CGCTCGCCAGGAGTGCGACATCGCGAGAGCTGTGGCCAGGGAGCTGTACCTGATTTCTACCAACCAGGC
 CCTGATTACGTCAAACAGAGATTTCAAGGAAGGTGTAGATGCTAAAGAAAATCCAGAAGAAAAGGTACCAG
 AAAAGCCACCTACACCAAAGGAGTCTCCTCATTTTATCGCAAAGGCACGACACCCCAAGATCTCCTGA
 GGCAAGTCCCAACACAGCCACTCTCCTGCTTCTCCCAAGCCCTGAAGAAGCAAAACCCAGCTCA
 GGGGCGAGACTCAACCAAGACAAAAGGAGTGTGGCTGATGAGCAGGTGACGGCCATTGTCAATAAGCCCT
 TGATGTCAAAGGCTCCACGAAGGAGGAGGAGCGGTTGTGCCCCAGTCCAAGTACTTGCCCGCCACCA
 CATCCCCAACCCAGTAACGGGGAGCTGCATTCTCAGTATCACGGCTACTACGTGAAGCTGAACGCCCCC
 CAGCACCTCCAGTAGACGTGGAGGACGGCGATGGATCCAGCCAGTCTTCTCAGCACTGGTGCACAAGC
 CATCCGCTAACAAAGTGGAGTCCCTCCAATCTGTGACAAAACAGTTGCCAAAGAAAGCAAAGCTGAGCC
 AAAAGCTAAGAAGTCTGAACTTGCTATACCAAAGAATCCAGCAAGCAACGATTATGCCCTGCTTTGGAA
 AAAGAAGCAATTCAGGCCCTAATTCAATCATGATTGTCTTGTGATGCTGTTGAATATCGGGTTGGCCA
 TTCTTTTGTTCACCTTTCTAACT

ACGCGTACGCGGCGGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG222984 representing NM_020647
 Red=Cloning site Green=Tags(s)


MTGGRFDFFDGGTYCGGWEEGKAHGHGICTGPKGQGEYSGSWSHGFVVGGYTWPSGNTYQGYWAQKGRH
 GLGVETKKGWYRGEWSHGFKGRYGVRSQSLCTPARYEGTWSNGLQDGYGVETYGDDGTYQGQWAGGMRHG
 YGVRQSVYPYGMATVIRSPRLTSLASLRSEQSNGSVLHDAADSPAGTRGGFVLNFHADAELAGKKKGG
 LFRGSLGSMKLRKSESKSSISSKRSSVRSDAAMSRISSSDANSTISFGDVDCDFCPVEDHVDATTET
 YMGWEKNDKRNFGVRSRNGMKYEGEWANNKRHGYGCTVFPDGSKEEGKYKNNILVRGIRKQLIPIRHT
 KTREKVDRAIEGAQRAAAMARTKVEIANSRTAHARAKADAADQAALAAEQECDIARAVARELSPDFYQPG
 PDYVKQRFQEGVDAKENPEEKVPEKPTPKESPHEFYRKGTTTPRSPEASPKHSHSPASSPKLKKQNPSS
 GARLNQDKRSVADEQVTAIVNKPLMSKAPTKEAGAVVPQSKYSGRHHIPNPSNGELHSQYHGYYVKNLAP
 QHPPVDVEDGDGSSQSSSALVHKPSANKWSPSKSVTKPVAKESKAEPKAKKSELAIPKNPASNDSCPALE
 KEANS GPNSIMIVLVMLLNI GLAILFVHFLT

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_020647

ORF Size: 1983 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020647.4](#)

RefSeq Size: 4373 bp

RefSeq ORF: 1986 bp

Locus ID: 56704

UniProt ID: [Q9HDC5](#)


Cytogenetics: 8q21.11

Domains: MORN

Protein Families: Transmembrane

Gene Summary: Junctional complexes between the plasma membrane and endoplasmic/sarcoplasmic reticulum are a common feature of all excitable cell types and mediate cross talk between cell surface and intracellular ion channels. The protein encoded by this gene is a component of junctional complexes and is composed of a C-terminal hydrophobic segment spanning the endoplasmic/sarcoplasmic reticulum membrane and a remaining cytoplasmic domain that shows specific affinity for the plasma membrane. This gene is a member of the junctophilin gene family. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG222984