

Product datasheet for **RG222987**

ZFYVE27 (NM_001002262) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ZFYVE27 (NM_001002262) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | ZFYVE27 |
| Synonyms: | PROTRUDIN; SPG33 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >RG222987 representing NM_001002262 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGACATCAGAACGTGAGGGGAGTGGGCCGGAGCTGAGCCCCAGCGTGATGCCCGAGGCTCCCCTGG
AGTCTCCACCTTTTCTACCAAGTCCCCAGCGTTTGACCTTTTCAACTTGGTTCTCTCTACAAGAGGCT
GGAGATCTACCTGGAACCCTTGAAGGATGCAGGTGATGGTGTTCGATACTTGCTCAGGTGGCAGATGCCT
TTGTGTTCTTGCTGACCTGCCTGGCCCTCAACGTCTTGTTCCTCACTTTGAATGAGGGTGCATGGTACT
CAGTAGGTGCCCTGATGATTTTCAGTGCCCGCCCTGCTGGGCTACCTTCAGGAGGTTTGCCGGGCACGGCT
GCCTGATTCGAGCTGATGCGGAGGAAGTATCATAGCGTGAGGCAGGAGGACCTGCAGAGAGGTCGCCTG
TCTCGTCCCGAGGCCGTGGCTGAGGTGAAGAGCTTCTTGATCCAGCTGGAGGCCCTTCTGAGCCGCCTGT
GCTGCACATGTGAAGCCGCTACCGCGTGTGCACTGGGAGAACCCCGTCGTGTCCTCACAGTTCTATGG
GGCTCTTCTGGGCACAGTCTGCATGCTGTATTTGCTGCCACTCTGCTGGGTTCTCACCTTTTAAACAGC
ACGCTCTTCTGGGGAATGTGGAGTTCTCCGAGTTGTGTCTGAGTACAGGGCATCTCTGCAGCAGAGGA
TGAACCCAAAGCAGGAAGAGCATGCCTTTGAGAGTCCCTCACCACCAGATGTTGGGGGAAGGATGGTCT
GATGGACAGCACGCTGCCCTCACACCCACGGAGACCTCACACCGGCAGCGTGGAGGAGGCTGAGGAG
GCTGAGCCAGATGAAGAGTTTAAAGATGCGATTGAGGAGGATGATGAGGGCCCCCGTCCCAGCAGAGG
ATGAGCTGGCCCTGCAGGACAACGGTTCTGAGCAAGAATGAGGTGCTGCGCAGCAAGGTGTCTCGGCT
CACGGAGCGGCTCCGCAAGCGCTACCCCAACAACCTTCGGAACCTGCACGGGCTGCTCGGCCACCTTC
TCAGTGTGAAGAAGAGGCGGAGCTGCAGTAATTGTGAAACAGCTTCTGCTCTCGATGCTGCTCTTCA
AGGTGCCCAAGTCTCCATGGGGCCACAGCCCTGAAGCCAGAGGGAGACTGTGTTTGTGTGCCTC
GTGTAACCAGACCTTGAGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG222987 representing NM_001002262
 Red=Cloning site Green=Tags(s)

MQTSEREGSGPELSPVMPPEAPLESPPFPPTKSPAFDLFNLVLSYKRLIYLEPLKDAGDGVRYLLRWQMP
 LCSLLTCLGLNLVFLTLNEGAWYSVGMALMISVPALLGYLQEVCRARLPDSELMRKYHYSVRQEDLQRGRL
 SRPEAVAIEVKSFLIQLEAFLSRLCCTCEAAYRVLHWENPVVSSQFYGALLGTVCMLYLLPLCWLTLLNS
 TLFLGNVEFFRVVSEYRASLQQRMPKQEEHAFESPPPPDVGGKGLMDSTPALTPTEDLTPGSVEEAE
 AEPDEEFKDAIEEDDEGAPCPAEDELALQDNGFLSKNEVLRSKVSRLTERLRKRYPTNNFGNCTGCSATF
 SVLKKRRSCSNCGNSFCRCCSFKVPKSSMGATAPEAQRETVFVCASCNQTLISK

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001002262

| | |
|-------------------------------|---|
| ORF Size: | 1212 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001002262.4 |
| RefSeq Size: | 3009 bp |
| RefSeq ORF: | 1215 bp |
| Locus ID: | 118813 |
| UniProt ID: | Q5T4F4 |
| Cytogenetics: | 10q24.2 |
| Protein Families: | Transmembrane |
| Gene Summary: | This gene encodes a protein with several transmembrane domains, a Rab11-binding domain and a lipid-binding FYVE finger domain. The encoded protein appears to promote neurite formation. A mutation in this gene has been reported to be associated with hereditary spastic paraplegia, however the pathogenicity of the mutation, which may simply represent a polymorphism, is unclear. [provided by RefSeq, Mar 2010] |