

Product datasheet for **RG223099**

WTAP (NM_152858) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: WTAP (NM_152858) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: WTAP
Synonyms: Mum2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG223099 representing NM_152858
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCAACGAAGAACCTCTTCCAAGAAGGTTTCGATTGAGTGAAACAGACTTCAAAGTTATGGCAAGAG
ATGAGTTAATTCTAAGATGGAACAATATGAAGCATATGTACAAGCTTTGGAGGGCAAGTACACAGATCT
TAAGTCTAATGATGTAAGTGGCCTAAGAGAGTCTGAAGAAAACTAAAGCAACAACAGCAGGAGTCTGCA
CGCAGGGAAAACATCCTTGTAATGCGACTAGCAACCAAGGAACAAGAGATGCAAGAGTGTACTACTCAA
TCCAGTACCTCAAGCAAGTCCAGCAGCCGAGCGTTGCCCAACTGAGATCAACAATGGTAGACCCAGCGAT
CAACTGTTTTTCTAAAATGAAAGGTGAAGTGAACAGACTAAAGACAACTGGAACAAGCCAAAAT
GAACTGAGTGCCTGGAAGTTTACGCCTGATAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG223099 representing NM_152858
Red=Cloning site Green=Tags(s)

MTNEEPLPKKVRLSETDFKVMARDELILRWKQYEAYVQALEGKYTDLNSNDVTGLRESEEKLNKQQQESA
RRENILVMRLATKEQEMQECTTQIQYLKQVQPSVAQLRSTMVDPAINLFFLKMKGELEQTKDKLEQAQN
ELSAWKFTPDR


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_152858

ORF Size: 453 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152858.2
RefSeq Size:	1742 bp
RefSeq ORF:	456 bp
Locus ID:	9589
UniProt ID:	Q15007
Cytogenetics:	6q25.3
Protein Families:	Druggable Genome
Gene Summary:	The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of certain cellular genes. This gene encodes a WT1-associating protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms. [provided by RefSeq, Jul 2012]