

Product datasheet for **RG223107**

Myelin oligodendrocyte glycoprotein (MOG) (NM_206810) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Myelin oligodendrocyte glycoprotein (MOG) (NM_206810) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MOG
Synonyms:	BTN6; BTNL11; MOGIG2; NRCLP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG223107 representing NM_206810 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAGCTTATCAAGACCCTCTGCCAGCTGCCTCTGCTCCTTCTCCTCCTCCTCCTCCAAG
TGTCTTCCAGCTATGCAGGGCAGTTCAGAGTGATAGGACCAAGACACCCTATCCGGGCTCTGGTCGGGA
TGAAGTGAATTGCCATGTCGCATATCTCTGGGAAGAACGCTACAGGCATGGAGGTGGGTGGTACCGC
CCCCCTTCTCTAGGGTGGTTCATCTCTACAGAAATGGCAAGGACCAAGATGGAGACCAGGCACCTGAAT
ATCGGGCCGGACAGAGCTGCTGAAAGATGCTATTGGTGAGGAAAGGTGACTCTCAGGATCCGGAATGT
AAGTTTCTCAGATGAAGGAGTTTACCTGCTTCTCCGAGATCATTCTTACCAAGAGGAGGCAGCAATG
GAATTGAAAGTAGAAGATCCTTTCTACTGGGTGAGCCCTGGAGTCTGTTCTCCTCGCGGTGCTGCCTG
TGCTCCTCCTGCAGATCACTGTTGGCCTCATCTTCTCTGCTGCACTACAGACTGAGAGGAAAACCTTCG
AGCAGAGATAGAGAATCTCCACCGGACTTTGGGCAATTCCTTGAAGAGCTACTTCCACCTGGAAGCC
CTCTCTGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG223107 representing NM_206810
Red=Cloning site Green=Tags(s)

MASLSRPSLPSCLCSFLLLLLLQVSSSYAGQFRVIGPRHPIRALVGDEVELPCRISPGKNATGMEVGWYR
 PPF SRVVHL YRNGKDQDGDQAPEYRGRTELLKDAIGEGKVTLRIRNVRF SDEGGF TCF FRDHSYQEEAAM
 ELKVEDPFYWVSPGVLVLLAVLPVLLLQITVGLIFLCLQYRLRGKLRAEIENLHRTFGQFLEELLFHLEA
 LSG


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_206810

ORF Size: 639 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_206810.4
RefSeq Size:	1803 bp
RefSeq ORF:	642 bp
Locus ID:	4340
UniProt ID:	Q16653
Cytogenetics:	6p22.1
Protein Families:	Transmembrane
Gene Summary:	<p>The product of this gene is a membrane protein expressed on the oligodendrocyte cell surface and the outermost surface of myelin sheaths. Due to this localization, it is a primary target antigen involved in immune-mediated demyelination. This protein may be involved in completion and maintenance of the myelin sheath and in cell-cell communication. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]</p>