

Product datasheet for **RG223216**

MBNL1 (NM_207296) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MBNL1 (NM_207296) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MBNL1
Synonyms:	EXP; MBNL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG223216 representing NM_207296 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGTTAGTGTACACCAATTCGGGACACAAAATGGCTAACACTGGAAGTATGTAGAGAGTTCCAGA
GGGGGACTTGCTCACGGCCAGACACGGAATGTAATTTGCACATCCTTCGAAAAGCTGCCAAGTTGAAAA
TGGACGAGTAATCGCCTGCTTTGATTTCATTGAAAGGCCGTTGCTCCAGGGAGAAGTCAAATATCTTCAT
CCACCCACATTTAAAAACGCAGTTGGAGATAAATGGACGCAATAACTTGATTCAGCAGAAGAATATGG
CCATGTTGGCCAGCAAATGCAACTAGCCAATGCCATGATGCCTGGTGGCCATTACAACCCGTGCCAAT
GTTTTAGTTGCACCAAGCTTAGCCACCAATGCATCAGCAGCCGCTTTAATCCCTATCTGGACCTGTT
TCTCCAAGCCTGGTCCCGCAGAGATCTTGCCGACTGCACCAATGTTGGTTACAGGGAATCCGGGTGCC
CTGTACCTGCAGCTGCTGCAGCTGCTGCACAGAAATTAATGCGAACAGACAGACTTGAGGTATGTCGAGA
GTACCAACGTGGCAATTGCAACCGAGGAGAAAAATGATTGTCGGTTTGCTCATCTGCTGACAGCACAATG
ATTGACACCAATGACAACACAGTCACTGTGTATGGATTACATCAAAGGGAGATGCTCTCGGAAAAAGT
GCAAATACTTTTCATCCCCTGCACATTTGCAAGCCAAGATCAAGGCTGCCAATACCAGGTCAACAGGC
TGCAGCTGCACAGGCTGCAGCCACCGCAGCTGCCATGTTCCATGGTGCACGGTGTACGCCAGCCACTG
TGTCGGCAGCAACAACATCTGCCACAAGTGTCCCTTCGCTGCAACAGCCACAGCCAACCCAGCCCATCC
TTGATGCTTCTACACTGTTGGGTGCAACATCCTGTCTGTCAGCAGGAAAAATGATACCCATAATATC
TGCCGAACATCTGACTAGCCACAAGTATGTTACCCAGATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG223216 representing NM_207296
 Red=Cloning site Green=Tags(s)


MAVSVTPIRDTKWLTLEVCREFQRGTC SRPDTECKFAHPSKSCQVENGRVIACFDSLKGRCSRENCKYLH
 PPPHLKTQLEINGRNLIQQKNMAMLAQQMQLANAMMPGAPLQPVPMFSVAPSLATNASAAAFNPYLGPV
 SPSLVPAEILPTAPMLVTGNPGVPVPAAAAAAQLMRTDRLEVCREYQRGNCNRGENDCRFAHPADSTM
 IDTNDNTVTVCMDYIKGRCSRECKYFHPPAHLQAKIKAAQYQVNQAAAAQAAATAAMFPWCTVLRQPL
 CPQQHLQVFPVSLQQPQPTSPILDASTLLGATSCPAAAGKMPIIISAHLTSHKYVTQM

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_207296

ORF Size: 1020 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_207296.2
RefSeq Size:	5246 bp
RefSeq ORF:	1023 bp
Locus ID:	4154
UniProt ID:	Q9NR56
Cytogenetics:	3q25.1-q25.2
Gene Summary:	This gene encodes a member of the muscleblind protein family which was initially described in <i>Drosophila melanogaster</i> . The encoded protein is a C3H-type zinc finger protein that modulates alternative splicing of pre-mRNAs. Muscleblind proteins bind specifically to expanded dsCUG RNA but not to normal size CUG repeats and may thereby play a role in the pathophysiology of myotonic dystrophy. Mice lacking this gene exhibited muscle abnormalities and cataracts. Several alternatively spliced transcript variants have been described but the full-length natures of only some have been determined. The different isoforms are thought to have different binding specificities and/or splicing activities. [provided by RefSeq, Sep 2015]