

Product datasheet for **RG223250**

ALDH1A2 (NM_003888) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALDH1A2 (NM_003888) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ALDH1A2
Synonyms:	RALDH(II); RALDH2; RALDH2-T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG223250 representing NM_003888
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACTTCCAGCAAGATAGAGATGCCCGCGAGGTGAAGGCCGACCCCGCCGCTCATGGCGTCGCTGC
 ACCTCCTGCCGTCGCCCACGCCCAATCTCGAAATTAAGTACACCAAGATCTTTATAAACAACGAGTGCCA
 GAACTCAGAGAGTGGGAGAGTGTTCCTGTCTATAATCCAGCCACAGGAGAACAGGTGTGTGAAGTTCAA
 GAAGCAGACAAGGCAGATATAGACAAAGCAGTGCAGGCAGCCGCTGGCTTTCTCTCTTGGTTCAGTGT
 GGAGAAGGATGGATGCTTCAGAAAGGGGACGTCTGTTGGATAAGCTTGACAGCTTGGTGAACGGGACAG
 GGCAGTTCTTGCAACCATGGAATCCCTAAATGGTGGCAAACCATTCCTGCAAGCTTTTTATGTGGATTTG
 CAGGGCGTCATCAAAACCTTTGATATTACGCAGGCTGGGCTGATAAAATTCATGGGATGACCATTCTG
 TAGATGGAGACTATTTTACCTTTACAAGACATGAACCCATTGGAGTGTGTGGACAGATCATCCCATGGAA
 CTTCCCCCTGCTGATGTTTGCCTGGAAATAGCTCCAGCTTTGTGCTGTGGCAATACAGTAGTTATTAAG
 CCAGCAGAGCAAAACCACTCAGTGCCTCTACATGGGAGCCCTCATCAAGGAGGCTGGCTTTCTCCCG
 GGGTCATCAATATTTTGGCAGGATATGGCCAACGGCTGGGGCAGCAATAGCTTCTCACATTGGCATAGA
 CAAGATTGCATTACAGGGTCTACTGAGGTTGGAAAGCTTATCCAAGAAGCAGCTGGAAGAAGTAATTTG
 AAGAGAGTAACCTGGAAGTGGAGGCAAAAGTCTTAATATTATTTTGTGATGCTGACTTGGACTATG
 CTGTGGAGCAGGCCACCAAGGTGTGTTCTTCAATCAAGGTCAGTGTGCACTGCAGGCTCTCGCATCTT
 CGTGGAGGAGTCCATCTATGAGGAGTTTGTGAGAAGAAGCGTGGAGCGGGCAAGAGGCGCGTAGTGGGG
 AGTCCCTTTGACCCCACTGAGCAGGTCCTCAGATTGATAAGAAACAGTACAACAAGATCTTGAAC
 TCATCCAGAGTGGTGTGGCTGAGGGCGCCAAGCTGGAATGTGGAGGCAAGGACTGGGCCGAAAGGGGT
 TTTCAATTGAGCCACAGTGTTCCTCAACGTCAGTGTGATGATGCGGATTGCCAAGGAGGAGATCTTTGGC
 CCTGTTACAGAAATTTTGAATTTAAGACGATGGATGAAGTTATCGAAAGAGCCAATAACTCAGACTTTG
 GACTCGTAGCAGCTGTCTTTACTAATGACATCAACAAGGCCCTCACAGTGTCTTCTGCAATGCAAGCTGG
 GACTGTTTGGATCAATTGTTACAATGCCTTAAATGCCAGAGCCCTTTGGGGATTCAAGATGTCTGGA
 AATGGGAGAGAAATGGGAGAATTTGGCTTGCAGGAGTACTCAGAAGTTAAGACGGTGACAGTAAAGATCC
 CCCAGAAGAACTCC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>RG223250 representing NM_003888
 Red=Cloning site Green=Tags(s)

MTSSKIEMPGEVKADPAALMASLHLLPSPTPNLEIKYTKIFINNEWQNSSESGRVFPVYNPATGEQVCEVQ
 EADKADIDKAVQAARLAFSLGSVWRRMDASERGRLLDKLADLVERDRAVLATMESLNGGKPFQAFYVDL
 QGVIKTFRYAGWADKIHGMITIPVDGDYFTFTRHEPIGVCGQIIPWNFPLLMFAWKIAPALCCGNTVVIIK
 PAEQTPLSALYMGALIKEAGFPVGVINILPGYGPTAGAAIASHIGIDKIAFTGSTEVGKLIQEAAGRSNL
 KRVTLLELGGKSPNIIFADADLDYAVEQAHQGVFFNQGCCTAGSRIFVEESIYEEFVRRSVERAKRRVVG
 SPFDPTTEQGPQIDKKQYNKILELIQSGVAEGAKLECGGKGLGRKGFFIEPTVFSNVTDDMRIAKEEIFG
 PVQEILRFKTMDEVIERANNSDFGLVAAVFTNDINKALTVSSAMQAGTVWINCYNALNAQSPFGGFKMSG
 NGREMGEFGLREYSEVKTVTVKIPQKNS

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_003888

ORF Size: 1554 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003888.4](#)

RefSeq Size: 3398 bp

RefSeq ORF: 1557 bp


Locus ID: 8854

UniProt ID: [O94788](#)

Cytogenetics: 15q21.3

Domains:	aldehyd
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Retinol metabolism
Gene Summary:	This protein belongs to the aldehyde dehydrogenase family of proteins. The product of this gene is an enzyme that catalyzes the synthesis of retinoic acid (RA) from retinaldehyde. Retinoic acid, the active derivative of vitamin A (retinol), is a hormonal signaling molecule that functions in developing and adult tissues. The studies of a similar mouse gene suggest that this enzyme and the cytochrome CYP26A1, concurrently establish local embryonic retinoic acid levels which facilitate posterior organ development and prevent spina bifida. Four transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, May 2011]

Product images:



Circular map for RG223250