

Product datasheet for **RG223403**

SLC2A13 (NM_052885) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC2A13 (NM_052885) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SLC2A13
Synonyms:	HMIT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG223403 representing NM_052885
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGGCGAGCGGCGCAGGAAGCAGCCGGAGCCGGACGCGGCGAGCGCGCGGGGAGTGCAGCCTCTGG
 CTGCCGCCGAATCGAGCACCAGCCTGCAGAGCGCGGGCGCGGGCGGCGGCGTCTGGGGACCTGGAGCG
 CGCGGCGCGGCGCAGTTCCAGCAGGACGAGACCCCGCCTTCGTGTACGTGGTGCCGCTTCTCCGCG
 CTGGGCGGCTTCTGTTGGCTATGACACCGGGGTGGTGTGAGGGCCATGCTGCTCAAGCGGCAGC
 TCAGTCTGGACGCGCTGTGGCAGGAGCTGCTGGTGTCCAGCACGGTGGGGCGGCTGCCGCTCGGCGCT
 GGCCGGAGGCGCCCTCAACGGCGTCTTCGGCCGCCGCGCTGCCATCCTCTGGCCAGTGCCCTCTTACC
 GCCGCTCCGCGGTGCTGGTGCGGCAACAACAAGGAGACACTGCTCGCCGGCCGCTGGTCTGGGAC
 TCGGCATCGGCATTGCTTCTATGACAGTGCCAGTGACATTGCGGAGGTCTCACCACCAATTTAAGAGG
 CCGATTAGTACCATTAAATACCCTCTTCATCACAGGAGGCGAGTCTTTGCAAGTGTGTTGATGGAGCC
 TTCAGTTATCTCCAGAAGGATGGATGGAGGTACATGTTGGGACTTGCARCAAGTCCGGCGGTATACAGT
 TTTTGGCTTTCTCTTTTGCCTGAAAGCCCTCGATGGCTTATTCAGAAAGGACAGACTCAGAAAGCCCG
 TAGAATTTTATCTCAGATGCGTGGTAACCAGACCATTGATGAGGAATATGATAGCATCAAAAACAACATT
 GAAGAGGAGGAAAAAGAGGTTGGCTCAGCTGGACCTGTGATCTGCAGAATGCTGAGTTATCCCCAACTC
 GCCGAGCTTTAATTGTGGGTGTGGCCTACAAATGTTCCAGCAGCTCTCAGGCATTAACACCATCATGTA
 CTACAGTGCAACCATTCTGCAGATGTCTGGTGTGAAGATGATAGACTTGCAATATGGCTGGCTTCAGTT
 ACAGCCTTCACAAATTTCAATTTTACACTTGTGGGAGTCTGGCTTGTGAGAAGGTGGGCGCAGAAAGC
 TTACCTTTGGTAGTTTAGCAGGTACCACGTAAGTCACTATTATCTTGCCTTGGGATTTGTGCTATCAGC
 CCAAGTTTCCCCACGCATCACTTTTAAGCCAATAGCTCCGTCAAGTCAGAACGCCACTTGACACAAGATAC
 AGTTACTGTAATGAATGTATGTTGGATCCAGACTGCGGTTTCTGCTWCAAGATGAACAAATCAACTGTCA
 TTGACTCCTCTGTGTTCCAGTTAATAAAGCATCTACAAATGAGGCAGCCTGGGCGAGGTGTGAAATGA
 AACCAAGTTCAAAACAGAAGATATTTTGGGCTTACAATTTCTGCCCTACTCCATACTCTGGACTGCA
 CTTCTGGGCTTATTTTATCTTGTCTTCTTGCACCTGGAATGGGACCAATGCCTTGGACTGTGAATT
 CTGAAATATATCCCTTTGGGCAAGAAGTACAGGAAATGCATGTTCTGGAATAAACTGGATTTTCAA
 TGTCTGGTTTCACTAACATTTTACACACAGCAGAGTATCTTACATACTATGGAGCTTTCTTCTCTAT
 GCTGGATTTGCTGCTGTGGGACTCCTTTTCTATCTATGGCTGTCTTCTGAGACCAAGGCAAAAAATTAG
 AGGAAATTTGAATCACTCTTTGACAAAGGCTATGTACATGTGGCACTTCAGATTCTGATGAAGGGAGATA
 TATTGAATATATTCGGGTAAGGGAAGTAACTATCATCTTTCTGACAATGATGCTTCTGATGTGGAA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG223403 representing NM_052885
 Red=Cloning site Green=Tags(s)


MGERRRKQPEPDAAASAAGECSLLAAAESSTLSQSAGAGGGVGDLEAARRQFQQDETFAFYVVVAVFSA
 LGGFLFGYDTGVVSGAMLLLRQLSLDALWQELLVSSTVGAAVVSALAGGALNGVFRRAAILLASALFT
 AGSAVLAANNKETLLAGRLVVLGIGIASMTVPVYIAEVSPNLRRLVTINTLFITGGQFFASVVDGA
 FSYLQKDGWRYMLGLAXVPAVIQFFGFLFLPESPRWLIQKGQTQKARRILSQMRGNQTIDEEYDSIKNNI
 EEEEEKEVGSAGPVICRMLSYPPTRRALIVGCGLQMFQQLSGINTIMYYSATILQMSGVEDDRLAIWLASV
 TAFTNFIIFTLVGWVLEKVGRRLTFGSLAGTTVALIILALGFVLSAQVSPRITFKPIAPSGQNATCTRY
 SYCNECMLDPDCGFCXKMNKSTVIDSSCPVNKASTNEAAWGRGENETKFKTEDIFWAYNFCPTPYSWTA
 LLGLILYLVFAPGMGPMPWTVNSEIYPLWARSTGNACSSGINWIFNVLSLTFLLHTAEYLYTYGAFFLY
 AGFAAVGLLIYGLPETKGKLEEIESLFDNRLCTCGTSDSDEGRYIEYIRVKGSNYHLSNDNASDVE

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_052885

ORF Size: 1887 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_052885.1](#), [NP_443117.1](#)


RefSeq Size: 3261 bp

RefSeq ORF: 1947 bp

Locus ID: 114134

UniProt ID:	<u>Q96QE2</u>
Cytogenetics:	12q12
Domains:	sugar_tr
Protein Families:	Transmembrane
Gene Summary:	H(+)-myo-inositol cotransporter (PubMed:11500374). Can also transport related stereoisomers (PubMed:11500374).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG223403