

## Product datasheet for **RG223411**

### Centlein (CNTLN) (NM\_017738) Human Tagged ORF Clone

#### Product data:

|                    |                               |
|--------------------|-------------------------------|
| Product Type:      | Expression Plasmids           |
| Tag:               | TurboGFP                      |
| Symbol:            | Centlein                      |
| Synonyms:          | bA340N12.1; C9orf39; C9orf101 |
| Mammalian Cell     | Neomycin                      |
| Selection:         |                               |
| Vector:            | pCMV6-AC-GFP (PS100010)       |
| E. coli Selection: | Ampicillin (100 ug/mL)        |

**ORF Nucleotide Sequence:** >RG223411 representing NM\_017738  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGCATCGCC

ATGGCGGCGCGTTTCGCTCCCTCACCGCACCCCTTCGCCCCAGCGCGACAGCTGGGCCCCAGGTCCCCAC  
GTGTTGGGCGGGGAGCTGAAGTACACGCAATGCGCAGCGAGGCCTCGGGTTTTGCCGCGCAGCGCGGGA  
GGTGGTCGCGGACGAAAGTGATAAAATCTGGGTGGGTGAAGAAGGGTCAGGGGGCCGGCAGGGCCTGGG  
GGGCGAGCTCCGGCTCATGCTCCCTCCTCAGCGCGCCCATGGGTCCAGACGGCTAGAGGGCATCTCGG  
TAGAGGAGGCGATGGTGACCCGACGACGCTGCTGGAGGAAGAGCTGAGCAGCCTAAAGGAGGAGTTGGC  
CCTGTGTGAGGCTGATAAAGAATTTGTATGGTCTTTGTGGAACGTCCTCAGGTTACAAACCCAGATCTC  
ACACAAGTGGTCAGTTTGGTTGTGGAAGAGAAAAACAGAAATCTGAAGCTAAAGACAGAAAAGTTCTAG  
AAATCTGCAAGTCAAGGATGCCAAAATACAAGAATTTGAACAGAGAGAGTCAGTACTGAAACAGGAAAT  
AAATGACCTTTGTAAGCGGAAATTTGAGTAGATGAAGAAATGCTTTCTTAAGGAAAGAATTCAGTGAC  
TTGGAGAAGAAATTTAAGATAAAAGTCAAGAAATTAAGGACACTAAGGAGTGTGTACAGAACAAGAAG  
AGCAAAACAGACTAGTTATAAAAAATCTGGAGGAGGAAACAGAAATTAAGTACCCGCTGCACTGACCT  
GCTAAATGACCTGGAGAAATTTGAGGAAGCAGGAAGCACATTTGAGAAAAGAAAAATATAGCACTGATGCA  
AAAATAAGACCTTTGAAGACAATTTAATTGAAGCAAGGAAAGAGTTGAAGTATCACAGAGTAAATACA  
ATGCTCTATCATTACAGTTGAGTAATAACAGACTGAAGTTATCCAGAAGGATATGGATATTACCTGGT  
CAGGAAGGAAGTGCAGGAGTGCAGAACTTTACAAACAGAACAGTACACATACAGCCAGCAAGCAGAG  
CTGATCCAGCAGCTTCAGGTTCTCAATATGGACACACAAAAAGTACTGAGAAATCAGGAAGATGTTACA  
CAGCTGAAAGTATATCATATCAAAAATTTACAATGAGTTACATATTTGTTTTGAAACCAAAAAATCAAA  
TGAAGCTATGCTCCGGCAAAGTGTTACTAATCTTCAGGATCAGCTATTACAAAAAGAGCAAGAAAATGCT  
AAGTTAAAAGAAAATTCAGGAATCACAGGGAGCACCTCTTCCTTTACCTCAAGAAAGTGATCCAGACT  
ACTCAGCACAGGTACCTCATCGCCCATCTTATCAAGCTTAGAAACGTTAATGGTTTCACAGAAGTCTGA  
AATTGAGTATTTACAGGAGAACTAAAGATAGCAAAATGAAAACTGTGAGAAAACATATCTGCCAACAAG  
GGTTTCTCCGAAAGAGCATCATGACAAGTGTGAAGGAAAACATAAGGAACCACTGTGAAACGTTCAA



GGTCTTTGTCCCCAAGAGCTCTTTCACAGACTCAGAAGAGCTACAGAAGCTGAGAAAAGCTGAAAGAAA  
GATTGAAAACCTTAGAGAAGGCACTACAACCTAAAGAGCCAAGAAAATGATGAGCTAAGAGATGCCCATGAA  
AAACGCAAGGAACGGCTACAGATGTTACAGACCAACTACAGAGCAGTAAAGAGCAATTAAACAGTGGG  
AAGAAGGCAGTGGCATGACTGAAATCAGGAAAAATAAGAGAGCAGATCCCCAACAACTTCGACAAGAAAGA  
TTCTGACGCTGTGTGGAATGAACTGGCATATTTCAAAAGGGAAAACCAGGAGCTAATGATTCAAAAGATG  
AATCTTGAAGAAGAATTAGATGAACCTTAAAGTACATATATCTATTGATAAGGCAGCAATACAAGAATTGA  
ATAGATGTGTGGCAGAGAGAAGAGAAGAACAGCTCTTTAGATCTGGTGAAGATGATGAGGTCAAGAGGAG  
TACTCCAGAGAAGAAATGAAAAAGAAATGTTGGAGCAGACATTACAGAAGGTCAGTGAAGTGGAAAATCGG  
CTGAAATCTTTTGAGAAAAGGTCGAGAAAATTAAGAAAGGGAATAAAAATTAATGAAAGAAAATGATT  
TTCTGAAATCCCTCTTAAACAGCAACAAGAAGATACAGAGACCAGAGAAAAAGAGCTAGAACAGATAAT  
AAAGGGGAGTAAAGATGTAGAAAAAGAAAACTGAACTTCAAGTAAAAATCAGTGAGCTGGAGACAGAA  
GTCACCTCCCTGAGGAGACAAGTGGCAGAAAGCTAATGCATTGAGAAATGAAAAAGAGCTGATCAACC  
CAATGGAGAATCACACCACTCAGCAGACAGAGCTAAATCCGAGATGGCCACCATGAAAGTGAGATCTGG  
ACGATATGATTGTAAGACAATATGACCAAGGTTAAATTTAAAGCTGCGAAGAAAAATTGCTCTGTGGGT  
CGTCACCACACTGTTCTCAATCATTCCATCAAGGTTATGAGCAATGTGTTTGAGAACCTCAGCAAGGACG  
GCTGGGAGGATGTGAGTGAAAGCAGCAGTATTCTGAAGCAGACCTCTCAAACCTTTGGGAACAATTAT  
TGTAAGAACATCCAGAAAAATAAGTCTACGGAAGATGGAAGAGCCAGAAAAGAAAGTATCCAACAGAA  
GACAGCCAAACACAAGGAAAAGAAATAGTACAGACATATTTAAATATAGATGGCAAGACCCCAAGGACT  
ATTTTCATGATAAGAATGCCAAAAACCACTTTTCAAAGAAGAAATGCAAGATGCAAAAGAGTTTACACA  
TACAGCAGTTCCCTACTAGAGTCAACAGAGAAAAAGTACAAAAATATAACTGCCAGAAATCAAGTAGCAAT  
ATTATTTTATTACGAGAACGGATTATATCCTTGCAACAACAAAAACAGTGTACTTCAGAATGCCAAGAAAA  
CAGCAGAATTGCTGTGTTAAAGAATATAAAGAAGTTAATGAAAAGCTCCTCCATCAACAGCAAGTATCCGA  
TCAACGATTTTCAGACAAGCAGGCAGACAATAAAGAAGCTAAATTTGGATTGGCTGGGCTTCGGAAGAA  
AAAGAAGATTTACTAAAGAAATTTGGAGTCCTCATCTGAAATCACAAAGTTTGGCAGAAAGAAAATCCCAGG  
TAACATTTCCACGGATACAAGTTACATCACTTAGTCCTTCAAGGAGCATGGATTGGAAATGAAGCAATT  
GCAGTATAAACTAAAGAATGCTACTAATGAACTCACTAAACAGTCATCAATGTGAAGACTTTGAAATTT  
GAACTCCTAGCAAAAGAAGAACACATAAAGGAAATGCATGAAAAGATATCTCGAATGGAGAGGGATATAA  
CTATGAAAAGACATTTGATAGAGGACTTGAAATTTGACAGAGAAAGTAAATTTGGAAAGTAACAAGAGTTT  
TAGTGAATGTTACAAAATCTTGATAAAAAAGGTAAGACATTAAGTGAAGAATGTTCCAACAAGAAGGTA  
TCAATTGATTCATAAGCAAGAGACTTAACGTTGCTGTAAAAGAAAAGTCACAGTATGAACAGATGTATC  
AGAAATCTAAAGAGGAGTTAGAAAAAAGGATCTCAAGCTTACTCTCTAGTATCAAGAATAAGTGAGAC  
TGAATCTGCAATGGCAGAAATGAAACAGCAGCATCTAAGCAGCTTCAAGAATTAGCATTGCAAAAGTGAA  
CAGGTCTAGAAAGGTGCACAGAAGACATTGCTGTTAGCCAATGAAAAAGTAGAAGAGTTCACCACATTTG  
TGAAGGCTTTGGCCAAAGAGTTGCAAAATGATGTCATGTGGTAAGGCGACAAATAAGAGAGCTTAAAAA  
AATGAAGAAAAACAGGGACGCTGTAAACCTCAACCCATAAAGCCAGACCTTGGCAGCTTCTATCCTG  
AACATTTACGGTCAGATTTAGAGGAAATATTAGACACAGAAGATCAAGTGGAATTTGAAAAACAAAAA  
TTGATGCTGAAATGACAAGGAATGGATGTTGTACATTCAGAACTTCTTGAAGGACAGCTTCTTTTGC  
CTCATATTTACTAGAAGCAGTACTGAAAAAATAATGAAAAAAGAACTAGTTGAAGGATATTTACACA  
ATTATGAAAGATATTAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG223411 representing NM\_017738

Red=Cloning site Green=Tags(s)


MAARSPSPSPSPARQLGPRSPRVGRGAIEVHAMRSEASGFAGAAREVVADESCKIIVGEEGSGGRRGPG  
GAAPAHAPLLSAPMGSRRLEGISVEEAMVTRTQLLEELSSLKEELALCQADKEFVWSLWKRLQVTNPDL  
TQVVSLLVEREKQKSEAKDRKVLIELQVKDAKIQEFQRESVLKQEIINDLVKRKIAVDEENAFRLKEFSD  
LEKKFKDKSQEIKDTKECVQNKKEQNRLVIKNLEENKKLSTRCTDLLNDLEKLRKQEAHLRKEKYSTDA  
KIKTFEDNLEARKEVEVSQSKYNALSLQLSNKQTEL IQKMDITLVRKELQELQNLQKQNSTHTAQQA  
LIQQLQVLNMDTQKVLNRQEDVHTAESISYQKLYNELHICFETTKSNEAMLRQSVTNLQDQLLQKEQENA  
KLKEKLQESQGAPLPLQESDPDYSAQVPHRPSLSLETLMVSQKSEIEYLQEKLKIANEKLSENISANK  
GFSRKSIMTSAEGKHKEPPVKRSRSLSPKSSFTDSEELQKLKRAERKIEKALQLKSQENDELDAHE  
KRKERLQMLQTNRAVKEQLKQWEEGSGMTEIRKIKRADPQQLRQEDSDAVWNELAYFKRENQELMIQKM  
NLEELDELKVHISIDKAAIQELNRCVAERREEQLFRSGEDDEVKRSTPEKNGKEMLEQTLQKVTELENR  
LKSFEKRSRKLKEGNNKMKENDFLKSLKQQQEDTETREKELEIIGSKDVEKENTELQVKISELETE  
VTSLRQVAEANALRNENEELINPMEKSHQSADRAKSEMATMKVRSGRYDCKTMTKVKFKAANKNSVG  
RHHTVLNHSIKVMSNVFENLSKDGWEDVSESSDSEAQTSQTLGTIIIVETSQKISPTEDGKDQKESDPTE  
DSQTQGKEIVQTYLNIDGKTPKDYFHDKNNAKPTFQKKNCKMQKSSHTAVPTRVNREKYKNITAQSSSN  
IILLRERIIISLQQQNSVLQNAKKTAE SVKEYEYNEKLLHQQVSDQRFQTSRQTIKKLNDLAGLRKE  
KEDLLKLESSEITSLAEENSQVTFPRIQVTSLSPSRSMDEMKQLQYKLNATNELTKQSSNVKTLKF  
ELLAKEEHKEMHEKISRMERDITMKRHLIEDLKFRQKVNLESNKSFSEMLQNLDKKVKTLTEECNNKKV  
SIDSLKQRLNVAVKEKSQYEQMYQKSKEELEKKDLKLTLLVSRISSETESAMAEIETAASKQLQELALQSE  
QVLEGAQKTLLLANEKVEEFTTFVKALAKELQNDVHVRRQIRELKMKKNRDACTSTHKAQTLAASIL  
NISRSDLLEEILDTEQVEIEKTKIDAENDKEWMLYIQKLLEGQLPFASYLLEAVLEKINEKKLVEGYFT  
IMKDIR

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:


NM\_017738

ORF Size:

4218 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_017738.4</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 5522 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 4221 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 54875                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <a href="#">Q9NXG0</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 9p22.2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Gene Summary:</b>          | Required for centrosome cohesion and recruitment of CEP68 to centrosomes.<br>[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                                                                                                                                                                            |

Product images:



Circular map for RG223411