

Product datasheet for **RG223434**

ATP6V1C2 (NM_001039362) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATP6V1C2 (NM_001039362) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ATP6V1C2
Synonyms:	ATP6C2; VMA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG223434 representing NM_001039362 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCGGAGTTTTGGTTAATTTCTGCCCTGGCGATAAGGAAAATTTGCAAGCTCTGGAGAGGATGAATA
CTGTAACCTCCAAGTCCAACCTGTCTTATAATACCAAATTCGCTATTCCTGACTTCAAGGTGGGACCTT
GGATTCCCTGGTTGGCCTCTCTGATGAGTTGGGAAACTCGACACCTTTGCTGAAAGCCTCATAAGGAGA
ATGGCTCAGAGCGTGGTGAAGTCATGGAGGACTCAAAGGGGAAGGTCCAGGAGCACCTCCTGGCAAACG
GAGTTGACTTAACATCCTTTGTGACCCACTTTGAATGGGACATGGCCAAATATCCTGTCAAGCAGCCGCT
CGTGAGTGTGGTGGACACAATAGCCAAGCAACTGGCGCAGATCGAGATGGACCTGAAGTCCCGAACGGCC
GCCTACAACACTCTGAAGACAAACCTGGAGAACCTGGAAGAAATCCATGGGGAACCTTCTCACCCGGA
CACTGAGTGATATTGTGAGCAAAGAGGACTTCGTGCTGGATTCTGAATATCTCGTCACACTTCTGGTCAT
CGTCCCCAAACCAAACTACTCACAATGGCAAAAACCTACGAATCTCTCAGACATGGTGGTCCCTCGA
TCAACCAAACTCATTACTGAGGACAAGGAAGGGGGCCTTTTCACTGTGACTCTGTTTCGAAAAGTGATTG
AAGATTTCAAACCAAGGCCAAAGAAAACAAGTTCCTGTTCTGTAAGAAAGCAACAGTATCAAACCTCCTGTGTT
GCTCTTAAAAAGGGATCATCCACCTTCCCGACCACAAGTTAAGTAACCCCGTAGGTAACCCCTGATA
GGCCTGCTGCGGGCAGACCGACAGAGAGAGAGAGTGAAGGCGAGGGTGAGGGCCCTGCTGCGCTG
GCTCAAGTGAACTTCAGTGAAGCCTTCATTGCCTGGATCCACATCAAGGCCCTGAGAGTGTGTTGGAG
TCCGTGCTCAGGTATGGACTACCAAGTGAACCTCCAGGCAGTGCTCCTGCAGCCGATAAGAAGTCATCCA
CCAAGCGTTTAAAGAGAGTTCTAAACTCTGTCTCCGACATCTGGATGAAGTAGCCGCTACAAGTACT
GGATGCATCTGTGGAGATCCCGGACTGCAACTCAATAACCAAGACTATTTTCTTATGTCTACTTCCAT
ATTGACCTTAGTCTTCTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG223434 representing NM_001039362
 Red=Cloning site Green=Tags(s)

MSEFWLISAPGDKENLQALERMNTVTSKSNLSYNTKFAIPDFKVGTLDSLVLGSDDELGKLDTF AESLIRR
 MAQSVVEVMEDSKGKVQEHL LANGVDLTSFVTHFEWDMAKYPVKQPLVSVVDTIAKQLAQIEMDLKSRTA
 AYNTLKTNLNLEKKS MG NLFTRT L S D I V S K E D F V L D S E Y L V T L L V I V P K P N Y S Q W Q K T Y E S L S D M V V P R
 S T K L I T E D K E G G L F T V T L F R K V I E D F K T K A K E N K F T V R E F Y Y D E K E I E R E R E E M A R L L S D K K Q Q Y Q T S C V
 A L K K G S S T F P D H K V K V T P L G N P D R P A A G Q T D R E R E S E G E G E P L L R W L K V N F S E A F I A W I H I K A L R V F V E
 S V L R Y G L P V N F Q A V L L Q P H K S S T K R L R E V L N S V F R H L D E V A A T S I L D A S V E I P G L Q L N N Q D Y F P Y Y F H
 I D L S L L D

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001039362

ORF Size: 1281 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039362.2](#)

RefSeq Size: 3250 bp

RefSeq ORF: 1284 bp

Locus ID: 245973


UniProt ID: [Q8NEY4](#)

Cytogenetics: 2p25.1

Protein Pathways: Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection

Gene Summary: This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A, three B, and two G subunits, as well as a C, D, E, F, and H subunit. The V1 domain contains the ATP catalytic site. This gene encodes alternate transcriptional splice variants, encoding different V1 domain C subunit isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG223434