

Product datasheet for **RG223524**

H3FI (HIST1H3F) (NM_021018) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H3FI (HIST1H3F) (NM_021018) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: H3FI
Synonyms: H3/i; H3C1; H3C2; H3C3; H3C4; H3C6; H3C8; H3C10; H3C11; H3C12; H3FI; HIST1H3F
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG223524 representing NM_021018
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTCGCACTAAGCAAACAGCTCGTAAGTCCACTGGCGGCAAAGCCCCGCGCAAGCAGCTGGCCACTA
 AGGCGGCCCGCAAAGCGCGCCAGCCACCGGTGGCGTGAAGAAGCCCCACCGCTACAGGCCTGGTACTGT
 CGCCCTCCGTGAAATCCGCCGTATCAGAAATCGACTGAGCTACTGATTCGCAAGCTACCATTCCAGCGT
 CTGGTACGTGAGATCGCGCAGGACTTCAAGACCGACCTGCGCTTCCAGAGCTCGGCTGTGATGGCGCTGC
 AGGAGGCCTGCGAGGCTTACCTGGTGGGCTCTTTGAGGACACCAACCTGTGTGCTATCCACGCCAAGCG
 AGTGACTATCATGCCCAAGGACATCCAGCTCGCTCGCCGATTCCGGGAGAGAGGGCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG223524 representing NM_021018
 Red=Cloning site Green=Tags(s)
 MARTKQTARKSTGGKAPRKQLATKAARKSAPATGGVKKPHRYRPGTVALREIRRYQKSTELLIRKLPFQR
 LVREIAQDFKTLRFQSSAVMALQEACEAYLVGLFEDTNLCAIHAKRVTIMPKDIQLARRIRGERA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_021018

ORF Size: 408 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_021018.2](#), [NP_066298.1](#)

RefSeq Size: 466 bp

RefSeq ORF: 411 bp

Locus ID: 8968
UniProt ID: [P68431](#)
Cytogenetics: 6p22.2
Protein Pathways: Systemic lupus erythematosus
Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3. [provided by RefSeq, Aug 2015]

Product images:



Circular map for RG223524