

Product datasheet for **RG223542**

CCDC96 (NM_153376) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CCDC96
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG223542 representing NM_153376
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACGTCAGCTCTGAGCACACTAAGGACCCCGCGGGGAAGGCGGAGACGGGAAAGCCTGGCCGCGC
 GGCCGTCCAAGATCAAGGCCAGCTCTGGCCCTCCGACCTCTCCGAACCCGGGAGCTGGAGTCGGAGCC
 GGAGGAGGAGGAGGAGGAGCAAGCGGCTTCGCGAGGAGGACCCGCCGCGGACGAGCAGGCCGAGGCCCGG
 AAAGGGCTAACGGCAGCCGAGGCTGCGGGCGAGGAGGGGCTGGAGAGCCGGGAGGCCGCGGAGCCCC
 AGCCCCGAGCCGAAGAGCCGGCCGAGGTCGGGGCTGAGGAGCCGCCAGCCGAGCCCGGAGCCGGGCC
 CGAGGAACTAGAGCGGAGGCGGGGCGGAGGAACTGGAGCAGGCGGAGAGGGGAAGGAAGTCAGGTTT
 CAGGCCTCTCTGCCGCTGACCAGGATCGATGAGGAAGAGGCTGCTGCAGCCCGGAGGCCGAGACAGAGC
 GGGTGGAGGGGAGGAAGAGGACAAGGAGGAGACGAGAGAGCGGCGCGGAGAGCAAGGAGAGGGATGG
 GGAAGGACGCCCGCCAAAAGCCAGGAAGAGGGGAAGCGCTCTATGGAAGGGACGAGTTTGAGGACCTC
 GAGTGGTCCGAGGAGGTCCAGAAGCTGCAGGAGCAGCAGCTGCGCAGCGACCTCTTGACCAAGTACCGTT
 CCCTGCTGGTGGAGCGGAACCGCTCCACGCGCTACAACCTATACCTGCAGCACAAGATCTTCGAGGCGCT
 GCGCAGAAAGAAGGGCTGGAGGCCGCTGAGGTGGTGACCGGGGCGCAGAGGCCGAGGCCCGGAGAAA
 GAGCAAGCGTACCTGCGCCATCTGGCATGCTGGAGGAGCTGAAGAAGCAGCAGGCAGACGACCTGCAGT
 GGTACCACAGGAGCTGGGCCAGTTGAAGCGGAGTGCAGGAGAGAAGCTCACCAGGTGGAGAAGGAGTG
 GCGACGCTTCCAGGCACTCAAGAAGCAGGTGGTATGCAGGCCATGGGAGCTGTGAGTGGAGGCGGG
 CGCCAGGCTGCTCTGCGAGAGGTGGAGCAGATCCAGGCGTTGGAGGATAAAAAGGAGAAGGAGATGAGCG
 CCGTGGGCTGGAGAACATTGAGTGAAGCAGAGCCTGGTGCATTTTGAACAGGATGAGGACCCAGGA
 GGACCTGACCCAGGCTGCTTCTTATTGACTTTGAACAGCTGAAGATTGAGAACCAACCTTCAATGAG
 AAAATTGAGGAACGAAATGAGGAACCTTTAACTTCGAGCAAGGTGACCAACAGTGTGCAAGTAATAA
 CCCACGTGAAGGAAAGCTGCATTCATGGACATGGAGAACGCGTGCAAAAAGACACAGCTGGCAGAAAT
 CGAGGCCAGGCCCGCTAGGAAGAGACATCCTGACCAAGACGAAGCAAGCCGAGAGGGGCTGCGGACG
 GACAACATCAGACTGAATCAGAAGTGGCGCTTCTAGGCAAGGACTCACTTCTTCGGGACTTGGAAGAAA
 AGGTGGACAAGACCAACTGCTTACCAGCGCCTGGAATCCCTGAAGCGCCACCACGCTAGCCTCACTTT
 GTCCTGCAGAGGCGTGAGGCGAGAAGATCAGGGAGGCCAAGCCTTTCTCCCTTCT

AGCGGACCGACGCGTACGCGCGCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG223542 representing NM_153376
 Red=Cloning site Green=Tags(s)

MDVSSEHTKDPGGEGGDGESLAARPSKIKASSGPPTSPEPGELESEPEEEEEQAAASQGGTAADEQAEAP
 KGLTAAEAAGEEGPGEPRPAEPQPEPEPAEVGAEEPAQPEPGAGPEELEAEAGAELEQAAEGKEVRF
 QASLPLTRIDEEEEAAPEATERVEGEEDKEETQRDGAESKERDGEGRPAKSQEEGKRLYGRDEFEDL
 EWSEEVQKLQEQLRSDLLDQYRSLVVERNRSQRYNLYLQHKIFEALRRKKGLEAAEVADRGAEEAPEK
 EQAYLRHLGMLEELKKQADDLQWYHQELGQLKRQCQEKLTRVEKEWRRFQALKKQVVMQAMGSCRMGG
 RQAALREVEQIQALEDKKEKEMSAVRLNIQLKQSLVHFETRMRTQEDLTQGLLLIDFEQLKIENQTFNE
 KIEERNEELLKRSKVTNSVQVITHVKEKLHFMDMENACKKTQLAEIEAQAALGRDILTKTKQAREGLRT
 DNIRLNQKCGLLGKDSLLRDLEEKVDKTELLHRRLESLKRHHASLTLSCRGVQRKIREAKAFLPS

SGPTRRRLE - GFP Tag - V

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_153376

ORF Size: 1665 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_153376.2](#), [NP_699207.1](#)

RefSeq Size: 1769 bp


RefSeq ORF: 1668 bp

Locus ID: 257236

UniProt ID: [Q2M329](#)

Cytogenetics: 4p16.1

Product images:



Circular map for RG223542