

Product datasheet for **RG223570**

EIF2B3 (NM_020365) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EIF2B3 (NM_020365) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EIF2B3
Synonyms:	EIF-2B; EIF2Bgamma
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG223570 representing NM_020365
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAATTTCAAGCAGTAGTGATGGCAGTAGTGGAGGATCTCGGATGACAGACCTAACTCCAGCATTC
 CCAAACCTCTGCTTCCAGTTGGGAACAAACCTTTAATTTGGTACCATTTGAACCTGCTTGAGCGTGTGG
 ATTTGAAGAAGTCATTGTGGTTACAACAGGGATGTTCAAAGGCTCTATGTGCAGAATCAAGATGAAA
 ATGAAGCCAGATATTGTGTGATTCTGATGACGCTGACATGGGAAGTGCAGATTCTTTGCGCTACATAT
 ATCCAAAACCTTAAGACAGATGTGCTGGTCTGAGCTGTGATCTGATAACAGACGTTGCCTTACATGAGGT
 TGTGGACCTGTTAGAGCTTATGATGCATCACTTGCTATGTTGATGAGAAAAGGCCAAGATAGCATAGAA
 CCTGTTCCCGGTCAAAGGGGAAAAAAGCAGTGGAGCAGCGTGACTTCATTGGAGTGGACAGCACAG
 GAAAGAGGCTGCTTTCATGGCTAATGAAGCAGACTTGGATGAAGAGCTGGTCATTAAGGGATCCATCCT
 ACAGAAGCATCCTAGAATACGTTTCCACACGGGCTTGTGGATGCCACCTCTACTGTTTAAAAAATAC
 ATCGTGGATTTCCTAATGGAAAATGGGTCAATAACTTCTATCCGGAGTGAAGTATCCATATTTAGTGA
 GAAAACAGTTTTCTCAGCTTCTCACAACAGGGACAAGAAGAAAAGAGGAGGATCTAAAGAAAAGGA
 GCTGAAGTCCTTAGATATCTACAGTTTTATAAAGAAGCCAATACACTGAACCTGGCTCCCTATGATGCC
 TGCTGGAATGCCTGTCGAGGAGACAGGTGGGAAGACTTGTCCAGATCACAGGTGCGCTGCTATGTCCACA
 TCATGAAAGAGGGGCTCTGCTCTCGAGTGAGCACACTGGGACTCTACATGGAAGCAAACAGACAGGTGCC
 CAAATTGCTGTCTGCTCTCTGTCAGAGAAGCAACACAGTCCATTTCGTAGCCAGATTGTGAGCAAAACAC
 CTGTTGGAGTTGACAGCCTCATTGGGCCAGAGACACAGATTGGAGAGAAGTCATCCATTAAGCGCTCAG
 TCATTGGCTCATCCTGTCTCATAAAGATAGAGTGACTATTACCAATTGCCTTCTCATGAACTCAGTCAC
 TGTGGAGGAAGGAAGCAATATCCAAGGCAGTGTCTATCTGCAACAATGCTGTGATCGAGAAGGGTGCAGAC
 ATCAAGGACTGCTTGATTGGAAGTGGCCAGAGGATTGAAGCCAAAGCTAAACGAGTGAATGAGGTGATCG
 TGGGAATGACCAGCTCATGGAGATC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>RG223570 representing NM_020365
 Red=Cloning site Green=Tags(s)

MEFQAVVMAVGGGSRMTDLTSSIPKLLPVGNKPLIWYPLNLLERVGFEEVIVVTTTRDVQKALCAEFKMK
 MKPDIVCIPDDADMGTADSLRYIYPKLTDLVLSCDLITDVALHEVVDLFRAYDASLAMLMRKQDSIE
 PVPQGKGKKKAVEQRDFIGVDSTGKRLLFMANEADLDEELVIKGSILQKHPRIRFHTGLVDAHL YCLKKY
 IVDFLMENGSI TSIRSELIPYLVRKQFSSASSQQGQEEKEEDLKKKELKSLDIYSFIKEANTLNLAPYDA
 CWNACRGDRWEDLSRSQVRCYVHIMKEGLCSRVTGLGLYMEANRQVPKLLSALCPEEPPVHSSAQIVSKH
 LVGVDSLIGPETQIGEKSSIKRSVIGSSCLIKDRVTITNCLLMNSVTVEEGSNIQGSVICNNAVIEKGAD
 IKDCLIGSGQRIEAKAKRYNEVIVGNDQLMEI

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_020365

ORF Size: 1356 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020365.5](#)

RefSeq Size: 1602 bp

RefSeq ORF: 1359 bp

Locus ID: 8891


UniProt ID: [Q9NR50](#)

Cytogenetics: 1p34.1

Domains: NTP_transferase

Gene Summary: The protein encoded by this gene is one of the subunits of initiation factor eIF2B, which catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP. It has also been found to function as a cofactor of hepatitis C virus internal ribosome entry site-mediated translation. Mutations in this gene have been associated with leukodystrophy with vanishing white matter. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]

Product images:



Circular map for RG223570