

## Product datasheet for **RG223657**

### HS2ST1 (NM\_012262) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	HS2ST1 (NM_012262) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HS2ST1
Synonyms:	dj604K5.2; NFSRA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG223657 representing NM_012262 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGGGCTCCTCAGGATTATGATGCCGCCAAGTTCAGCTGCTGGCGGTGGTGGCCTTCGCGGTGGCGA  
TGCTCTTCTTGGAAAACAGATCCAGAACTGGAGGAGTCCCGCTCGAAGCTAGAAAGGGCTATTGCAAG  
ACACGAAGTCCGAGAAATTGAGCAGCGACATAACAATGGATGGCCCTCGGCAAGATGCCACTTTAGATGAG  
GAAGAGGACATGGTGATCATTTATAACAGAGTTCCTCAAAACGGCAAGCACTTCATTTACCAATATCGCT  
ATGACCTGTGTCAAAGAATAAATACCATGTCCTTCATATCAACTACCAAAAATAATCCAGTGATGTC  
ATTGCAAGATCAGGTGCGCTTTGTAAGAATATAACTTCCTGGAAAGAGATGAAACCAGGATTTTATCAT  
GGACACGTTTTCTACTTGGATTTTGCAAAATTTGGTGTGAAGAAGAAACCAATTTACATTAATGTCAAA  
GGGATCCTATTGAGAGGCTAGTTTCTTATTACTTTCTGAGATTTGGAGATGATTATAGACCAGGGTT  
ACGGAGACGAAAACAAGGAGACAAAAGACCTTTGATGAATGTGTAGCAGAAGGTGGCTCAGACTGTGCT  
CCAGAGAAGCTCTGGCTTCAATCCCGTTCTTCTGTGGCCATAGCTCCGAATGCTGGAATGTGGGAAGCA  
GGTGGGCTATGGATCAAGCCAAGTATAACCTAATTAATGAATATTTCTGGTGGGAGTTACTGAAGAAT  
TGAAGATTTTATCATGTTATTGGAGGCAGCATTGCCCGGTTTTTCAGGGGTGCTACTGAACTCTATCGC  
ACAGGAAAGAAATCTCATCTTAGGAAAACACAGAGAAGAACTCCCACTAAACAAACCTTGGCAAAAC  
TACAGCAATCTGATATTTGGAAAATGGAGAATGAGTTCTATGAATTTGCACTAGAGCAGTTCCAATTCAT  
CAGAGCCCATGCCGTTCCGAGAAAAGATGGAGACCTCTACATCCTCGCACAAAACTTTTCTATGAAAAG  
ATTTACCTAAGTCGAAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG223657 representing NM\_012262  
 Red=Cloning site Green=Tags(s)

MGLLRIMMPPKLQLLAVVAFVAVMLFLENQIQKLEESRSKLERAIARHEVREIEQRHTMDGPRQDATLDE  
 EEDMVIIYNRVPKTASTSFTNIAIDLCAKNKYHVLHINTTKNNPVMSLQDQVRFVKNITSWKEMKPGFYH  
 GHVSYLDFAKFGVKKKPIYINVIKPIERLVSYFFLRFGDDYRPLRRRQGDKKTFDECVAEGGSDCA  
 PEKWLQIPFFCGHSSECVNVSRAWMDQAKYNLINEYFLVGVTEELDFIMLLEAALPRFFRGATELYR  
 TGKSHLRKTTTEKKLPTKQTIKLLQSDIWKMENEFYEFALQFQFIRAHAVREKDGDLIYLAQNFFYEK  
 IYPKSN

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_012262

**ORF Size:** 1068 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_012262.4](#)

**RefSeq Size:** 6708 bp

**RefSeq ORF:** 1071 bp

**Locus ID:** 9653

**UniProt ID:** [Q7LGA3](#)


**Cytogenetics:** 1p22.3

**Domains:** HS2ST

**Protein Pathways:** Heparan sulfate biosynthesis

**Gene Summary:** Heparan sulfate biosynthetic enzymes are key components in generating a myriad of distinct heparan sulfate fine structures that carry out multiple biologic activities. This gene encodes a member of the heparan sulfate biosynthetic enzyme family that transfers sulfate to the 2 position of the iduronic acid residue of heparan sulfate. The disruption of this gene resulted in no kidney formation in knockout embryonic mice, indicating that the absence of this enzyme may interfere with the signaling required for kidney formation. Two alternatively spliced transcript variants that encode different proteins have been found for this gene. [provided by RefSeq, Aug 2008]

### Product images:



Circular map for RG223657