

Product datasheet for **RG223667**

GLYAT (NM_201648) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	GLYAT
Synonyms:	ACGNAT; CAT; GAT
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG223667 representing NM_201648
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGTTACCATTGCAAGGTGCCAGATGCTGCAGATGCTGGAGAAATCCTTGAGGAAGAGCCTCCAG
CATCCTTAAAGGTTTATGGAAGTGTCTTTCACATAAACCATGGAATCCATTCAATCTGAAGGCTGTGGT
GGACAAGTGGCCTGATTTTAAACAGTGGTTGTCTGCCCTCAGGAGCAGGATATGACAGATGACCTTGAT
CACTATACCAATACTTACCAATCTACTCCAAAGATCCCCAACTGTGAGGAATTCCTTGGATCACCAG
AACTCATCAACTGGAACAGCATTACAGATTCAAAGTTCACAGCCTAGCCTGAATGAGGCTATACAAAA
TCTTGAGCCATTAAGTCTTCAAAGTCAAACAAACACACGATTCTCTATATGGCAGCTGAAACAGCC
AAGGAAGTACTCCTTCTGCTGAAATCAAAGATTTTATCTCCAGTGGTGGCAAACCAAGGCCATCA
ACCAAGAGATGTTTAACTCTCATCCATGGATGTTACCCATGCTCACTTGGTGAATAAATTCTGGCATT
TGGTGGTAATGAGAGGAGCCAGAGATTCATTGAGCGCTGCATTGACACCTTCCACCTGCTGTCTCCTG
GGGCCTGAGGGGACCCCTGTGTGCTGGGATCTAATGGACCAGACTGGAGAGATGAGAATGGCAGGCACCT
TGCCGGAATACCGGCTCCATGGCCTTGTGACGTATGTCATCTATTTCCACGCCAGAAATTGGGCAAACT
TGGGTTTCCTGTCTATTCTCATGTAGACTACAGCAATGAAGCTATGCAAAAAATGAGTTACACACTGCAA
CATGTTCCCATCCAGAAGCTGGAACCACTGGAAGTGTGTACCTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG223667 representing NM_201648
Red=Cloning site Green=Tags(s)

MMLPLQGAQMLQMLEKSLRKSLPASLKVYGTVFHINHGNPFLKAVVDKWPDFNTVVVCPQEQDMTDDL
HYTNTYQIYSKDPQNCQEFLLGSPELINWKHLQIQSSQPSLNEAIQNLAIAKSFVKVQTQRILYMAAETA
KELTPFLKSKILSPSGGPKAINQEMFKLSSMDVTHAHLVNKFHFGGNNERSQRFIERCIQTFTPCCLL
GPEGTPVCWDLMDQTGEMRMAGTLPEYRLHGLVTYVIYSHAQKLGLGFPVYSHVDYSNEAMQKMSYTLQ
HVPIPRSWNQWNCVPL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGTCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGCCGCGGAGCT

- - E E R V Stop

ACCN: NM_201648

ORF Size: 888 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_201648.1](#), [NP_964011.1](#)

RefSeq Size: 2068 bp

RefSeq ORF: 891 bp


Locus ID: 10249

UniProt ID: [Q6IB77](#)

Cytogenetics: 11q12.1

Gene Summary: The glycine-N-acyltransferase protein conjugates glycine with acyl-CoA substrates in the mitochondria. The protein is thought to be important in the detoxification of endogenous and xenobiotic acyl-CoA's. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG223667