

Product datasheet for **RG223735**

MAGEA11 (NM_001011544) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAGEA11 (NM_001011544) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAGEA11
Synonyms:	CT1.11; MAGE-11; MAGE11; MAGEA-11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG223735 representing NM_001011544 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAAGGTGAGCACTATGTTCTCAGAGGACGACTTCCAGTCAACAGAAAGAGCCCCATATGGTCCAC
AACTACAGTGGTCCCAGGATCTGCCAAGAGTCCAGGTTTTAGAGAACAGGCCAACCTGGAGGACAGGAG
TCCAGGAGAACCAGAGGATCACTGGAGGAGAACAAGTCTGTGGGGCCCCATACCCAGATATTTCCC
ACAGTTCGGCCTGCTGACCTAACAGAGTCATCATGCCTCTTGAGCAAAGAAGTCAGCACTGCAAGCCTG
AGGAAGGCCTTCAGGCCAAGAAGAAGACCTGGGCCTGGTGGGTGCACAGGCTCTCCAAGCTGAGGAGCA
GGAGGCTGCCTTCTCTCTACTCTGAATGTGGGCACTCTAGAGGAGTTGCCTGCTGCTGAGTACCA
AGTCTCTCCAGAGTCTCAGGAAGAGTCTTCTCTCCACTGCCATGGATGCCATCTTTGGGAGCCTAT
CTGATGAGGGCTCTGGCAGCCAAGAAAAGGAGGGGCCAAGTACCTCGCCTGACCTGATAGACCCTGAGTC
CTTTTCCAAGATATACTACATGACAAGATAATTGATTTGGTTCATTTATTGCTCCGCAAGTATCGAGTC
AAGGGCTGATCACAAGGCAGAAATGCTGGGGAGTGTCAAAAATTATGAGGACTACTTTCCTGAGA
TATTTAGGGAAGCCTCTGTATGCATGCAACTGCTCTTTGGCATTGATGTGAAGGAAGTGGACCCACTAG
CCACTCCTATGCTCTGTACCTCCCTCAACCTCTTTATGATGGCATAACAGTGTAAAGCAGAGCATG
CCCAAGTCTGGCCTCTGATAATAGTCTGGGTGAATCTTCATGGAGGGGAAGTGCATCCCTGAAGAGG
TTATGTGGGAAGTCTGAGCATTATGGGGTGTATGCTGGAAGGAGCACTTCTCTTTGGGGAGCCCAA
GAGGCTCCTTACCCAAAATTGGGTGCAGAAAAGTACCTGGTGTACCGGCAGGTGCCCGGCACTGATCCT
GCATGCTATGAGTTCCTGTGGGGTCCAAGGGCCACGCTGAGACCAGCAAGATGAAAGTTCTTGAGTACA
TAGCCAATGCCAATGGGAGGGATCCCACTTCTTACCCATCCCTGTATGAAGATGCTTTGAGAGAGGAGG
AGAGGGAGTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG223735 representing NM_001011544
 Red=Cloning site Green=Tags(s)

MSKVSTMFSEDDFQSTERAPYGPQLQWSQDLPRVQVFREQANLEDRSPRRTQRITGGEQVLWGPITQIFP
 TVRPADLTRVIMPLEQRSQHCKPEEGLQAQEEDLGLVGAQALQAEEQEAFFSSTLNVTLEELPAAESP
 SPPQSPQEESFSPTAMDAIFGSLSDDESGSQEKEGPSTSPDLIDPESFSQDILHDKIIDLVHLLLRKYRV
 KGLITKAEMLGSVIKNYEDYFPEIFREASVCMQLLFGIDVKEVDPTSHSYVLVTSNLNSYDGIQCNEQSM
 PKSGLLIIVLGVIFMEGNCIPEEVMWEVLSIMGVYAGREHFLFGEPKRLLTQNWVQEKYL VYRQVPGTDP
 ACYEFLWGPRAHAETSKMKVLEYIANANGRDPTSYPSTYEDALREEGEGV

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001011544

ORF Size:	1200 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001011544.2
RefSeq Size:	1893 bp
RefSeq ORF:	1203 bp
Locus ID:	4110
UniProt ID:	P43364
Cytogenetics:	Xq28
Gene Summary:	<p>This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p>