

Product datasheet for **RG223760**

C17orf49 (NM_174893) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: C17orf49
Synonyms: BAP18; HEPIS
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG223760 representing NM_174893
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGACGTCAGCGTCCACAAAGGTCGGAGAGATCTTCTCGGCGGCCGCGCCGCTTCACGAAGCTCGGGG
AGCTGACGATGCAGCTGCATCCCGTGGCCGACTCTTCTCCTGCGGGCGCGAAGTGGACGGAGACGGAAAT
AGAGATGCTGAGGGCTGCTGTGAAGCGATTTGGGACGATCTTAATCACATCAGCTGTGTCATCAAGGAA
CGGACAGTGGCCAGATAAAGGCCACTGTGAAACGCAAGGTATATGAAGATTCTGGCATCCCCCTTCCAG
CTGAGTCACCCAAGAAAGGGCCCAAGAAAGGTGGCATCTGGTGTCTTGTACCTCTCCAGCTGCCCTCC
TCCCAGCAGCTCCAGTGTCCCTGAGGCCGGGGTCCCCCATAAAGAAACAGAAGGCTGATGTGACACTC
AGTGCTCTGAACGACTCCGATGCCAACGTGACGTGGTGGATATTGAAGGGCTAGGAGAACTCCTCCAG
CTAAGAACTCAACTTCGACCAGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG223760 representing NM_174893
Red=Cloning site Green=Tags(s)


MTSASTKVGEIFSAAGAAFTKLGLT MQLHPVADSSPAGAKWTETETIEMLRRAAVKRFGDDLNHISCVIKE
RTVAQIKATVKRKVYEDSGIPLPAESPKKGPVKVASGLVSPPPAAPPSSSSVPEAGGPPIKKQKQADVTL
SALNDSANDVDVIEGLGETPPAKKLNFDQA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_174893

ORF Size: 516 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_174893.1](#), [NP_777553.1](#)

RefSeq Size:	841 bp
RefSeq ORF:	519 bp
Locus ID:	124944
UniProt ID:	Q8IXM2
Cytogenetics:	17p13.1
Gene Summary:	Component of chromatin complexes such as the MLL1/MLL and NURF complexes. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG223760