

## Product datasheet for **RG223825**

### CD36 (NM\_001001547) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CD36 (NM_001001547) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CD36
Synonyms:	BDPLT10; CHDS7; FAT; GP3B; GP4; GPIV; PASIV; SCARB3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG223825 representing NM\_001001547  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGGCTGTGACCGAACTGTGGGCTCATCGTGGGGCTGTCATTGGTGTCTCCTGGCTGTGTTTGGAG  
 GTATTCTAATGCCAGTTGGAGACCTGCTTATCCAGAAGACAATAAAAAGCAAGTTGTCCTCGAAGAAGG  
 TACAATTGCTTTTAAAAATTGGGTTAAAAACAGGCACAGAAGTTTACAGACAGTTTGGATCTTTGATGTG  
 CAAAATCCACAGGAAGTGATGATGAACAGCAGCAACATTCAAGTTAAGCAAAGAGTTTCTTATACGTACA  
 GAGTTCGTTTTCTAGCCAAGGAAAATGTAACCCAGGACGCTGAGGACAACACAGTCTCTTTCCTGCAGCC  
 CAATGGTGCCATCTCGAACCTTCACTATCAGTTGGAACAGAGGCTGACAACCTTACAGTTCTCAATCTG  
 GCTGTGGCAGCTGCATCCATATCTATCAAAATCAATTTGTTCAAATGATCCTCAATCACTTATTAACA  
 AGTCAAAATCTTCTATGTTCCAAGTCAGAACTTTGAGAGAAGTGTATGGGGCTATAGGGATCCATTTTT  
 GAGTTTGGTTCGGTACCCTGTTACTACCACAGTTGGTCTGTTTTATCCTTACAACAATACTGCAGATGGA  
 GTTTATAAAGTTTTCAATGGAAAAGATAACATAAGTAAAGTTGCCATAATCGACACATATAAAGGTA  
 GGAATCTGTCCTATTGGGAAAGTCACTGCGACATGATTAAATGGTACAGATGCAGCCTCATTTCCACCTTT  
 TGTGAGAAAAGCCAGGTATTGCAGTCTTTTCTTCTGATATTTGCAGGTCAATCTATGCTGATTTGAA  
 TCCGACGTTAATCTGAAAGGAATCCCTGTGTATAGATTTGTTCTTCCATCCAAGGCCCTTTCCTCTCCAG  
 TTGAAAACCCAGACAACATTGTTTCTGCACAGAAAAAATATCTCAAAAAATTGACATCATATGGTGT  
 GCTAGACATCAGCAAATGCAAAGAAGGGAGACCTGTGTACATTTCACTTCCTCATTTTCTGTATGCAAGT  
 CCTGATGTTTCAAGACCTATTGATGGATTAACCCAAATGAAGAAGAACATAGGACATACTGGATATTG  
 AACCTATAACTGGATTCACTTTACAATTTGCAAAACGGCTGCAGGTCAACCTATTGGTCAAGCCATCAGA  
 AAAAAATCAAGTATTAAGAATCTGAAGAGGAACATATTTGTCCTATTCTTTGGCTTAATGAGACTGGG  
 ACCATTGGTGATGAGAAGCAAACATGTTCAAGTCAAGTAACTGGAAAAATAAACCTCCTTGGCCTGA  
 TAGAAATGATCTTACTCAGTGTTGGTGTGGTGTGTTGTTGCTTTTATGATTTCAATTGTGCATGCAG  
 ATCGAAAAACAATAAAA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:**

>RG223825 representing NM\_001001547  
 Red=Cloning site Green=Tags(s)

MGCDRNCGLIAGAVIGAVLAVFGGILMPVGDLLIQKTIKKQVVLEEGTIAFKNWKVKGTEVYRQFWIFDV  
 QNPQEVMMNSSNIQVKQRGSYTYRVRFLAKENVTQDAEDNTVSFLQPNGAIFEPSLSVGTEADNFTVLNL  
 AVAAASHIYQNFVQMILNSLINKSKSSMFQVRTLRELLWGYRDPFLSLVPYPVTTTGLFYYPYNTADG  
 VYKVFNGKDNISKVAIIDTYKGRNLSYWESHCDMINTDAASFPPFVEKSQLQFFSSDICRSIYAVFE  
 SDVNLKGI PVYRFVLPKAFASPVENPDNYCFTEKIIISKNCTSYGVLDISKCKEGRPVYISLPHFLYAS  
 PDVSEPIDGLNPNEEEHRTYLDIEPITGFTLQFAKRLQVNLLVKPSEKIQLKLNKRNYIVPILWLN  
 ETGTIGDEKANMFRSQVTGKINLLGLIEMILLSVGVVMFVAFMISYCACRSKTIK

**TRTRPLE** - GFP Tag - V


**Restriction Sites:**

Sgfl-Mlul

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001001547

ORF Size: 1416 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001001547.1</a> , <a href="#">NP_001001547.1</a>
<b>RefSeq Size:</b>	2050 bp
<b>RefSeq ORF:</b>	1419 bp
<b>Locus ID:</b>	948
<b>UniProt ID:</b>	<a href="#">P16671</a>
<b>Cytogenetics:</b>	7q21.11
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Adipocytokine signaling pathway, ECM-receptor interaction, Hematopoietic cell lineage, PPAR signaling pathway
<b>Gene Summary:</b>	<p>The protein encoded by this gene is the fourth major glycoprotein of the platelet surface and serves as a receptor for thrombospondin in platelets and various cell lines. Since thrombospondins are widely distributed proteins involved in a variety of adhesive processes, this protein may have important functions as a cell adhesion molecule. It binds to collagen, thrombospondin, anionic phospholipids and oxidized LDL. It directly mediates cytoadherence of Plasmodium falciparum parasitized erythrocytes and it binds long chain fatty acids and may function in the transport and/or as a regulator of fatty acid transport. Mutations in this gene cause platelet glycoprotein deficiency. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2014]</p>