

Product datasheet for **RG223885**

RHOBTB1 (NM_001032380) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	RHOBTB1
Synonyms:	1700008H16Rik; 3110048G13Rik; AI173445; AI573858; AV350930; KIAA0740, MGC33059, MGC33841; MGC107538; Rho-related BTB domain containing 1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG223885 representing NM_001032380
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACGCTGACATGGACTACGAAAGACCAACGTTGAACTATCAAATGTGTGGTCGTGGGTGACAATG
 CCGTGGGGAAGACGCGCTTGATCTGTGCCAGGGCGTGCAACACCACACTCACGAGTATCAGCTGCTGGC
 CACCCACGTGCCAACAGTGTGGGCGATTGACCACTACCGCGTGTGCCAGGAGGTCTTGGAGCGTTCTCGG
 GATGTTGTTGATGAAGTGAGTGTCTCTCAGGCTTTGGGATACCTTTGGTGATCATCACAAGACAGAC
 GCTTTGCATATGGCAGGTCTGATGTTGTGGTCTCTGTTTTTCGATTGCTAATCCCAATTCCTAAATCA
 TGTGAAAAGCATGTGGTATCCAGAAATCAAGCACTTTTGCCCTCGAACACCCGTTATCCTTGTGGGTGC
 CAGCTTGATCTCCGCTATGCCGACCTGGAAGCTGTTAATCGAGCCAGGCGCCGTTAGCAAGGCCATAA
 AGAGAGGGGATATTTGCCCCAGAAAAGGCCGAGAGGTAGCAAGGAAGTGGCTTACCATACTATGA
 AACAGCGTGTGGACAGTTTGGTATCAAGGATGTGTTGACAATGCAATCCGAGCAGCGCTGATTTC
 CGCAGGCACCTGCAATTCTGGAATCCACCTAAAGAAAGTCCAGAAACCTTTACTTCAGGCACCTTCC
 TACCTCCAAAAGCCCTCCACCGGTATCAAAATCCAGAGTGTCTTCCATGGGGACAAATGAAGCTGC
 CTGTTTACTGGACAATCTCTATGTGCCGATGTTCTGTTATCCTTCAGGACCAGGAACACATCTTGCA
 CATCGAATTTACCTCGCTACCTCTTCTTCCAAATTTATGATCTGTTTTAATGGAATGTGAAGAATCCC
 CAAATGGGAGTGAAGGAGCCTGTGAGAAAAGAGAAGCAGAGCAGAGATTTCCAGGGGCGGATATTGAGTGT
 CGACCCAGAGGAAGAAAAGGGAGGAGGGCCCGCTAGGATTCCTCAGGCCGACCAAGTGAAGTCTTCAAAC
 AAGAGCCTGGTGGAGGCTCTGGGGCTGGAAGCCGAGGGTGCAGTTCCTGAGACACAGACTTTGACCGGAT
 GGAGTAAGGGGTTCAATGGCATGCACAGGAAATGCAAGTCAACCCCATTTCAAAGCGGATGGGGCCAT
 GACTGTGGTCAGGATGGACGCTTCAGTCCAGCCAGGCCCTTTTCGACCCTGCTCCAGTTTCTTTATACG
 GGACAATGGATGAAAAGGAAAAGGATTTGGTGGGCCTGGCTCAGATCGCAGAGGTCTCGAGATGTTTCG
 ATTTGAGGATGATGGTGAAAACATCATGAACAAGGAAGCCTTCATGAACCAGGAGATTACGAAAGCCTT
 TCACGTAAGGAAAGCCAAATCGGATAAAGAGTGTCTCAGCAAGGGAACGTTCTCGGACGTGACATTTAA
 TTGGACGATGGAGCCATCAGTGCCCAAGCCGCTGCTGATCTGTAGCTGTGAGTGGATGGCAGCCATGT
 TCGGGGGGTCATTTGTGAAAGTGCCAACAGTGAAGTGTATCTCCGAACATAAACAAGATATCAATGCA
 AGCAGTATTGGATTATCTATACCAAGCAGTTGTCTCCTAAGTGGATCTGGACCCGCTGGAATTAATT
 GCCTTGGCAAACAGATTTGCCTGCCACACTTGGTTGCACTTGCAAGACAGCATGCCGTTCCAGGAGTTGA
 CCAAAGCCGCCACGAGTGGCGTGGGCATTGACGAGAAAGTGTCTCTTACTTGGAAATGGCTCAGTTTCA
 CAATGCCACCAAGTTGGCCGCTGGTGTGTCACCACATCTGCACCAACTACAACAGTGTATGCTCCAAG
 TTCCGTAAGGAAATCAAATCAAATCTGCAGACAACAGGAATACTTCGAGCGGCACCGTGGCCCCCTG
 TGTGGTACCTGAAGGAAGAAGTCACTACCAGCGTGTGAAAAGGGAACGAGAGAAGGAAGATATTGCACT
 AAATAAGCATCGCTCAAGACGAAAGTGGTGCTTCTGGAATTCATCTCCAGCAGTGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG223885 representing NM_001032380
 Red=Cloning site Green=Tags(s)

MDADM DYERPNVETIKCVVVDNAV GKTRLICARACNTTLTQYQLLATHVPTVVAIDQYRVQEVLEERSR
 DVVDEVS VSLRLWDTFGDHHKDRRFAYGRSDVVVLCFSIANPNSLNHVKSMWYPEIKHFCPRTPVILVGC
 QLDLRYADLEAVNRARRPLARPIKRGDILPPEKGREVAKEGLPPYETSVFDQFGIKDVF DNAI RAALIS
 RRHLQFWKSHLKKVQKPLLQAPFLPPKAPPPVIKIPECPSMGTNEAACLDNPLCADVLFILQDQEHIFA
 HRIYLATSSSKFYDLFLMECEESPNGSEGACEKEKQSRDFQGRILSVDPEEERE EGPPRIPQADQWKSSN
 KSLVEALGLEAEGAVPETQTLTGWSKGFIMHREM QVNPISKRMGPMTVVRMDASVQPGPFRTLLQFLYT
 GQLDEKEKDLVGLAQIAEVLEMFDLRMMVENIMNKEAFMNQEITKAFHVRKANRIKECLSKGTFSDVTFK
 LDDGAISAHKPLLICSCWMAAMFGGSFVESANSEVYLPNINKISMQAVLDYLYTKQLSPNLDLPLELI
 ALANRFCLPHLVALAEQHAVQELTKAATSGVGIDGEVLSYLELAQFHNAHLAAWCLHHICTNYSVCSK
 FRKEIKSKSADNQEFERHRWPPVWYLKEEDHYQVRKEREKEDIALNKHSRRKWCFFWNSSPAVA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001032380

ORF Size: 2088 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001032380.1](#), [NP_001027552.1](#)

RefSeq Size: 4381 bp


RefSeq ORF: 2090 bp

Locus ID: 9886

Cytogenetics: 10q21.2

Gene Summary: The protein encoded by this gene belongs to the Rho family of the small GTPase superfamily. It contains a GTPase domain, a proline-rich region, a tandem of 2 BTB (broad complex, tramtrack, and bric-a-brac) domains, and a conserved C-terminal region. The protein plays a role in small GTPase-mediated signal transduction and the organization of the actin filament system. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Dec 2008]

Product images:



Circular map for RG223885