

Product datasheet for **RG223895**

SPRR3 (NM_005416) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: SPRR3

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG223895 representing NM_005416
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGTTCTTACCAGCAGAAGCAGACCTTTACCCACCACCTCAGCTTCAACAGCAGCAGGTGAAACAAC
CCAGCCAGCCTCCACCTCAGGAAATATTTGTTCCACAACCAAGGAGCCATGCCACTCAAAGTTCCACA
ACCTGGAAACACAAAGATTCCAGAGCCAGGCTGTACCAAGGTCCTGAGCCAGGCTGTACCAAGGTCCT
GAGCCAGGTTGTACCAAGGTCCTGAGCCAGGATGTACCAAGGTCCTGAGCCAGGTTGTACCAAGGTC
CTGAGCCAGGCTACCAAGGTCCTGAGCCAGGAGCATCAAGGTCCTGACCAAGGCTTCATCAAGTT
TCCTGAGCCAGGTCATCAAGTTCTGAGCAAGGATACCAAAAGTTCTGTGCCAGGCTACACAAAG
GTACCAGAGCCATGTCTTCAACGGTCACTCCAGGCCAGCTCAGCAGAAGACCAAGCAGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAA

Protein Sequence: >RG223895 representing NM_005416
Red=Cloning site Green=Tags(s)

MSSYQQKQTFPPQLQQQVKQPSQPPQEIFVPTTKEPCHSKVPQPGNTKIPEPGCTKVPEPGCTKVP
EPGCTKVPEPGCTKVPEPGCTKVPEPGYTKVPEPGSIKVPDQGF IKFPEPGA IKVPEQGYTKVPVPGYTK
VPEPCPSTVTPGPAQQKTKQK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_005416

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_005416.1, NP_005407.1

RefSeq Size: 984 bp

RefSeq ORF: 510 bp


Locus ID: 6707

UniProt ID: Q9UBC9

Cytogenetics: 1q21.3

Gene Summary: Cross-linked envelope protein of keratinocytes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG223895