

Product datasheet for **RG223986**

FBXO25 (NM_183420) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	FBXO25
Synonyms:	FBX25
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG223986 representing NM_183420
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCATTTTGGGTCAGGACTGGAGATCTCCTGGATGGAGTTGGATTAAGACAGAAGATGGCTGGAAGA
GATGTGAATCTTGTAGTCAGAACTTGAAAGAGAGAATAACCGTTGTAACATCAGTCACAGCATTATCTT
AAATAGTGAAGATGGAGAAATATTCATAATGAAGAGCATGAATATGCATCGAAAAAAGGAAAAAGGAC
CATTTTAGAAATGACACAAATACTCAAAGTTTTATCGTGAAAAATGGATCTATGTCCATAAAGAAAGCA
CAAAGGAAAGGCATGGCTATTGCACCTTGGGAGAAGCCTTTAATCGGTTAGACTTCTCAAGTGCAATTCA
AGATATCCGAAGGTTCAATTATGTGGTCAAACCTGTTGCAGCTAATTGCAAAATCCCAGTTAACTTCATTG
AGTGGCGTGGCACAGAAGAATTACTTCAACATTTGGATAAAATCGTTCAAAGGTTCTTGATGACCACC
ACAATCCTCGCTTAATCAAAGATCTTCTGCAAGACCTAAGCTCTACCTCTGCATTCTTATTAGAGGAGT
AGGGAAGTCTGTATTAGTGGGAAACATCAATATTTGGATTTGCCGATTAGAACTATTCTCGCCTGGCAA
CAACAGCTACAGGATCTTCAGATGACTAAGCAAGTGAACAATGGCCTCACCTCAGTGACCTTCCTCTGC
ACATGCTGAACAACATCCTATACCGGTTCTCAGACGGATGGGACATCATCACCTTAGGCCAGGTGACCCC
CACGTTGTATATGCTTAGTGAAGACAGACAGCTGTGGAAGAAGCTTTGTCACTACCATTTTGTGAAAAAG
CAGTTTTGTAGACATTTGATCCTTTTCAGAAAAAGGTCATATTGAATGGAAGTTGATGTACTTTGCATTC
AGAAACATTACCCAGCGAAGGAGCAGTACGGAGACACACTGCATTCTGTGCGCACTGCAGCATTCTCTT
TTGGAAGGACTCAGGACACCCCTGCACGGCGGCCGACCCTGACAGCTGCTTACGCCTGTGTCTCCGCAG
CACTTCATCGACCTCTTCAAGTTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG223986 representing NM_183420
Red=Cloning site Green=Tags(s)

MPFLGQDWRSPGWSWIKTEDGWKRCESCSQKLERENNRCNISHSIILNSEDGEIFNNEEHEYASKKRKKD
HFRNDTNTQSFYREKWIYVHKESTKERHGYCTLGEAFNRLDFSSAIQDIRRFNYVVKLLQLIAKSQLTSL
SGVAQKNYFNILDKIVQKVLDDHNPRLIKDLLQDLSSSLCILIRGVGKSVLVGNINIWICRLETILAWQ
QQLQDLQMTKQVNGLTSLDLPLHMLNNILYRFDGWDIITLGQVTPTLYMLSEDRQLWKKLCQYHFAEK
QFCRHLILSEKGHIEWKLMYFALQKHYPKEQYGDTLHFCRHCSILFWKDSGHPCTAADPDSCFTPVSPQ
HFIDLFKF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_183420

ORF Size: 1074 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_183420.2</u>
RefSeq Size:	2502 bp
RefSeq ORF:	1077 bp
Locus ID:	26260
UniProt ID:	<u>Q8TCJ0</u>
Cytogenetics:	8p23.3
Gene Summary:	This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG223986