

Product datasheet for **RG224045**

PRAMEF15 (NM_001098376) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	PRAMEF15
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG224045 representing NM_001098376
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGAAGATGAGCATCCGGACTCCACCCAGACTCCTGGAGCTTGCAGGGCGGAGCCTGCTGAGGGACCAAG
CTTTGGCCATGTCCACCCTGGAGGAGCTGCCACAGAACTTTCCCCCACTGTTTCATGGAGGCCTTCAG
CAGGAGACGCTGTGAGGCCCTGAAGCTGATGGTGCAGGCCCTGGCCCTCCGCCGCCTCCCTCTGAGGCCT
CTGATAAAGATGCCTTGTCTGGAGGCCCTCCAAGCTGTGCTCGATGGGCTTGATGCACTGCTTACCCAAG
GGGTTTCGTCACAGAGGTGAAACTCCAAGTCTGGATTTACAGGATGTCTGTGAGAACTTCTGGATGGT
TTGGTCTGAAGCTATGGCCATGGGTGCTTCTCAATGCCAAGAGGAACAAAAACCAGTGCAGGACTGT
CCAAGGATGAGAGGACGGCAGCCCTTGACTGTGTTTCGTAGAAGCTTTGGCTCAAGAACAGGACTCTGGATG
AATACCTCACCTACCTCCTTCTATGGGTCAAGCAGAGGAAAGATTTACTACACCTGTGCTGTAAGAAGCT
GAAAAATTTGGGAATGCCCTTCCGCAATATCAGAAGCATCCTGAAAATGGTGAACCTAGACTGTATCCAG
GAGGTGGAAGTGAATTGCAAGTGGTACTGCCATCCTGACACAGTTTACCCCATACCTGGGCCACATGA
GGAATCTTCAGAAGCTCGTTCTCTCCACATGGATGTCTCTCGCTACGTTTCCCAGAGCAGAAGAAGGA
GATTGTTACCCAGTTCACCACTCAGTTCCTCAAGCTGCGCTGCCTCCAAAAGCTTTATATGAACTCTGTT
TCTTTCTCGAAGGCCACCTGGACCAGCTGCTCAGCTGTCTGAAGACCTCGTTAAAAGTCTCACAATAA
CTAACTGTGTGCTTTTGAATCAGACTTGAAGCATCTATCCAGTGCCCGAGTATCAGTCAACTAAAGAC
CCTGGACCTGAGTGGCATCAGACTGACCAATTATAGTCTTGTGCCCTCCAAATTCTCTAGAAAAAGTT
GCAGCCACCCTTGAGTACCTGGATTTAGATGACTGTGGCATCATAGACTCCCAAGTCAACGCCATCCTGC
CTGCCCTGAGCCGCTGCTTTGAGCTCAACACCTTCAGCTTCTGTGGAAATCCCATCTGCATGGCCACCCT
GGAGAACCTGCTGAGCCACACAATCATACTCAAAAACCTATGTGTGGAGCTGTATCCTGCCCCCGAGAG
AGTTATGGTCTGATGGTACTCTCTGCTGGAGCAGATTGCTCAAATTAGGGCTGAGCTGATGAACAGAG
TGAGGGACTTAAGGCACCCCAAGAGGATCTTGTCTGTACTGACTACTGCCCTGACTGTGGCAACAGGTC
ATTTTATGACCTGGAGGCAGATCAATACTGCTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG224045 representing NM_001098376
 Red=Cloning site Green=Tags(s)

MKMSIRTPPRLLELAGRSLLRDQALAMSTLEELPTLEFPPLFMEAFSRRRCEALKLMVQAWPFRRLPLRP
 LIKMPCLEAFQAVLDGLDALLTQGVPRRWKLQVLDLQDVCENFMVWSEAMAHGCFNNAKRNNKPVQDC
 PRMRGRQPLTVFVELWLNRTLDEYLYLLLWVKQRKDLLHLCKKILGMPFRNIRSILKMVNLDCIQ
 EVEVNCKWVLPILTQFTPYLGHMRLQKLVLSHMDVSRYSPEQKKEIVTQFTTQFLKRLCLQKLYMNSV
 SFLEGLDQLLSCLKTSKVLITNCVLLSGLKHLSCPSISQLKTLDSGIRLTNYSLVPLQILLEKV
 AATLEYLDLDDCGIIDSQVNAILPALSRFELNTFSFCGNPICMATLENLLSHTIILKNLCVELYPAPRE
 SYGADGTLCSRFQIRAEMLNRVRLRHPKRILFCTDYCPDCGNRSFYDLEADQYCC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001098376

ORF Size: 1434 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001098376.3](#)

RefSeq Size: 1847 bp


RefSeq ORF: 1437 bp

Locus ID: 653619

UniProt ID: [Q5VWM5](#)

Cytogenetics: 1p36.21

Product images:



Circular map for RG224045