

Product datasheet for **RG224091**

CD164L2 (NM_207397) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD164L2 (NM_207397) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CD164L2
Synonyms: EAPG6122; UNQ6122
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG224091 representing NM_207397
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGCTCCGGGACCCCGCGCCTTGCAGACTGCGCTCTGTGGCGGCTGTTGCTGCCTCCTCCTATGTG
CCCAGCTGGCTGTGGCTGGTAAAGGAGCTCGAGGCTTTGGGAGGGGAGCCCTGATCCGCCTGAATATCTG
GCCGGCGTCCAAGGGGCTGCAACAGCTGGAGGCTGTGAGCACTGCGTGGAGGGAGACAGAGCGCGC
AATCTCTCCAGCTGCATGTGGGAGCAGTCCGGCCAGAGGAGCCAGGACTGTGTGCCCAATCTGAGG
TGGTCAAGGAAGTTGCTCCATCTACAACCGCTCAGAGGCATGTCCAGCTGCTCACCACCACCCACCTA
TGAAACGAAGACAGTCACAACAGGGAGCCCCAGTCCCTGAGGCCACAGCCCTGGATTTGACGGGGCC
AGCTTTATCGGAGGTGTGTGCTGGTGTGAGCCTACAGCGGTTGGCTTTCTTTGTGCTGCACTTCTCTCA
AGGCCAAGGACAGCACCTACCAGACGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG224091 representing NM_207397
Red=Cloning site Green=Tags(s)

MEAPGPRALRTALCGGCCLLLLCAQLAVAGKARGFGRGALIRLNIWPAVQGACKQLEVCEHCVEGDRAR
NLSSCMWEQCRPEEPGHCAQSEVVKEGCSIYNRSEACPAHHHTYEPKTVTTGSPVPEAHSPPGFDGA
SFIGGVVLVLSLQAVAFFVLHFLKAKDSTYQTL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207397.2</u>
RefSeq Size:	1517 bp
RefSeq ORF:	522 bp
Locus ID:	388611
UniProt ID:	<u>Q6UWJ8</u>
Cytogenetics:	1p36.11
Protein Families:	Transmembrane