

Product datasheet for **RG224502**

hnRNP R (HNRNPR) (NM_005826) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	hnRNP R (HNRNPR) (NM_005826) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	HNRNPR
Synonyms:	hnRNP-R; HNRPR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG224502 representing NM_005826
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTAATCAGGTGAATGGTAATGCGGTACAGTTAAAAGAAGAGGAAGAACCAATGGATACTTCCAGTG
 TAACTCACACAGAACTACAAGACACTGATAGAGGCAGGCCCTCCACAGAAGGTGGCAGAAAGACTTGA
 TGAATATTTTCAGACAGGATTGGTAGCTTATGTCGATCTTGATGAAAGAGCAATTGATGCTCTCAGGGAA
 TTTAATGAAGAAGGAGCTCTGTCTGACTACAGCAGTTCAAGGAAAGTACTTATCACATGTTTCAGAACA
 AAAGTGCATTTTTATGTGGAGTTATGAAGACCTACAGGCAGAGAGAGAAACAGGGGAGCAAGGTGCAAGA
 GTCCACAAAGGGACCTGATGAAGCGAAGATCAAGGCCTTGCTTGAGAGAAGTGGTATACTCTGGATGTA
 ACCACAGGACAGAGGAAGTATGGTGGTCTCCACCAGACAGTGTGTACTCTGGCGTCAACCTGGAATTG
 GAACGGAGGTATTTGTAGGCAAAATACCAAGGGATTTATATGAGGATGAGTTGGTCCCTTTTTGAGAA
 GGCCGGACCCATTTGGGATCTACGTCTTATGATGGATCCACTGTCCGGTCAAGAATAGAGGGTATGCATTT
 ATCACCTTCTGTGAAAGGAAGCTGCACAGGAAGCCGTGAAACTGTGTGACAGCTATGAAATTCGCCTG
 GTAAACACCTTGGAGTGTGCATTTCTGTGGCAAACAACAGACTTTTTGTTGGATCCATTCGAAGAATAA
 GACTAAAGAAAACATTTTGAAGAATTCAGTAAAGTCACAGAGGGTTTGGTGGACGTTATTTCTCTATCAT
 CAACCCGATGACAAAAGAAGAAATCGGGGGTCTGCTTCTTGAATATGAGGATCACAAGTCAGCAGCAC
 AAGCCAGACGCCGGCTGATGAGTGGAAAAGTAAAAGTGTGGGAAATGTAGTTACAGTTGAATGGCTGA
 CCCTGTGGAAGAACCAGATCCAGAAGTCATGGCTAAGGTAAGTAAAAGTTTTGTTGTGAGAAACTGGCTACT
 ACGGTGACAGAAGAAATTTGAAAAGTCAATTTCTGAATTTGAAAAGTCAAGAGTAAAGAAGTTGA
 AAGATTATGCATTTGTTCAATTTGAAGACAGAGGAGCAGCTGTTAAGGCTATGGATGAAATGAATGGCAA
 AGAAATAGAAGGGGAAGAAATGAAATAGTCTTAGCCAAGCCACCAGACAAGAAAAGGAAAGAGCGCCAA
 GCTGCTAGACAGGCTCCAGAAGCACTGCGTATGAAGATTACTACCACCTCCTCCTCGCATGCCAC
 CTCCAATTAGAGGTGGGGTCTGGTGGGGGAGAGGTGGATATGGCTACCCTCCAGATTACTACGGCTA
 TGAAGATTACTATGATGATTACTATGGTTATGATTATCACGACTATCGTGGAGGCTATGAAGATCCCTAC
 TACGGCTATGATGATGGCTATGCAGTAAGAGGAAGAGGAGGAGGAAGGGGAGGGCAGGTGCTCCACCAC
 CACCAAGGGGAGGGGAGCACCCTCCAAGAGGTAGAGCTGGCTATTCACAGAGGGGGCACCTTTGGG
 ACCACCAAGAGGCTCTAGGGTGGCAGAGGGGTCTGCTCAACAGCAGAGAGCCGTGGTCCCCTGGA
 TCTCGGGCAATCGTGGGGCAATGTAGGAGGCAAGAGAAAGGCAGATGGGTACAACCAGCTGATTCCA
 AGCGTCTGCAGACCAACAACCAACAGAAGTGGGGTCCCAACCCATCGCTCAGCAGCCGCTTCAGCAAGG
 TGGTGACTATTCTGGTAAGTATGTTACAATAATGACAACCAGGAATTTATCAGGATACTTATGGGCAA
 CAGTGGAAAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG224502 representing NM_005826
 Red=Cloning site Green=Tags(s)


MANQVNGNAVQLKEEEPMDTSSVTHTEHYKTLIEAGLPQKVAERLDEIFQTGLVAYVDLDERAIDALRE
 FNEEGALSVLQQFKESDLSHVQNKSAFLCGVMKTYRQREKQGSKVQESTKGPDEAIKALLERTGYTLDV
 TTGQRKYGGPPDSVYSGVQPGIGTEVFVGIKIPRDLYEDELVPLFEKAGPIWDLRLMMDPLSGQNRGYAF
 ITFCGKEAAQEAVKLSDSYEIRPGKHLGVCISVANNRLFVGSIPKNKTKENILEEFKSVTEGLVDVILYH
 QPDDKKKNRGFCFLEYEDHKSAAQARRRLMSGKVKVWGNVVTVEWADPVEEPDPEVMAKVKVLFVRNLAT
 TVTTEEILEKSFSEFGKLERVKKLDYAFVHFEDRGAAVKAMDEMNGKEIEGEEIEIVLAKPPDKRKRERQ
 AARQASRSTAYEDYYYHPPRMPPIRGRGRGGGGYPPDYGYEDYYDDYYGYDHYDRGGYEDPY
 YGYDDGYAVRGRGGGRGAPPPRGRGAPPRGRAGYSQRGAPLPPRGRSRRGGPAQQQRGRGRSRG
 SRGNRGGNVGGKRAKADGYNQPSKRRQTNNQNWGSQPIAQQPLQQGGDYSYGNYYNNDNQEFYQDITYGQ
 QWK

TRTRPLE – GFP Tag – V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_005826

ORF Size: 1899 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005826.5
RefSeq Size:	2663 bp
RefSeq ORF:	1902 bp
Locus ID:	10236
UniProt ID:	O43390
Cytogenetics:	1p36.12
Domains:	RRM
Gene Summary:	This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10. [provided by RefSeq, Jul 2014]