

Product datasheet for **RG224559**

WFDC8 (NM_181510) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WFDC8 (NM_181510) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	WFDC8
Synonyms:	C20orf170; dj461P17.1; HEL-S-292; WAP8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG224559 representing NM_181510 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGACTGTCCGAAGTGAAGGAGGGCACTTTCTCTCCATAGCCCCACCTTCTCCTGGAGGAATGTAG
CTTTCTGCTGCTTCTCTCCCTTGCTTTGGAGTGGACTTCTGCAATGCTGACCAAGAAGATCAAACACAA
ACCAGGGTTATGTCCAAAGAGAGGCTCACCTGTACCACTGAAGTCCGGACTCATGTAACACAGATTTT
GACTGCAAGGAATACCAGAAGTGTCTTTTTGCTGTCAGAAGAAGTGCATGGATCCCTTTCAAGAAC
CCTGCACGCTACCTGTGAGGCATGAAACTGTAATCATGAGGCACAGCGCTGGCATTCTGACTTTAAAAA
TTACCGCTGCACACCCTTCAAATACAGGGGCTGCGAAGGGAATGCCAACAACTTCTTAAATGAGGATGCC
TGCAGAACGCGCTGCATGTTAATTGTTAAGGATGGACAATGCCCACTCTTCCCTTCACTGAACGTAAGG
AGTGTCCACCTTCATGTCACAGTGACATCGATTGTCCCCAGACAGACAAATGTTGTGAATCCAGGTGTGG
CTTTGTTGTGCCAGGGCCTGGACAGTCAAAAAAGTTTCTGCCACGCAAGCCCTTGCTATGTACCAAG
ATTGATAAACCAAGTGCCTGCAGGATGAGGAGTGCCATTGGTGGAAAAGTGTCTCACATTGTGGAC
TGAATGTATGGACCCAGACGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG224559 representing NM_181510
Red=Cloning site Green=Tags(s)

MWTVRTEGGHFP LHSPTFSWRNVAFL LLLSLALEWTSAMLTKKIKHKPGLCPKERLTCTTELPDSCNTDF
 DCKEYQKCCFFACQKKCMDFPQEPCTLPVRHGNCNHEAQRWHDFKNYRCTPFKYRGCEGNANFLNEDA
 CRTACMLIVKDGQCPLFPF TERKECPPSCHSDIDCPQTDKCCESRCGFVCARAWTVKKGFCPRKPLLCTK
 IDKPKCLQDEECLVEKCCSHCGLKCMDPRR

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_181510

ORF Size: 723 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_181510.1 , NP_852611.1
RefSeq Size:	1269 bp
RefSeq ORF:	726 bp
Locus ID:	90199
UniProt ID:	Q8IUAA0
Cytogenetics:	20q13.12
Protein Families:	Secreted Protein
Gene Summary:	This gene encodes a member of the WAP-type four-disulfide core (WFDC) domain family. The WFDC domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor. The encoded protein contains a Kunitz-inhibitor domain, in addition to three WFDC domains. Most WFDC genes are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene belongs to the telomeric cluster. Two alternatively spliced transcript variants have been found for this gene, and they encode the same protein. [provided by RefSeq, Jul 2008]