

Product datasheet for **RG224635**

CCNP (NM_024877) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: CCNP (NM_024877) Human Tagged ORF Clone
 Tag: TurboGFP
 Symbol: CCNP
 Synonyms: CNTD2
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >RG224635 representing NM_024877
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGGTGAGAGGCAGGGACCAGGGTCCGGCTCCCGGCTCGGGCCTATCGTTAGGCGCTGGGCCCA
 GGCCCTCTCCTTTCAGAGTCTCGCTGCCTCCCTCGACGCAGAGCCTCAAGCGCCGAGTCCCGACGG
 CTTCCCGCGGGCCCACTGTCTCCCAAGACGCCTGGCAGGCCCGGGGCTGGAGGAGGCGCTGAGC
 GCGCTGGGGCTGCAGGGAGAACGCGAGTACGCCGGGACATCTTCGCCAAGTCATGGTGTGCCCGTGC
 TGCCCTGAGAGCCCTGCCCGCTGTGACCCCGGAGATGCGCGCCCTGGTGGTAGACTGGCTGGTCCA
 GGTGCACGAGTACCTGGTCTGGCTGGTACACACTTTATCTGGCGGTTACCTGCTTGATTCTACCTG
 AGCGCTGGCCGCTGCGTCTACATCGCCTGCAGCTGTGGGCGTGGCTTGCTTGTGGCGTGCAAAA
 TGAAGAGTGCCTGCTTCCCGAGCCCGCTTCTCTGCCTCCTGAGCGCGGACTCCTTCTCACGGCGGA
 GCTGCTGCGCGCCGAGCGTGCATCCTGAGCCGCTGGATTTCCGGCTGCACCACCCCGGCCGCTGCTG
 TGCTCGGGTGTGGCCGCTGGCAGGGAGCAGCCCCAGGTGATGTTACTTGCCACCTACTTCTGG
 AGCTGTCTTGTGGAGCCGAGCGCGGGATGGGAGCCGGTCTGTCGCGGCTGCGGCTCTGAGCCT
 GGCGCACCGCTTGTCTGACGGGGCGGGCTCCAGGCTCCAGCCAGAAGTTACAGGTGTAGTCTTGGCGGA
 GGAAGTGTATGGGTACCCGAGCTTACGGGACTTACCTTCTGGTCATTTTTACGGTCTCGGAGAATGA
 GAGACAATTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG224635 representing NM_024877
 Red=Cloning site Green=Tags(s)

MLVVRGRDQGSRSRLGPIVRRWAPRPSPLQSLAASLDAEPSSAAVPDGFAPGPTVSPRRLARPPGLEEALS
 ALGLQGEREYAGDIFAEVMVCRVLPRLRALPRAVTPEMRALVVDWLQVHEYLGLAGDTLYLAVHLLDSYL
 SAGRVRLHRLQLLGVACLFVACKMEECVLPPEFAFLCLLSADSFSAELLRAERRILSRDLFRLHHPGPLL
 CLGLLAALAGSSPQVMLLATYFLEL SLLAEAAAGWEPGRRAAAALSLAHRLLDGAGSRLQPELYRCSLGG
 GSVWGHRSFRDLPSWSFLRSRRMRDNY


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_024877

ORF Size: 921 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024877.4](#)

RefSeq Size: 1817 bp

RefSeq ORF: 924 bp

Locus ID: 79935

UniProt ID: [Q9H8S5](#)

Cytogenetics: 19q13.2

Gene Summary: Seems to be involved in the regulation of proliferation and migration.[UniProtKB/Swiss-Prot Function]