

Product datasheet for **RG224774**

MACC1 (NM_182762) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MACC1 (NM_182762) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MACC1
Synonyms:	7A5; SH3BP4L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG224774 representing NM_182762
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCTAATCACTGAAAGAAAACATTTTCGGTCAAGAAATTGCACAAAGTATGTCTGAAGCAAATTTGA
 TTGACATGGAAGCTGGAAAAGTCTCAAAAAGTTGCAATATTACAGAATGCCAGGACCCAGACTTGCTTCA
 CAATTGGCCGGATGCTTTCCACCCTTCGTGGTAATAATGCTTCCAAAAGTTGCAAAATCCATTCTGGAATCAA
 CTGTCTGCTTCTAACCCATTTTTGGATGACATAACTCAACTAAGAAAATACAGGAAGAGAAAATAATATTT
 CCATCTTAAAGGAAGATCCTTTTCTTTCTGTAGAGAAAATAGAAAATGGAACTCTTTTGATTCTCCGG
 TGATGAACTTGATGTGCATCAGTTACTTAGGCAGACTTCCCTAAGAAATCTGGAAGATCTAAAAGTGT
 TCAGAACTCTGGACATTTAGACGACACAGCACATGCCATCAGAGTACATAACTCTGACCAGATCC
 TACTACACGACTTAGAGTGGCCTAAAAATGATCGGGAGGCTTATAAAATGGCTTGGTTAAGTCAACGCCA
 GCTGGCCCGCTCCTGCCTTGATTTGAATACAATTAGTCAGAGCCCTGGATGGGCCAGACACAATTGCG
 GAGGTCACCATAGCTTGCAAAGTAAACCATCAAGGAGGGTCAGTACAATTACCTGAATCAGACATCACTG
 TTCATGTGCCCAAGGCATGTGGCTGTGGGAGAATCCAAGAGGTGTCTCTAAGGGCTTTCCTTGATCC
 GCCACACATGCTTAACCATAAATCTTTCTGTGACTGTGAGCCCGTTGTTGGAAATCATGTTAGGCAACCTC
 AATACAATGGAAGCCCTTTTCTGGAGATGAAAATGGGGCTGAAGTAAAGAAAGGATCCTTTAGCCAAAG
 TCATGACAGAAATGGTGTGTTTACACAGCTTGGGTAAAGAAGGCCCTTTTAAAGTTTTAAGCAACTGCTA
 CATTTATAAAGACACCATCCAAGTCAAGCTAATCGACTGAGTCAGGTAATGTATCTAGTGGTTGCTGCA
 CAAGCTAAAGCTCTTCCGTCACCACTGCCACCATTTGGGATTATCCACAAAACCACTCAATTGGAA
 TTTATGGACCCAAATATATCCATCCCAATTTTACTGTTGTTTTAACAGTTTGTGGACACAATTATATGCC
 AGGACAGCTTACAATTTCTGATATTAAGAAGGGTGGAAAAACATATCTCCAGTTGTGTTTCAGCTCTGG
 GGAAGCAGTCATTTTTACTTGACAAGCCACAAGATTTAAGTATTTCTATTTTTCTGTGATCCTGATT
 TTGAAGTAAAGACAGAAGGAGAAAAGAAAATTAACAAAAGCAGTTGGAAGCAGGTGAAGTAGTTCA
 TCAACAATTTTTATTTCTTTAGTTGAGCACAGAGAGATGCACTTGTGTTGATTTTTGTGTTCAAGTGGAG
 CCTCCCAATGGTGAACCAGTTGCACAGTTCTCTATCACTACTCCTGATCCAACCCCAAACCTAAAAAGAC
 TCTTGAATCTGCCAGGCTATTTGCAGAAGAAGGAGGAAAATCAAGTCTGCTCCTTATCACAAAAATTTCT
 TGTTAAATATCCTACATTTCAAGATAAAACATTGAACTTTAGCAACTATGGGGTAACCTGAAGGCAGTG
 CTAAGACAAAAGCAAGATTGATTACTTCTTGAATATTTCAAAGGGGACACAATAGCTCTCCTCGGGGAAG
 GTAAGGTAAAAGCTATTGGTCAGTCCAAAGTGAAGAATGGTATGTAGGAGTCTCAGAGGTAAGATTGG
 ACTTGTACTGCAAAAAATGTCAAGGTGATTTCAAAGGAGCAAGTAAATGTTTATGTCAGATAGTGTCTTT
 ACAACCAGAAATCTTCTTGAACAGATTGTCTGCCTTTAAAAAATTTGACTTATATCTACTCAGTTGTAT
 TAACCTTGGTGTGAGAAAAGTTTATGATTGAAAAGTTTTAGCTGATGTCTGGGTTACTCACATCTGTC
 CCTGGAAGATTTTATCAAAATCAAGCAGACAAAGAATCAGAGAAAAGTTTCTTATGTTATAAAGAAGTTA
 AAGGAAGATTGCCACACAGAGAGAAAATACAAGGAAGTTTCTGTATGAACTATTGTGGCTCTTCTGAAA
 TGGATTGCCAAGAGTTAGTCGCACGTCTCATCAAGAAGCTGCTGTTCTGACTTCAGCTGCAAGCTTGG
 AAAAGGCTGGAGGAACTAGCTGAAAAGTTAGTACGACTCACAAAGCAACAAATGGAGGCATATGAAAT
 CCTCATCGAGGAAACACTGGAGATGTTGCTGTTGAGATGATGTGAAAACCTGCCTATGATTTTCTGTATA
 CCTGGAGTGCTCACTATGAAAATAACTACACAGATGTGTTACAAGACCTTCAGTCAGCTTTGGACAGAAT
 GAAAAACCTGTGACTAAACACTGGAGAGAATTAAGTGGAGTTTTAATACTAGTAAATCTTTGGAGGTT
 TTGAGAGTAACTGCATTCTCCACTTCTGAGGAAGTA

ACGGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG224774 representing NM_182762
Red=Cloning site Green=Tags(s)

```
MLITERKHFRSGRIAQSMSEANLIDMEAGKVS KSCNITECQDPDLLHNWPAFTLRGNNASKVANPFWNQ  
LSASNPFLDDITQLRNNRKRNNISILKEDPFLFCREIENGNSFDSSGDEL DVHQLLRQTSSRNSGRSKSV  
SELLDILDDTAHAHQSIHNSDQILLHDLEWPKNDREAYKMAWLSQRQLARSCLDLNTISQSPGWAQTQLA  
EVTIACKVNHQGGSVQLPESDITVHVPQGHVAVGEFQEVSLRAFLDPPHMLNHNLSCTVSPILLEIMLGNL  
NTMEALLLEMKIGAEVRKDPFSQVMTMVCLHSLGKEGPFKVL SNCYIYKDTIQVKLIDL SQVMYLVVAA  
QAKALPSPAATIWDYIHKTTSIGIYGPKYIHPNFTVVL TVCGHNYMPGQLTISDIKKGKNI SPVVFQLW  
GKQSFLLDKQPQLSISIFSCDPDFEVKTEGERKEIKQKQLEAGEVVHQQFLFSLVEHREMHLDFCVQVE  
PPNGEPVAQFSITTPDPTPNLKRLLNLPGLYQKKEEIKSAPLSPKILVKYPTFQDKTLNFSNYGVTLKAV  
LRQSKIDYFLEYFKGDTIALLGEGKVKAIGQSKVKEWYVGVLRGKIGLVHCKNKVVISKEQVMFMSDSVF  
TTRNLLEQIVLPLKCLTYIYSVVLTLVSEKVDWKVLADVLGYSHLSLEDFDQIQADKESEKVSYVIKLL  
KEDCHTERNTRKFLYELIVALLKMDCQELVARLIQEAAVLTSVAVKLGKGWRELAEKLVR LTKQQMEAYEI  
PHRGNTGDVAVEMMKPAYDFLYTWSAHYGNNYTDVLDLQDLSALDRMKNPVTKHWREL TGVLLVNSLEV  
LRVTAFSTSEEV
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_182762


ORF Size: 2556 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_182762.2</u> , <u>NP_877439.2</u>
RefSeq Size:	3188 bp
RefSeq ORF:	2559 bp
Locus ID:	346389
UniProt ID:	<u>Q6ZN28</u>
Cytogenetics:	7p21.1
Gene Summary:	MACC1 is a key regulator of the hepatocyte growth factor (HGF; MIM 142409)-HGF receptor (HGFR, or MET; MIM 164860) pathway, which is involved in cellular growth, epithelial-mesenchymal transition, angiogenesis, cell motility, invasiveness, and metastasis. Expression of MACC1 in colon cancer (MIM 114500) specimens is an independent prognostic indicator for metastasis formation and metastasis-free survival (Stein et al., 2009 [PubMed 19098908]). [supplied by OMIM, Mar 2009]

Product images:



Circular map for RG224774