

Product datasheet for **RG224892**

NSRP1 (NM_032141) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NSRP1 (NM_032141) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	NSRP1
Synonyms:	CCDC55; HSPC095; NSrp70
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG224892 representing NM_032141
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGATTCCGGGCAGGCAGTATGGGCTATTTTGC~~AAAGAAAACACAGCAGTTGCACCTGTTTTC~~
~~AAAAACCATCAGTGTGGGAATGATTCTGATGATGATGATGAGACCTCTGTGAGTGAAAGCCTTCAGAG~~
~~GGAAGCTGCTAAGAAGCAGGCCATGAAACAGACCAAACCTGGAATCCAGAAGGCCCTGCAGAAGATGCT~~
~~ACTGTGTATGAATATGACAGTATTTATGATGAAATGCAGAAAAAAGGAGGAAAATAATCCCAAATTCG~~
~~TTTTGGGAAAGACAGAAAGCCCAAGTATATTCACAACCTGCTAAAAGCAGTTGAGATCAGAAAAAGGA~~
~~ACAGGAAAAAGAATGAAAAAGAAAATACAGAGAGAACGAGAAAATGAAAAAGGGGAGTTTGATGATAAA~~
~~GAAGCATTGTGACATCTGCATAAAGAAAAACTGCAAGAGAGAGCTGAAGAAGAAGAAAGAGAAAAAGA~~
~~GGGCTGCTGCACTGGAAGCATGTTGGATGTAACCAAGCAGAAAGATCTCAGTGGATTTTATAGGCACCT~~
~~ATTAATCAAGCAGTTGGTGAAGAGGAAGTACCTAAATGCAGCTTCGTGAAGCCAGATCTGGTATAAAG~~
~~GAAGAAAAATCAAGGGGCTTCTCCAATGAAGTAAGTTCAAAAAACAGAATACCACAAGAGAAATGCATTC~~
~~TTCAAACTGATGTGAAAGTAGAGGAAAACCCAGATGCAGACAGTACTTCGATGCTAAGAGCAGTGCCGA~~
~~TGATGAAATAGAAGAACTAGAGTGAATGCAGAAAGGAAAAGGTCATAGAGACCCCTGAGAATGACTTC~~
~~AAGCACCACAGGAGTCAAACCACTCTCGGTCACCTAGTGAAGAAAGAGGGCACAGTACCAGGCACCACA~~
~~CGAAAGGATCACGAACGTCGAGAGGACATGAGAAAAGGGAAGATCAGCACCAGCAGAAGCAATCCAGAGA~~
~~CCAAGAGAACCATTACTGACCGTGATTACCGGAAAGAAAGGGATTCTCATAGGCACAGAGAGGCCAGT~~
~~CATAGAGATTCCCATTGGAAGAGGCATGAACAGGAAGATAAACAAGGGCGAGGGACCAAGAGAAAAGAA~~
~~GTGACAGAGTATGAAAAAGGGAAGAAAGATAGGGAGAAAATATTCCAAAGAGAACAAGAAAGATAGACA~~
~~ACAAAATGATCAGAACCACCCAGTGAGAAAGGAGAGAAAGGAAGAGAAAAGCAAAGCAAAGGAAGGCAT~~
~~ATGAAAGTAAGGAAGGAAAGATATGAAAATAATGATAAATACAGAGATAGAGAAAAACGAGAGGTAGGTG~~
~~TTCAGTCTTCAGAAAGAAATCAAGACAGAAAGGAAAGCAGCCAAATTCAGGGCAAAGGATAAATTTCT~~
~~TGACCAAGAAAGATCCAACAAAATGAGAAACATGGCAAAGGACAAAGAAAGAAACCAAGAGAAACCTCT~~
~~AATTCTGAATCATCACTGGGAGCAAAACACAGACTCACAGAGGAAGGGCAAGAGAAGGGTAAAGAACAAG~~
~~AGAGACCACCTGAGGCAGTGAGCAAGTTTGCAAAGCGGAACAATGAAGAACTGTATGTCAGCTAGAGA~~
~~CAGGTAATGGCCAGGCAGATGGCGCGGTTAATGCAAAGACCTATATTGAGAAAGAAGATGAT~~

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG224892 representing NM_032141
 Red=Cloning site Green=Tags(s)

MAIPGRQYGLILPKKTQQLHPVLQKPSVFGNDSDDDDETS~~SVSESLQREAAKKQAMKQTKLEIQKALAEDA~~
~~TVYEYDSIYDEMQKKKEENPKLLL GKDRPKYIHNLLKAVEIRKKEQEKRMKKIQREREMEKGEFDDK~~
~~EAFVTSAYKKKLQERAEEREKRAAALEACL~~DVTKQKDL~~SGFYRHLLNQAVGEEEVPKCSFREARSGIK~~
~~EEKSRGFSNEVSSKNRIPQEKCILQTDVKEENPDADSD~~FD~~AKSSADDEIEETRVNCRREKVIETPENDF~~
~~KHHRQNHSRSPSEERGHSTRHHTKGSRTSRGHEKREDQH~~QKQSRDQENHYTDRDRYKERD~~SHRHREAS~~
~~HRDSHWKRHEQEDKPRARDQRERSDRVVKREKREKYS~~QREQERDRQNDQNR~~SEKGEKEEKSKAKEEH~~
~~MKVRKERYENNDKYRDREKREVGQSSERNQDRKESS~~PN~~SRAKDKFLDQERSNKMNRMAKDKERNQEKPS~~
~~NSESSLGAKHRLTEEGQEKGKEQERPP~~EA~~VSKFAKRNN~~EE~~TVMSARDRYLARQMARVNAKTYIEKEDD~~

TRTRPLE - GFP Tag - V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN:	NM_032141
ORF Size:	1674 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032141.4
RefSeq Size:	2547 bp
RefSeq ORF:	1677 bp
Locus ID:	84081
UniProt ID:	Q9H0G5
Cytogenetics:	17q11.2
Gene Summary:	RNA-binding protein that mediates pre-mRNA alternative splicing regulation. [UniProtKB/Swiss-Prot Function]