

## Product datasheet for **RG225041**

### PLAC8 (NM\_001130715) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PLAC8 (NM\_001130715) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** PLAC8  
**Synonyms:** C15; DGIC; onzin; PNAS-144  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG225041 representing NM\_001130715  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAAGCTCAGGCGCCGGTGGTCGTTGTGACCCAACCTGGAGTCGGTCCCGGTCCGGCCCCCAGA  
ACTCCAAGTGGCAGACAGGCATGTGTGACTGTTTCAGCGACTGCGGAGTCTGTCTCTGTGGC  
ACATTTTGTTCCTCGCCTTGGGTGTCAAGTTGCAGCTGATATGAATGAATGCTGTCTGTGTG  
GAACAAGCGTCGCAATGAGGACTCTCTACAGACCCGATATGGCATCCCTGGATCTATTTGTG  
ATGACTATATGGCAACTCTTTGCTGTCTCATTGTA

CTCTTTGCCAAATCAAGAGAGATATCAACAGAAGGAGGCCATGCGTACTTTC  
**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA**

**Protein Sequence:** >RG225041 representing NM\_001130715  
Red=Cloning site Green=Tags(s)  
MQAQAPVVVVVTPQGVGPGPAPQNSNWQTMCDCFSDCGVCLCGTFCFPCLGCQVAADMNECCLCGTSVAM  
RTLYRTRYGIPGSICDDYMATLCCPHCTLCQIKRDINRRRAMRTF

**TRTRPLE - GFP Tag - V**

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2369\\_f05.zip](https://cdn.origene.com/chromatograms/ja2369_f05.zip)

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001130715

**ORF Size:** 345 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001130715.2](#)

**RefSeq Size:** 1003 bp


**RefSeq ORF:** 348 bp

**Locus ID:** 51316

UniProt ID: [Q9NZF1](#)

Cytogenetics: 4q21.22

**Product images:**



Circular map for RG225041