

Product datasheet for **RG225079**

CD59 (NM_001127227) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CD59 (NM_001127227) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CD59
Synonyms: 1F5; 16.3A5; EJ16; EJ30; EL32; G344; HRF-20; HRF20; MAC-IP; MACIF; MEM43; MIC11; MIN1; MIN2; MIN3; MIRL; MSK21; p18-20
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG225079 representing NM_001127227
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGAATCCAAGGAGGGTCTGTCTGTTCTGGGCTGCTGCTCGTCTGGTGTCTTCTGCCATTCAGGTC
 ATAGCCTGCAGTGTACAACCTGCTAACCCTGCTGACTGCAAAACAGCCGTCAATTGTTTCATCTGA
 TTTTGATGCGTGTCTCATTACCAAAGCTGGGTTACAAGTGTATAACAAGTGTGGAAGTTTGAGCATTGC
 AATTTCAACGACGTCACAACCCGCTTGAGGAAAATGAGCTAACGTAAGTACTGCTGCAAGAAGGACCTGT
 GTAACCTTAACGAACAGCTTAAAATGGTGGGACATCCTATCAGAGAAAACAGTCTTCTGCTGGTGAC
 TCCATTTCTGGCAGCAGCCTGGAGCCTTCATCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG225079 representing NM_001127227
 Red=Cloning site Green=Tags(s)
 MGIQGGSVLFGLLLVLAVFCHSGHSLQCYNCPNPTADCKTAVNCSSDFDACLITKAGLQVYNKCKWKFHC
 NFNDVTTTRLRENELTYCCKKDLNCFNEQLENGGTSLSKTVLLLLVTPFLAAAWSLHP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001127227.1, NP_001120699.1</u>
RefSeq Size:	7573 bp
RefSeq ORF:	387 bp
Locus ID:	966
UniProt ID:	<u>P13987</u>
Cytogenetics:	11p13
Protein Families:	Druggable Genome
Protein Pathways:	Complement and coagulation cascades, Hematopoietic cell lineage
Gene Summary:	This gene encodes a cell surface glycoprotein that regulates complement-mediated cell lysis, and it is involved in lymphocyte signal transduction. This protein is a potent inhibitor of the complement membrane attack complex, whereby it binds complement C8 and/or C9 during the assembly of this complex, thereby inhibiting the incorporation of multiple copies of C9 into the complex, which is necessary for osmolytic pore formation. This protein also plays a role in signal transduction pathways in the activation of T cells. Mutations in this gene cause CD59 deficiency, a disease resulting in hemolytic anemia and thrombosis, and which causes cerebral infarction. Multiple alternatively spliced transcript variants, which encode the same protein, have been identified for this gene. [provided by RefSeq, Jul 2008]