

Product datasheet for **RG225174**

CNBP (NM_001127194) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CNBP (NM_001127194) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CNBP
Synonyms: CNBP1; DM2; PROMM; RNF163; ZCCHC22; ZNF9
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG225174 representing NM_001127194
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCAGCAATGAGTGTCTCAAGTGTGGACGATCTGGCCACTGGGCCCGGGAATGTCCTACTGGTGGAG
 GCCGTGGTCGTGGAATGAGAAGCCGTGGCAGAGGTTCCAGTTTGTTCCTCGTCTTCCAGACATTTG
 TTATCGCTGTGGTGTGAGTCTGGTCATCTTGCCAAGGATTGTGATCTTCAGGAGGATGTTGAAGCCTGCTAT
 AACTGCGGTAGAGGTGGCCACATTGCCAAGGACTGCAAGGAGCCCAAGAGAGCGAGAGCAATGCTGCT
 ACAACTGTGGCAAACCAGGCCATCTGGCTCGTGACTGCGACCATGCAGATGAGCAGAAATGCTATTCTTG
 TGGAGAATTCGGACACATTCAAAAGACTGCACCAAAGTGAAGTCTATAGGTGTGGTAAAAGTGGTCAT
 GTAGCCATCAACTGCAGCAAGACAAGTGAAGTCAACTGTTACCGCTGTGGCGAGTCAGGGCACCTTGAC
 GGAATGCACAATTGAGGCTACAGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG225174 representing NM_001127194
 Red=Cloning site Green=Tags(s)

MSSNECFKCGRSGHWARECPTGGGRGRGMRSRGRGFQVSSSLPDICYRCGESGHLAKDCDLQEDVEACY
 NCGRGGHIAKDCKEPREREQCCYNCGKPGHLARDCDHADEQKCYSCGEFGHIQKDCTKVKCYRCGETGH
 VAINCSKTSEVNCYRCGESGHLARECTIEATA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001127194.2</u>
RefSeq Size:	3378 bp
RefSeq ORF:	519 bp
Locus ID:	7555
UniProt ID:	<u>P62633</u>
Cytogenetics:	3q21.3
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	This gene encodes a nucleic-acid binding protein with seven zinc-finger domains. The protein has a preference for binding single stranded DNA and RNA. The protein functions in cap-independent translation of ornithine decarboxylase mRNA, and may also function in sterol-mediated transcriptional regulation. A CCTG expansion from <30 repeats to 75-11000 repeats in the first intron of this gene results in myotonic dystrophy type 2. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2016]