

## Product datasheet for **RG225178**

### **TMEM240 (NM\_001114748) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM240 (NM\_001114748) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** TMEM240  
**Synonyms:** C1orf70; SCA21  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG225178 representing NM\_001114748  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTCCATGAGCGCAACACCATGATCTTCATGATTCTGGGGCGTCGGTCGTGATGGCCATCGCGTGCT  
 TGATGGACATGAACGCGCTGCTGGACCGATTCCACAACACATCTCCCGCACCTGCGGGCGAGGACCG  
 CGTCTGCCACTGCAACTGTGGCCGGCACCATATCCACTACGTGATCCCGTACGACGGGGACCAGTCGGTG  
 GTGGACGCTCCGAGAATACTTTGTGACGACAGTGTGACCAAGCAGGAGATCGACCTCATGCTGGGGC  
 TGCTGCTGGGCTTTTGCATCAGCTGGTTCCTGGTGTGGATGGACGGCGTCCTGCACTGCGCCGTGCGCGC  
 CTGGAGAGCCGGACGGCGCTACGATGGCTCGTGGACCTGGCTGCCAAGCTGTGCAGCCTGCGGGAGCTG  
 GGCCGGCGCCGACAGGCCCTTCGAGGAGGCCCGCGGAACATGGTACACGTGAAGCAGAACTCTACC  
 ACAATGGCCACCCAGCCCGCGGCACCTG

**ACGCGT**ACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG225178 representing NM\_001114748  
 Red=Cloning site Green=Tags(s)

MSMSANTMIFMILGASVVMIAIACLMMDNALLDRFHNYILPHLRGEDRVCHCNCGRHHIHYVIPPYDGDQSV  
 VDASENYFVTDSVTKQEIDLMLGLLLGFCISWFLVWMDGVLHCAVRAWRAGRRYDGSWTWLPKLCSLREL  
 GRRPHRPFEEAAGNMVHVQKLYHNGHPSRHL

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001114748

**ORF Size:** 519 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001114748.2](#)

**RefSeq Size:** 1128 bp


**RefSeq ORF:** 522 bp

**Locus ID:** 339453

**UniProt ID:** [Q5SV17](#)

Cytogenetics:	1p36.33
MW:	19.9 kDa
Gene Summary:	This gene encodes a transmembrane-domain containing protein found in the brain and cerebellum. Mutations in this gene result in spinocerebellar ataxia 21. [provided by RefSeq, Dec 2014]

### Product images:



Circular map for RG225178