

Product datasheet for **RG225204**

PHR1 (PLEKHB1) (NM_001130035) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PHR1 (PLEKHB1) (NM_001130035) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PHR1
Synonyms: KPL1; PHR1; PHRET1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG225204 representing NM_001130035
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCCTGGTGAGGGGCGGCTGGCTGTGGAGACAGAGCTCCATCCTCCGCCGCTGGAAGCGGAAGTGGT
 TTGCCCTGTGGCTGGACGGGACCCTGGGATACTACCACGATGAGACAGCGCAGGACGAGGAGGACCGTGT
 GCTCATCCACTTCAATGTCGCTGACATAAAGATCGGCCAGAGTGCCATGATGTGCAGCCCCAGAGGGC
 CGGAGCCGAGATGGCCTGCTGACTGTGAACCTACGGGAAGGCGGCCGCTGCACCTCTGTGCGGAGACCA
 AGGATGATGCCCTAGCATGGAAGACAGCACTGCTGGAGGCAAATCCACCCGGTGCGCTCTACAGCCC
 GTACCAAGACTACTACGAGGTGGTGCCCCCAATGCACACGAGGCCACGTATGTCCGCAGCTACTACGGA
 CCGCCCTACGCAGGCCCTGGCGTGACGCACGTGATAGTGCGGGAGGATCCCTGCTACAGCGCCGGCGCCC
 CTCTGGCCATGGGCATGCTTGCGGGAGCCGCACTGGGGCGGCGCTGGGCTCGCTCATGTGGTCGCCCTG
 CTGGTTC

A**CGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG225204 representing NM_001130035
 Red=Cloning site Green=Tags(s)


MALVRGGWLWRQSSILRRWKRNFALWLDGTLGYHDETAQDEEDRVL IHFNVVDIKIGPECHDVQPPEG
 RSRDGLLTVNLREGGRLHLCAETKDDALAWKTALLEANSTPVRVYSPYQDYEVVPPNAHEATYVRSYYG
 PPYAGPGVTHVIVREDPCYSAGAPLAMGMLAGAATGAALGSLMWSPCWF

T**R**T**R**P**L**E - GFP Tag - V

Restriction Sites: SgfI-MluI



View online »

Cloning Scheme:


ACCN: NM_001130035

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130035.2](#)

RefSeq Size: 2027 bp


RefSeq ORF: 570 bp

Locus ID: 58473

UniProt ID: [Q9UF11](#)

Cytogenetics: 11q13.4

Product images:



Circular map for RG225204