

Product datasheet for **RG225242**

Vgl4 (VGLL4) (NM_001128220) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vgl4 (VGLL4) (NM_001128220) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: VGLL4
Synonyms: VGL-4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG225242 representing NM_001128220
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATTAAGTGAGGAACAAGACTGCCAATGGAGACTGCCGCAGAGACCCCGGGAGCGGAGCCGCAGCC
 CCATCGAGCGCGCTGTGGCCCCACCATGAGCCTGCACGGCAGCCACCTGTACACCTCCCTCCCAAGCCT
 TGGCCTGGAGCAGCCCTCGCACTGACCAAGAACAGCCTGGACGCCAGCAGGCCAGCCGGCTCTCGCCC
 AACTGACCCCGGGGAGCGGCAGCAGAACCGGCCCTCCGTGATCACCTGTGCCTCGGCTGGCGCCCGCA
 ACTGCAACCTCTCGCACTGCCCATCGGCACAGCGGCTGTGCCGCGCCCGGGCCTGCCAGCTACCGGAG
 GCCACCGAGCGCTGCCACCACCTGTGACCCCGTGGTGGAGGAGCATTTCGCAGGAGCCTGGGCAAGAAT
 TACAAGGAGCCCGAGCCGGCACCAACTCCGTGTCCATCACGGGCTCCGTGGACGACCACTTTGCCAAAG
 CTCTGGGTGACACGTGGCTCCAGATCAAAGCGGCCAAGGACGGAGCATCCAGCAGCCCTGAGTCCGCCTC
 TCGCAGGGGCCAGCCCGCAGCCCTCTGCCACATGGTCAGCCACAGTCACTCCCTCTGTGGTCTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG225242 representing NM_001128220
 Red=Cloning site Green=Tags(s)

MIKVRNKTANGDCRRDPRERSRSPIERAVPTMSLHGSHLYTSLPSLGLAQPLALTKNSLDASRPAGLSP
 TLTPGERQQNRPSVITCASAGARNLNSHCPIAHSGCAAPGPASYRRPPSAATTCDPVVEEHFRRSLGKN
 YKEPEPAPNSVITGSVDDHFAKALGDTWLQIKAAKDGASSPESASRRQPASPSAHMVSHSHSPSVVS

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001128220

ORF Size: 630 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001128220.3](#)

RefSeq Size: 3182 bp

RefSeq ORF: 633 bp

Locus ID: 9686

UniProt ID: [Q14135](#)

Cytogenetics: 3p25.3-p25.2

Protein Families: Transcription Factors

Gene Summary: May act as a specific coactivator for the mammalian TEFs.[UniProtKB/Swiss-Prot Function]