

Product datasheet for **RG225254**

CDV3 (NM_001134422) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDV3 (NM_001134422) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CDV3
Synonyms:	H41
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG225254 representing NM_001134422 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGAGACGGAGGAGCGGAGCCTGGACAACCTCTTTGCCAAGAGGGACAAGAAGAAGAAGGAGC
GGAGCAACCGGGCGGCGAGTGCCGCGGGCGCAGCGGGCAGCGCCGGCGGAAGCAGTGAGCCGGGTGC
GGCGGGCGGGCGGGCGGGGACCCGGCCGGGTGACGGCGGGACCGCCAGCGCGGGGCTCGGGC
CCAGGGGCGCCACCAAGGCTGTGACGAAGGACGAAGATGAATGGAAAGAATTGGAGCAAAAAGAGTTG
ATTACAGCGGCCCTCAGGGTTCAGGCAATGCAAATAAGCAGTGAAAAGGAAGAAGACGATAATGAAAAGAG
ACAAGATCCAGGTGATAACTGGGAAGAAGTGGAGGTGGTGGTGGAGGTATGGAAAAATCTTCAGGTCCC
TGGAAATAAACAGCTCCAGTACAAGCACCTCCTGCTCCAGTAATTGTTACAGAAACCCAGAACAGCGA
TGACTAGTGGTGTGTATAGGCCTCCTGGGGCCAGGTAAACCACAACAAGGAAAACACCACAAGGACCACC
AGAAATCTACAGTGATACACAGTTCCTATCCCTGCAGTCAACTGCCAAGCATGTAGAAAGCCGGAACAGG
TACTTAAA

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG225254 representing NM_001134422
Red=Cloning site Green=Tags(s)

MAETEERSLDNFFAKRDKKKKKERSNRAASAAGAAGSAGGSSGAAGAAGGGAGAGTRPGDGGTASAGAAG
 PGAATKAVTKDEDEWKELEQKEVDYSGLRVQAMQISSEKEEDDNEKRQDPGDNWEEGGGGGMEKSSGP
 WNK TAPVQAPPAPVIVTETPEPAMTSGVYRPPGARLTTTRKTPQGPPEIYSDTQFPSLQSTAKHVESRNR
 YLK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001134422

ORF Size: 639 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134422.2](#)

RefSeq Size: 4213 bp


RefSeq ORF: 642 bp

Locus ID: 55573

UniProt ID: [Q9UKY7](#)

Cytogenetics: 3q22.1

Product images:



Circular map for RG225254