

Product datasheet for **RG225290**

AZI2 (NM_001134433) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: AZI2 (NM_001134433) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: AZI2
Synonyms: AZ2; NAP1; TILP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG225290 representing NM_001134433
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATGCACTGGTAGAAGATGATATCTGTATTCTGAATCATGAAAAAGCCATAAGAGAGATACAGTGA
CTCCAGTTTCAATATATTCAGGAGATGAATCTGTTGCTTCCCATTTTCTTGTCTACTGCATATGAAGA
CATCAAAAACGACTTAAGGATTCAGAGAAAGAGAAGCTTTTGTAAAGAAGAGAATAAGATTTTGGAA
GAAAAGCTAATAGCTCGATTTGAAGAAGAAACAAGTCCGTGGGACGAGAACAAGTAAATAAGGCCTATC
ATGCATATCGAGAGGTTTGCATTGATAGAGATAATTTGAAGAGCAAAGTGGACAAAATGAATAAAGACAA
CTCTGAATCTTTGAAAGTATTGAATGAGCAGCTACAATCTAAAGAAGTAGAACTCCTCCAGCTGAGGACA
GAGGTGGAAGCTCAGCAGGTGATGAGGAATTTAAATCCACCTTCATCAAAGTGGGAGGTGAAAAAGTTGA
GCTGTGACCTGAAGATCCATGGTTTGAACAAGAGCTGGAAGTATGAGGAAAGAATGTAGCGATCTCAA
AATAGAACTACAGAAAGCCAAACAAACGGATCCATATCAGGAAGACAATCTGAAGAGCAGAGATCTCCAA
AACTAAGCATTCAAGAGTATCTGTAGCACAGCCTGGAGTGCAGTGGCGCAATCATGGGATAGCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG225290 representing NM_001134433
 Red=Cloning site Green=Tags(s)

MDALVEDDICILNHEKAHKRDTVTPVSIYSGDESVAHFALVTAYEDIKKRLKDSEKENSLLKKRIRFLE
 EKLIARFEEETSSVGREQNKAYHAYREVCIDRDNLKSKLDKMKNKDNSESLKVLNEQLQSKEVELLQLRT
 EVETQQVMRNLNPPSSNWEVEKLSCDLKIHGLEQELELMRKECSDLKIELQKAKQTDPTYQEDNLKSRDLQ
 KLSISRVSVAQPGVQWRNHGIA


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001134433

ORF Size: 696 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001134433.2
RefSeq Size:	1456 bp
RefSeq ORF:	699 bp
Locus ID:	64343
UniProt ID:	Q9H6S1
Cytogenetics:	3p24.1
Protein Pathways:	RIG-I-like receptor signaling pathway
Gene Summary:	AZI2, or NAP1, contributes to the activation of NFKB (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK (TBK1; MIM 604834) (Fujita et al., 2003 [PubMed 14560022]).[supplied by OMIM, Mar 2008]