

Product datasheet for **RG225309**

Myeloid leukemia factor 1 (MLF1) (NM_001130157) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Myeloid leukemia factor 1 (MLF1) (NM_001130157) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MLF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG225309 representing NM_001130157 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCCAGGATGCTGAACAGCAGTTTTGAGGATGACCCCTTCTTCTCTGAGTCCATTCTTGCACACCGAG
AAAATATGCGACAGATGATAAGAAGTTTTCTGAACCCTTTGAAGAGACTTGCTCAGTATCTCTGATGG
TAGAGGGAGAGCTCATAATCGTAGAGGACATAATGATGGTGAAGATTCTTGACTCATACAGATGTCAGC
TCTTTCCAGACAATGGACCAAATGGTGTCAAATATGAGAACTATATGCAGAAATTAGAAAGAAACTTCG
GTCAACTTTCAGTGGATCCAAATGGACATTCATTTTGTCTTCCCTCAGTTATGACTTATTCCAAAATAGG
AGATGAACCGCCAAAGGTTTTTCAGGCCTCAACTCAAACTCGTCGAGCTCCAGGAGGAATAAAGGAAACC
AGGAAAGCAATGAGAGATTCTGACAGTGGACTAGAAAAATGGCTATTGGTCATCATATCCATGACCGAG
CTCATGTCATTAAGTCAAAGAACAAGAAGACTGGAGATGAAGAGGTCAACCAGGAGTTCATCAATAT
GAATGAAAGTGTGCTCATGCTTTTGTGAGGAGTGGCAAAGTGAAGTTTGAAGTACAAACCAGGACGA
CACAATCTAGGAAACACTAGAATGAGAAGTGTGGCCATGAGAATCCTGGCTCCCAGAGACTTAAAGAA
GGGAGAAACCTCAACAAAGTCCAGCCATTGAACATGGAAGGAGATCAAATGTTTTGGGGACAACTCCA
CATCAAAGGCTCATCTGTGAAAAGCAACAAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG225309 representing NM_001130157
Red=Cloning site Green=Tags(s)

MFRMLNSSFEDDPFFSEILAHRENMRQMIRSFSEPFGRDLLSISDGRGRAHNRGRHNDGEDSLTHTDVS
 SFQTMQMVSNMARNYMQKLERNFQQLSVDPNHSHFCSSSVMTYSKIGDEPPKVFQASTQTRRAPGGIKET
 RKAMRSDSDGLEKMAIGHHHIHDRAHVIKSKNKKTGDEEVNQEFINMNESDAHAFDEEWQSEVLKYKPGR
 HNLGNTRMRSVGHENPGSRELKRREKPPQSPAIEHGRRSNVLGDKLHIKGSVSKSNKK


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001130157

ORF Size: 807 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130157.2 , NP_001123629.1
RefSeq Size:	2336 bp
RefSeq ORF:	732 bp
Locus ID:	4291
UniProt ID:	P58340
Cytogenetics:	3q25.32
Protein Families:	Druggable Genome
Gene Summary:	This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoetic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010]