

Product datasheet for **RG225349**

AIF (AIFM1) (NM_001130846) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AIF (AIFM1) (NM_001130846) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	AIFM1
Synonyms:	AIF; AUNX1; CMT2D; CMTX4; COWCK; COXPD6; DFNX5; NADMR; NAMSD; PDCD8; SEMDHL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG225349 representing NM_001130846
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTCCGGTGTGGAGGCTGGCGCGGGTCTTTGAAGCAGAAGCTGGTGCCCTTGGTGCGGACCGTGT
 GCGTCCGAAGCCGAGGCAGAGGAACCGCTCCAGTTGTGCAGTCTCATCACCTAGGATCCCCTCTAG
 ATCACTAGCATCTACAGGTGCTTCTGGGAAAGATGGCAGCAACCTAGTGTACTTCTTAATTGTAGGAGCA
 ACAGTCACTGGGGCAGGAGTTTATTATGCCTACAAGACTATGAAAGAGGACGAAAAAGATAACAATGAAA
 GAATTTACAGGTTAGGGCTGACACCAGAACAGAAAAGGCGCGTTATCTGCTTCAAGGAGAGA
 GGAAGTTCCTCAAGACAAGGCGCAAGTCAATGTTCTTTCTGCTAATTGGTGGAGGCACAGCTGCTTTT
 GCTGCAGCCAGATCCATCCGGGCTCGGGATCCTGGGGCCAGGGTACTGATTGTATCTGAAGATCCTGAGC
 TGCCGTACATGCGACCTCTCTTCAAAGAAGTGTGGTTTTTCAGATGACCCAAATGTCACAAAGACT
 GCGATTCAAACAGTGAATGGAAAAGAGAGAAGCATATATTTCCAGCCACCTTCTTCTATGTCTCTGCT
 CAGGACCTGCCTCATATTGAGAAATGGTGGTGTGGCTGTCCTCACTGGGAAGAAGGTAGTACAGCTGGATG
 TGAGAGACAACATGGTGAACCTAATGATGGCTCTCAAATAACCTATGAAAAGTGCTTGATTGCAACAGG
 AGGTACTCCAAGAAGTCTGTCTGCCATTGATAGGGCTGGAGCAGAGGTGAAGAGTAGAACACGCTTTTC
 AGAAAAGATTGGAGACTTTAGAAGCTTGGAGAAGATTTACGGGAAGTCAAATCAATTACGATTATCGGTG
 GGGGCTTCTTGGTAGCGAAGTGGCCTGTGCTCTTGGCAGAAAGGCTCGAGCCTTGGGCACAGAAGTGT
 TCAACTCTTCCCCGAGAAAGGAAATATGGGAAAGATCCTCCCGAATACCTCAGCAACTGGACATGGAA
 AAAGTCAGACGAGAGGGGGTTAAGGTGATGCCCAATGCTATTGTGCAATCCGTTGGAGTCAGCAGTGGCA
 AGTTACTTATCAAGCTGAAAGACGGCAGGAAGGTAGAAACTGACCACATAGTGGCAGCTGTGGCCTGGA
 GCCCAATGTTGAGTTGGCCAAGACTGGTGGCTGGAAATAGACTCAGATTTTGGTGGCTTCCGGGTAAT
 GCAGAGCTACAAGCACGCTCTAACATCTGGGTGGCAGGAGATGCTGCATGCTTCTACGATATAAAGTTGG
 GAAGGAGGCGGGTAGAGCACCATGATCACGCTGTTGTGAGTGGAAAGATTGGCTGGAGAAAAATGACTGG
 AGCTGCTAAGCCGACTGGCATCAGTCAATGTTCTGGAGTGAATTTGGGCCCGATGTTGGCTATGAAGCT
 ATTGGTCTTGTGGACAGTAGTTTCCCACAGTTGGTGTGTTTTGCAAAGCAACTGCACAAGACAACCCCA
 AATCTGCCACAGAGCAGTCAGGAAGTATCCGATCAGAGAGTGAGACAGAGTCCGAGGCCCTCAGAAAT
 TACTATTCCTCCCAGCACCCGGCAGTTCACAGGCTCCCGTCCAGGGGGAGGACTACGGCAAAGGTGTC
 ATCTTCTACCTCAGGGACAAAGTGGTCTGTTGGGATTGTGCTATGGAACATCTTAAACCGAATGCCAATAG
 CAAGGAAGATCATTAAAGACGGTGGCAGCATGAAGATCTCAATGAAGTAGCCAAACTATTCAACATTCA
 TGAAGAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG225349 representing NM_001130846
 Red=Cloning site Green=Tags(s)

MFRCGGLAAGALKQKLVPLVRTVCVRSRQRNRLPVVQSHHLGSPSRSLASTGASGKDGSNLVYFLIVGA
 TVTGAGVYYAYKTMKEDEKRYNERISGLGLTPEQKQKKAALSASEGEEVPQDKAPSHVPLLIGGGTAAF
 AAARSIRARDPGARVLIVSEPELPMRPLSKELWFSDPNVTKLRFKQWNGKERSIYFQPPSFYVSA
 QDLPHIENGGVAVLTGKKVVQLDVRDNMVKLNDGSQITYEKCLIAATGGTPRSLAIDRAGAEVKSRTTLF
 RKIGDFRSLEKISREVKSIITIGGGFLGSELACALGRKARALGTEVIQLFPEKGNMGKILPEYLSNWTME
 KVRREGVKVMPNAIVQSVGVSSGKLLIKLKDGRKVETDHIVAAVGLEPNVELAKTGGLEIDSDFGGFRVN
 AELQARSNIWVAGDAACFYDIKLRRRRVEHHDHAVVSGRLAGENMTGAAKPYWHQSMFWSDLGPDVGYEA
 IGLVDSLPTVGVFAKATAQDNPKSATEQSGTGIRSESETESEASEITIPPSTPAVPQAPVQGEDYKGV
 IFYLDRKVVVGVILWNIFNRMPIARKIIKDGEQHEDLNEVAKLFNIHED

TRTRPLE – GFP Tag – V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001130846
 ORF Size: 1830 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130846.2 , NP_001124318.1
RefSeq Size:	1946 bp
RefSeq ORF:	825 bp
Locus ID:	9131
UniProt ID:	O95831
Cytogenetics:	Xq26.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis
Gene Summary:	This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells, and it is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it affects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Mutations in this gene cause combined oxidative phosphorylation deficiency 6 (COXPD6), a severe mitochondrial encephalomyopathy, as well as Cowchock syndrome, also known as X-linked recessive Charcot-Marie-Tooth disease-4 (CMTX-4), a disorder resulting in neuropathy, and axonal and motor-sensory defects with deafness and cognitive disability. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 10. [provided by RefSeq, Aug 2015]