

Product datasheet for **RG225375**

LDHA (NM_001135239) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | LDHA (NM_001135239) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | LDHA |
| Synonyms: | GSD11; HEL-S-133P; LDHM; PIG19 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >RG225375 representing NM_001135239 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAACTCTAAAGGATCAGCTGATTTATAATCTTCTAAAGGAAGAACAGACCCCCAGAATAAGATTA
CAGTTGTGGGGTTGGTGTCTTGGCATGGCCTGTGCCATCAGTATCTTAATGAAGGACTTGGCAGATGA
ACTTGCTCTTGTGATGTCATCGAAGACAAATTGAAGGGAGAGATGATGGATCTCCAACATGGCAGCCTT
TTCCTTAGAACACCAAAGATTGTCTCTGGCAAAGTGGATATCTTGACCTACGTGGCTTGAAGATAAGTG
GTTTTCCAAAAACCGTGTATTGGAAGTGGTTGCAATCTGGATTCAGCCCGATTCCGTTACCTGATGGG
GGAAAGGCTGGGAGTTCACCCATTAAGCTGTATGGGTGGTCCCTGGGGAACATGGAGATTCCAGTGTG
CCTGTATGGAGTGAATGAATGTTGCTGGTGTCTCTCTGAAGACTCTGCACCCAGATTTAGGGACTGATA
AAGATAAGGAACAGTGGAAAGAGGTTCAACAAGCAGGTGGTTGAGAGTGCTTATGAGGTGATCAAACCTCAA
AGGCTACACATCCTGGGCTATTGGACTCTCTGTAGCAGATTTGGCAGAGAGTATAATGAAGAATCTTAGG
CGGGTGCACCCAGTTCCACCATGATTAAGGGTCTTTACGGAATAAAGGATGATGTCTTCCTTAGTGTT
CTTGCAATTTGGGACAGAATGGAATCTCAGACCTGTGAAGGTGACTCTGACTTCTGAGGAAGAGGCCCG
TTTGAAGAAGAGTGCAGATACACTTTGGGGATCCAAAAGGAGCTGCAATTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG225375 representing NM_001135239
 Red=Cloning site Green=Tags(s)

MATLKDQLIYNLLKKEEQTPQNKITVVGVGAVGMACAISILMKDLADELALVDVIEDKLGEMMDLQHGSL
 FLRTPKIVSGKVDILTYVAWKISGFPKNRVIGSGCNLDSARFRYLMGERLGVHPLSCHGWVLGEHGDSSV
 PVWSGMNVAGVSLKTLHPDLGTDKDKQKQKVEVHKQVVESAYEVIKLGKGYTSWAIGLSVADLAESIMKNLR
 RVHPVSTMIKGLYGIKDDVFLSVPCILGQNGISDLVKVTLTSEEEARLKKADTLWGIQKELQF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001135239

ORF Size: 822 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135239.1](#), [NP_001128711.1](#)

RefSeq Size: 2052 bp

RefSeq ORF: 825 bp

Locus ID: 3939

UniProt ID: [P00338](#)


Cytogenetics: 11p15.1

Protein Families: Druggable Genome

Protein Pathways: Cysteine and methionine metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism

Gene Summary: The protein encoded by this gene catalyzes the conversion of L-lactate and NAD to pyruvate and NADH in the final step of anaerobic glycolysis. The protein is found predominantly in muscle tissue and belongs to the lactate dehydrogenase family. Mutations in this gene have been linked to exertional myoglobinuria. Multiple transcript variants encoding different isoforms have been found for this gene. The human genome contains several non-transcribed pseudogenes of this gene. [provided by RefSeq, Sep 2008]

Product images:



Circular map for RG225375