

Product datasheet for **RG225508**

PPP6C (NM_001123355) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP6C (NM_001123355) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PPP6C
Synonyms:	PP6; PP6C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG225508 representing NM_001123355 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCCGCTAGACCTGGACAAGTATGTGAAATAGCGCGGCTGTGCAAGTACCTGCCAGAGAACGACC
TGAAGGTGAGCCCTATTTGCGGGCTCGCCCTTCCGGCTGTGGGGCGCCCGCGGCCCTTCCCTTC
TCCAGGACCCCGCCGGTATTTCAATTTCTTCGTTTTCTAAAGGAGCGGCTATGTGACTACGTTTGTGAC
CTCCTCTTAGAAGAGTCAAATGTTCAAGCAGTATCAACACCAGTAAACAGTGTGGAGATATCCATGGAC
AGTTTTATGACCTTTGTGAAGTTCAGAACTGGAGGTCAGGTTCTGACACAACTACATATTTATGGG
TGATTTTGTAGACAGAGTTACTATAGTTTGGAGACCTTCACTTACCTTCTTGCATTAAGGCTAAATGG
CCTGATCGTATTACACTTTTGCAGGAAATCATGAGAGTAGACAGATAACACAGGCTATGGATTTTATG
ATGAGTGCCAAACCAAATATGGAATGCTAATGCCTGGAGATACTGTACCAAAGTTTTTGACATGCTCAC
AGTAGCAGCTTAAATAGATGAGCAGATTTGTGTGCCATGGTGGTTTATCTCCTGATATCAAAACACTG
GATCAAATTCGAACCATCGAACGGAATCAGGAAATTCCTATAAAGGAGCATTGTGTGATCTGGTTTGGT
CAGATCCTGAAGATGTGGATACCTGGGCTATCAGTCCCGAGGAGCAGGTTGGCTTTTGGAGCAAAGGT
CACAAATGAGTTTGTTCATATCAACAACCTAAAACCTCATCTGCAGAGCACATCAACTAGTGCACGAAGGC
TATAAATTTATGTTTGTGAGAAGCTGGTGCAGTATGGTCTGCTCCTAATTACTGCTATCGTTGTGGAA
ATATTGCTTCGATCATGGTCTTCAAAGATGTAATAACAAGAGAACCAAAGTTATTCCGGGCAGTTCAGAA
TTCAGAACGTGTTATTCTCCAGAACGACAACGCCATATTTCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG225508 representing NM_001123355
 Red=Cloning site Green=Tags(s)

MAPLDDLDKYVEIARLCKYLPENDLKVSPICGLAPSGCGAPAGRPFLSPGPPPVFHFRLFLKERLCDYVCD
 LLLLEESNVQPVSTPVTVCGLDIHQFYDLCELFRGTGGQVPDTNYIFMGDFVDRGYYSLETFTYLLALKAKW
 PDRITLLRGNHESRQITQVYGFYDECQTKYGNANAWRYCTKVFDMLTVAALIDEQILCVHGGGLSPDIKTL
 DQIRTIERNQEIPHKGAFCDLVWSDPEDVDTWAI SPRGAGWLFGAKVTNEFVHINNKLICRAHQLVHEG
 YKFMFDEKLVTVWSAPNYCYRCGNIASIMVFKDVNTREP KLFRAVPDSERVIPRRTTTPYFL


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001123355

ORF Size: 1026 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001123355.1 , NP_001116827.1
RefSeq Size:	4367 bp
RefSeq ORF:	1029 bp
Locus ID:	5537
UniProt ID:	O00743
Cytogenetics:	9q33.3
Protein Families:	Druggable Genome, Phosphatase
Gene Summary:	This gene encodes the catalytic subunit of protein phosphatase, a component of a signaling pathway regulating cell cycle progression. Splice variants encoding different protein isoforms exist. The pseudogene of this gene is located on chromosome X. [provided by RefSeq, Jul 2008]