

Product datasheet for **RG225657**

alpha 1 Antitrypsin (SERPINA1) (NM_001127702) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	alpha 1 Antitrypsin (SERPINA1) (NM_001127702) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SERPINA1
Synonyms:	A1A; A1AT; AAT; alpha1AT; nNIF; PI; PI1; PRO2275
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG225657 representing NM_001127702 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGCTTCTGTCTCGTGGGGCATCCTCCTGCTGGCAGGCCTGTGCTGCCTGGTCCCTGTCTCCCTGG
CTGAGGATCCCCAGGGAGATGCTGCCAGAAAGACAGATACATCCCACCATGATCAGGATCACCCAACCTT
CAACAAGATCACCCCAACCTGGCTGAGTTCGCCTTCAGCCTATACCGCCAGCTGGCACACCAGTCCAAC
AGCACCAATATCTTCTTCCCCAGTGAGCATCGCTACAGCCTTTGCAATGCTCTCCCTGGGGACCAAGG
CTGACACTCACGATGAAATCCTGGAGGGCCTGAATTTCAACCTCACGGAGATTCCGGAGGCTCAGATCCA
TGAAGGCTTCCAGGAACTCCTCCGTACCTCAACCAGCCAGACAGCCAGCTCCAGCTGACCACCGCAAT
GGCCTGTTCCCTCAGCGAGGGCCTGAAGCTAGTGGATAAGTTTTGGAGGATGTTAAAAAGTTGTACCACT
CAGAAGCCTTCACTGTCAACTTCGGGGACACCGAAGAGGCCAAGAAACAGATCAACGATTACGTGGAGAA
GGTACTCAAGGGAAAAATTGTGGATTTGGTCAAGGAGCTTGACAGAGACACAGTTTTGTCTGGTGAAT
TACATCTTCTTTAAAGGCAATGGGAGAGACCTTTGAAGTCAAGGACACCGAGGAAGAGGACTTCCACG
TGGACCAGGTGACCACCGTGAAGGTGCCTATGATGAAGCGTTTAGCATGTTAACATCCAGCACTGTAA
GAAGCTGTCCAGCTGGGTGCTGCTGATGAAATACCTGGCAATGCCACCGCCATCTTCTCTGCCTGAT
GAGGGGAAACTACAGCACCTGGAAAATGAACTCACCCACGATATCATACCAAGTTCCTGGAAAATGAAG
ACAGAAGGTCTGCCAGTTACATTTACCCAAACTGTCCATTACTGGAACCTATGATCTGAAGAGCGTCTT
GGGTCAACTGGGCATCACTAAGGTCTTACGCAATGGGGCTGACCTCTCCGGGGTACAGAGGAGGCCACCC
CTGAAGCTCTCCAAGGCCGTGCATAAGGCTGTGCTGACCATCGACGAGAAAGGGACTGAAGCTGCTGGG
CCATGTTTTAGAGGCCATACCCATGTCTATCCCCCGAGGTCAAGTTCAACAAACCTTTGTCTTCTT
AATGATTGAACAAAATACCAAGTCTCCCCTCTTCATGGGAAAAGTGGTGAATCCCACCAAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG225657 representing NM_001127702
 Red=Cloning site Green=Tags(s)

MPSSVSWGILLLAGLCLLPVSLAEDPQGDAAQKTDTSHHDDHPTFNKITPNLAEFASFSLYRQLAHQSN
 STNIFFSPVSIATAFAMLSLGTKADTHDEILEGLNFNLTEIPEAQIHEGFQELLRTLNPDSQLQLTTGN
 GLFLSEGLKLVDFLEDVKKLYHSEAFVNFVGDTEEAKKQINDYVEKGTQGKIVDLVKELDRDVFALVN
 YIFFKKGWERPFVVKDTEEDDFHVDQVTTVKVPMMKRLGMFNIQHCKLSSWVLLMKYLGNATAIFFLPD
 EGKLQHLENELTHDIITKFLNEDRRSASLHLPKLSITGTYDLKSVLGQLGITKVFVNGADLSGVTEEAP
 LKLSKAVHKAVLTIDEKGTAAAGAMFLEAIPMSIPPEVKFNKPFVFLMIEQNTKSPLFMGKVVNPTQK


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001127702

ORF Size:	1254 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001127702.2
RefSeq Size:	3340 bp
RefSeq ORF:	1257 bp
Locus ID:	5265
UniProt ID:	P01009
Cytogenetics:	14q32.13
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	Complement and coagulation cascades
Gene Summary:	The protein encoded by this gene is a serine protease inhibitor belonging to the serpin superfamily whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. This protein is produced in the liver, the bone marrow, by lymphocytic and monocytic cells in lymphoid tissue, and by the Paneth cells of the gut. Defects in this gene are associated with chronic obstructive pulmonary disease, emphysema, and chronic liver disease. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Aug 2020]